

# Validation of Gamma Spectroscopy, Chemical Separation Alpha Spectrometry, Gas Proportional Counting, and Liquid Scintillation Analytical Data

Effective Date: 4/24/2017Next Review Date: 4/24/2020

**Hazard Class:**  Low  Moderate  High/Complex  
**Usage Mode:**  Reference  UET  Both UET & Reference

The Responsible Manager has determined that the following organizations' review is required for initial procedure release as well as subsequent major revisions. Review documentation is contained in the Document History File.

Technical Leads

Quality Assurance

Classification Review:  Unclassified  UCNI  Classified

Diana Hollis / 111125 / /s/ Diana Hollis / 4/18/2017  
 Name (print) Z# Signature Date

Responsible Manager, Division and Title

Nita Patel / 153003 / /s/ Nita Patel / 4/20/2017  
 Name (print) Z# Signature Date

*Reference*

---

**REVISION HISTORY**

Document No./Revision No.	Issue Date	Action	Description
OIO-TP-5166, Rev. 0	8/2/2016	New	New Document, changed Doc # from SOP-5166 to OIO-TP-5166.
ER-AP-20314, R0	4/24/2017	Major Revision	Revised to remove NNSA Model Validation

Reference

---

**TABLE OF CONTENTS**

<u>Section</u>	<u>Page</u>
TITLE PAGE.....	1
REVISION HISTORY .....	2
TABLE OF CONTENTS .....	3
1. PURPOSE .....	4
2. SCOPE.....	4
3. BACKGROUND.....	4
4. PRECAUTIONS.....	4
5. PREREQUISITE ACTIONS.....	5
6. PERFORMANCE.....	5
6.1 Validation Process .....	5
6.2 Analyte Quantitation .....	6
6.4 Holding Times and Sample Preservation .....	7
6.5 Initial and Continuing Calibration.....	8
6.6 Carrier.....	8
6.7 Tracer.....	8
6.8 Blanks .....	9
7. RECORDS.....	10
8. REFERENCES .....	10
9. ATTACHMENTS .....	10
<u>Attachments</u>	
Attachment 1, Data Validation Cover Sheet .....	11
Attachment 2, Radiochemistry Analytical Data Validation Checklist .....	12

*Reference*

---

**1. PURPOSE**

This procedure establishes guidance for the qualification of radionuclide analytical data. This document is intended to assist in the technical review of analytical data generated by environmental laboratories.

**2. SCOPE**

Qualification of data is the product of data validation, analytical laboratory analysis, and focused validation that describe validation anomalies and their consequences. Radiochemistry measurements include:

- gamma-emitting isotopes by gamma spectroscopy;
- alpha-emitting isotopes (americium-241; uranium-234, -235, and -238; thorium-230, -232, and -234; and plutonium-238 and -239/-240) by chemical separation alpha spectrometry;
- strontium-90 by gas-proportional counting (GPC);
- gross-alpha and -beta analyses by GPC; and
- tritium by liquid scintillation.

**3. BACKGROUND**

Data qualifiers and reason codes are assigned to analytical results from radiochemical analyses according to the specifications in this method-specific procedure. This procedure conforms to the requirements of U.S. Environmental Protection Agency (EPA) methodologies. Data qualifiers and reason codes are assigned according to the specifications in this method specific procedure.

**4. PRECAUTIONS**

Nothing in this procedure precludes the data validator from going beyond the minimum requirements specified herein. If additional directions are required, the data validator shall reference EPA method-specific guidelines. Implementation of this procedure may be followed by a more focused and data use-specific evaluation of the data by the project chemist, especially if the implementation of this procedure indicates the data may contain technical deficiencies.

*Reference*

---

**5. PREREQUISITE ACTIONS**

Data Validators must:

- Possess a minimum of a Bachelor's degree in chemistry or one of the physical sciences and either two (2) years of experience in generating analytical data in an environmental analytical laboratory or two (2) years of experience in data validation.
- Complete Attachment 1, Data Validation Cover Sheet, and Attachment 2, Radiochemistry Analytical Data Validation Checklist, during data validation.

**6. PERFORMANCE**

**6.1 Validation Process**

EIM applies a subset of qualifiers described in this procedure to analytical data using auto-validation subroutines. EIM auto-validation applies qualification to analytical records using tests listed in Attachment 2 that have a Valid Reason Description containing "(AV)". When the project leader requests a focused validation the assigned data validator completes the following steps to assess all potential analytical data qualification:

- [1] **REVIEW** the qualifiers assigned during EIM auto-validation to verify that qualifiers were assigned consistently with this procedure. If auto-validation qualification is found to be inconsistent with this procedure then the validator initiates a change request using ER-AP-20304, Change Control for Data in the Environmental Information Management (EIM) Database.
- [2] **PRINT** Attachment 1 and **REVIEW** the data package for potential qualification using Attachment 2.
- [3] **NOTE** conditions causing recommendation for qualification and options for qualification.
- [4] **COMPLETE** Attachment 1 and **FORWARD** to the project leader with conditions and options.

*Reference*

---

**6.1 Validation Process (continued)**

The project leader is the responsible party for making the decision of record if validation qualifiers should be assigned and EIM validation records updated. This record of decision is added to comments section of Attachment 1.

Once the decision of record has been made, Attachment 1 is sent to the Sample Management Office (SMO) staff. The SMO staff re-print the data validation record from EIM and add Attachment 1 that includes the record of decision to the final records package.

**6.2 Analyte Quantitation**

The assignment of the detection status to analytical measurements is the first step of analytical data validation. Most validation qualifiers and validation reason codes are applied based on the measurement's initial detection status. Results that are less than the report minimum detectable activity (MDA) are qualified as nondetect with the U validation qualifier and R5 validation reason code. Results greater than or equal to the MDL are qualified as detected with the NQ validation qualifier.

<b>Criteria</b>	<b>Validation Qualifier</b>	<b>Validation Reason Code</b>
Target analyte result is < MDA; a nondetect	U	R5
Target analyte result is $\geq$ MDA; a detect	NQ	NQ

Since a result can have only one validation qualifier and one validation reason code the sequencing of validation steps is important. Analyte quantitation occurs first, then analyte identification, because most other validation functions depend on the correct identification and quantitation of the analytical parameter. When two or more qualifiers can be applied to a record, the qualifier representing the more severe consequence to data usability supersedes the qualifier with less severe consequence. The R validation qualifier has the greatest impact on data usability and supersedes other validation qualifiers.

*Reference*

---

**6.2 Analyte Quantitation (continued)**

<b>Order Of Severity</b>	<b>Validation Qualifier</b>	<b>Description</b>
1	R	The reported sample result is classified as rejected due to serious non-compliances regarding quality control acceptance criteria. The presence or absence of the analyte cannot be verified.
2	UJ	The analyte is classified as not detected, with an expectation that the reported result is more uncertain than usual.
3	U	The analyte is classified as not detected.
4	J	The analyte is classified as detected but the reported concentration value is expected to be more uncertain than usual.
5	NQ	No validation qualifier flag is associated with this result, and the analyte is classified as detected.

LANL project chemists may identify quality deficiencies in analytical results affecting analyte quantitation. These deficiencies can include analytical results with detection limits elevated above project data-quality objectives, concentrations above the calibration range of the instrument or method, results exhibiting carryover or detector contamination, spectral interference from another radioisotope, and other quality deficiencies. The reason code of R19 is applied to affected records by the project chemist to identify these quality deficiencies when they are identified.

**6.3 Analyte Identification**

The identification of an analytical parameter is the second step of analytical data validation. Identification of radioisotopes depends upon the energy signature of radioactive decay, isotope luminescence, and/or chemical separation. When counting uncertainty, peak-width, or abundance do not meet criteria the R5a reason code is applied to affected parameters.

**6.4 Holding Times and Sample Preservation**

Sample handling requirements are specified to ensure integrity and defensibility of analytical measurements. Samples are to be prepared and analyzed within specified time limits. Samples are also preserved chemically with the addition of acids and physically by controlling temperature. When sample handling requirements are not met the R9 series of reason codes are applied to affected samples.

*Reference*

---

**6.5 Initial and Continuing Calibration**

Calibration is performed to ensure that the instrument is performing within specifications. The calibration establishes the energy/channel relationship, counting efficiency, and energy resolution. Instrument performance checks are performed periodically to confirm the instrument remains within operating specifications. When a calibration has expired the R7 series of reason codes are applied to affected analytes in all samples analyzed since the expiration of the calibration for that instrument. When an instrument performance check is not performed at the analytical laboratory's established frequency or is not within limits the R7 series of reason codes is applied to affected analytes in all samples analyzed after the unacceptable instrument performance check to the next acceptable instrument performance check for that instrument.

**6.6 Carrier**

Carriers are compounds not normally found in the environment, but which have quantitation limits and energies similar to the radioisotopes of interest in a sample. Carriers are added to samples, standards, and QC samples to determine the effectiveness of analyte quantitation and chemical separation. Sample results are not adjusted based on carrier recoveries. When carrier recovery criteria are not met the R1 series of reason codes are applied to affected samples.

**6.7 Tracer**

Tracers are compounds not normally found in the environment, but which are easily measurable. They are added to samples, standards, and QC samples to compensate for fluctuations in the analytical system. Sample results are quantitated or adjusted by the relative response of associated tracers. When tracer criteria are not met the R3 series of reason codes are applied to the affected sample.



*Reference*

---

**6.8 Blanks**

The Method Blank is an analyte-free matrix that is prepared and analyzed in the laboratory with the samples. The method blank determines contamination from the analytical processes. Method blanks are prepared with every preparation batch. If more than one method blank is associated with a given sample, qualification is based upon a comparison with the associated blank having the highest concentration of the parameter. When method blank criteria are not met the R4 series of reason codes are applied affected samples.

The Field Blank is an analyte-free matrix opened to the atmosphere at the time of sample collection. Field blanks are used to determine if atmospheric conditions resulted in contamination of samples during sample collection. Samples collected the same day as the field blank that does not meet blank criteria are qualified with the R4 series of reason codes.

The Equipment Blank is an analyte-free matrix poured over or through sample collection equipment. Equipment blanks are used to determine the cleaning effectiveness of sampling equipment between samples. Samples collected using the same tools as the equipment blank that does not meet blank criteria are qualified with the R4 series of reason codes.

**6.9 Laboratory Control and Matrix Spike Samples**

The laboratory control sample is created by adding known amounts of parameters of interest to an aliquot of a blank matrix. The laboratory control sample is used to evaluate the effect of the analytical process of the recovery of analytes. When laboratory control sample criteria are not met the R12 series of reason codes are applied to all associated samples.

The matrix spike is created by adding known amounts of parameters of interest to an aliquot of a sample matrix. The matrix spike is used to evaluate the effect of the sample matrix on the recovery of analytes. When matrix spike criteria are not met the R6 series of reason codes are applied to all associated samples.

**6.10 Sample Duplicate**

Field duplicate samples are collected from the same material at the same time as the primary sample. The relative percent difference between the results of the parent sample and the field duplicate sample is used to determine the field and laboratory precision of the analytical measurement. When field duplicate precision criteria are not met the appropriate R10 reason code is applied to the parent sample.

*Reference*

**7. RECORDS**

Records generated by this procedure will be submitted to the Environmental Protection Records Management Office for document management in accordance with Institutional Records Management Procedure, P1020-1 and EP-AP-10003, Records Management.

- Completed Data Validation Cover Sheets
- Completed Radiochemistry Analytical Data Validation Checklists

**8. REFERENCES**

EP-AP-10003, Records Management

ER-AP-20304, Change Control for Data in the Environmental Information Management (EIM) Database

P1020-1, Laboratory Records Management

**9. ATTACHMENTS**

Attachment 1: Data Validation Cover Sheet

Attachment 2: Radiochemistry Analytical Data Validation Checklist

Reference

**ATTACHMENT 1**

Page 1 of 1

**Data Validation Cover Sheet**

<b>Section I.</b>							
Request Number: _____		Validation Date: _____		Lab Code: _____			
Contract Laboratory Name: _____							
Validator: _____				Organization: _____			
Analytical Suite (Check All That Apply):							
<input type="checkbox"/> TPH-GRO	<input type="checkbox"/> High Explosives	<input type="checkbox"/> Dioxin Furans	<input type="checkbox"/> LCMSMS Perchlorates				
<input type="checkbox"/> TPH-DRO	<input type="checkbox"/> Metals & Cyanide	<input type="checkbox"/> PCB Congeners	<input type="checkbox"/> Organochlorine Pesticides/Polychlorinated Biphenyls				
<input type="checkbox"/> General Chemistry	<input type="checkbox"/> Radiochemistry	<input type="checkbox"/> LCMSMS High Explosives					
<input type="checkbox"/> Other (Describe): _____							
<b>Section II. Completeness Check</b>							
YES	NO	N/A	(check one)	YES	NO	N/A	(check one)
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	1. Chain-Of-Custody Form(S)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	6. Raw/BSS Data
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	2. Case Narrative	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	7. Quality Control Forms
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	3. Sample Result Forms	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	8. Quantitation Reports
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	4. Sample Chromatograms	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	9. TICS Forms
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	5. Standard Chromatograms	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	10. TICS Mass Spectra
Comments/problems noted (include information about requests for further information submitted to the contract laboratory and agreed-upon date of resolution and contract laboratory point of contact):							
Validator's Signature: _____ Date: _____							
ER-AP-20314, R0				Los Alamos Environmental Safety & Health			
				(Attach additional comment sheets as necessary)			

Reference

**ATTACHMENT 2**

Page 1 of 4

**Radiochemistry Analytical Data Validation Checklist**

Yes No N/A				Assign Qualifier Listed Below If Criterion = Yes	
(Check One)				Non-detected Analyte	Detected Analyte
<b>Holding Times and Sample Preservation</b>					
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	1. The affected radionuclide is regarded as rejected because the analytical holding time was exceeded, 28 days C-14, I-129, I-131 in water, 180 days all others. (AV)	R, R9b	J-, R9b
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	2. C-14, I-129, I-131 temperature > 10°C upon receipt at the laboratory.	UJ, R9c	J-, R9c
<b>Instrument Calibration</b>					
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	3. Instrument calibration or instrument performance check is expired or deficient.	UJ, R7	J, R7
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	4. Instrument calibration not performed.	R, R7a	R, R7a
<b>Blanks</b>					
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	5. The sample result is $\leq 5$ times the concentration of the related analyte in the method blank. (AV)	N/A	U, R4
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	6. The affected analytes are considered estimated and biased high because this analyte was identified in the method blank but was >5 times. (AV)	N/A	J+, R4a
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	7. The sample result is $\leq 5$ times the concentration of the related analyte in the trip blank, rinsate blank, or equipment blank.	N/A	U, R4d
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	8. Required method blank information is missing. Data may not be acceptable for use. Contact the SMO or external laboratory for information.	R, R4e	R, R4e

Reference

**ATTACHMENT 2**

Page 2 of 4

Yes No N/A (Check One)				Assign Qualifier Listed Below If Criterion = Yes	
				Non-detected Analyte	Detected Analyte
<b>Tracer/Carrier</b>					
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	9 . The tracer/carrier is <10%R. Follow the external laboratory limits located within the associated data package. Tracer %R is not applicable for gamma spectroscopy.	R, R3	R, R3
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	10. The tracer/carrier is less than the lower acceptance level (LAL) but ≥10%R. Follow the external laboratory limits located within the associated data package. tracer/carrier %R is not applicable for gamma spectroscopy.	UJ, R3a	J-, R3a
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	11. The tracer/carrier %R value is greater than the upper acceptance limit (UAL). Follow the external laboratory limits located within the associated data package. tracer/carrier %R is not applicable for gamma spectroscopy.	N/A	J+, R3b
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12. Required tracer/carrier information is missing. Data may not be acceptable for use. Contact the SMO or external laboratory for information. tracer/carrier %R is not applicable for gamma spectroscopy.	R, R3d	R, R3d
<b>Laboratory Control Samples</b>					
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	13. The LCS %R was <10%. Follow the external laboratory limits located within the associated data package. (AV)	R, R12	R, R12
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	14. The LCS %R was < the LAL but >10%. Follow the external laboratory limits located within the associated data package. (AV)	UJ, R12a	J-, R12a
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15. The LCS %R was > the UAL. Follow the external laboratory limits located within the associated data package. (AV)	N/A	J+, R12b
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16. The LCS documentation is missing. Data may not be acceptable for use. Contact the SMO or external laboratory for information.	R, R12c	R, R12c

Reference

**ATTACHMENT 2**

Page 3 of 4

Yes No N/A  (Check One)				Assign Qualifier Listed Below If Criterion = Yes	
				Non-detected Analyte	Detected Analyte
<b>Matrix Spike</b>					
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	17. The associated matrix spike recovery was <10%. Follow the external laboratory limits. MS/MSD is not applicable to gamma spectroscopy.	UJ, R6a	R, R6a
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	18. The associated matrix spike recovery was less than the LAL but greater than 10%. Follow the external laboratory limits. MS/MSD is not applicable to gamma spectroscopy.	UJ, R6b	J-, R6b
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	19. Required matrix spike information is missing. Data may not be acceptable for use. Contact the SMO or external laboratory for information. If LCS information is present, do not reject. Qualify data based on LCS information. MS/MSD is not applicable to gamma spectroscopy.	R, R6c	R, R6c
<b>Duplicate</b>					
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20. Associated duplicate sample has DER or RER greater than the analytical laboratory's acceptance limits. (AV)	R, R10	J, R10
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	21. The duplicate sample was not prepared and/or analyzed with the samples for unspecified reasons. The duplicate information is missing. Data may not be acceptable for use. Contact the SMO or external laboratory for information.	R, R10d	R, R10d

Reference

**ATTACHMENT 2**

Page 4 of 4

<b>Analyte Quantitation</b>					
			22. The non-detected analytes have elevated detection limits and may not meet project data-quality objectives because the sample was diluted without any target analytes identified as a result of matrix interference. Reject non-detected results if the analytical laboratory cannot provide proof for matrix interference.	UJ, R, R15	NA
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	23. The results for the affected analytes are considered not detected (U) because the associated sample concentration was less than or equal to the MDC. (AV)	N/A	U, R5
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24. The analytical laboratory qualified the result with UI. The result is greater than MDA but is a non-identified nuclide.	U, R5a	U, R5a
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	25. Interferences prevent positive identification of the detected analyte or misidentification of the non-detected analyte.	R, R5a	R, R5a
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	26. The MDC and/or TPU documentation is missing. Data may not be acceptable for use. Contact the SMO or external laboratory for information.	R, R5b	R, R5b
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	27. The LANL project chemist identified quality deficiencies in the reported data that require further qualification. This code can be used ONLY under advisement of the LANL project chemist.	UJ, R, R19	J, R, R19
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	28. Qualification of data via data validation did occur, however no data quality control requirements in this procedure were applicable. Adhere to the external laboratory qualifiers found within the Form 1 analytical data summary sheets generated by the external laboratory. (AV)	U, U_LAB	J, J_LAB NQ, NQ