

LA-UR-18-29563

Approved for public release; distribution is unlimited.

Title: MC-15 Users Manual

Author(s): Moss, Calvin Elroy
Nelson, Mark Andrew
Rothrock, Richard Brian
Sorensen, Eric Byron

Intended for: Report

Issued: 2018-10-09

Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by the Los Alamos National Security, LLC for the National Nuclear Security Administration of the U.S. Department of Energy under contract DE-AC52-06NA25396. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

MC-15 Users Manual



1. Introduction

The current iteration of field-capable neutron multiplicity counters is the Multiplicity Counter 15 tube detector (MC-15) shown in Figure 1. This was a tri-lab effort between Los Alamos National Laboratory (LANL), Lawrence Livermore National Laboratory (LLNL), and Sandia National Laboratories (SNL). The MC-15 incorporates significant improvements over previous man-portable neutron multiplicity counters of the LANL NPOD and the Ortec Fission Meter™. These improvements include finer time resolution of 100 ns, more advanced on-board computing power for near-real-time onboard analysis, a sunlight readable color touchscreen for ease of use, a USB 2.0 interface for data transfer, and hot-swappable rechargeable batteries that provide 14 hours of run time.

The layout of the MC-15 has a total of fifteen ^3He tubes embedded in high density polyethylene (HDPE). Seven of the ^3He tubes are aligned in a removable front panel. The remaining eight tubes are in the main body in two rows with the second, or middle, row containing six ^3He tubes and the third, or last, row containing two ^3He tubes.

The MC-15 records when a neutron is detected and which tube detected the neutron. The time resolution is 100 ns and this raw data is stored in a file that can be downloaded later. This list-mode file for a five min run is on the order of 100 Mbytes. Additionally, the MC-15 can analyze the data internally and display the results on the touchscreen. Data and results are stored in the internal non-volatile memory and can be exported to a USB storage device. Data collection can be started by turning the unit on, touching [Collect Data](#), touching [Start](#) on the next screen, confirming or changing the parameters on the next screen, and touching [Start](#). While the data collection is in process, analysis occurs, and real time summary data is displayed. More capabilities are discussed in this document.



Figure 1. MC-15 with canvas cover.

1.1. Purpose

This purpose of this document is to aid users of the MC-15 to achieve a better understanding of the capabilities and limitations this instrument.

1.2. Scope

This document is written to provide information and guidance to users of the MC-15.

1.3. Precautions and Limitations

1.3.1. Electrical Safety

None of the procedures described in this manual require disassembly of the MC-15. However, inside each unit is a high-voltage power supply capable of 2 kV and 7 μ A of current. There is 0.002 μ F of capacitance, which if charged to 2 kV, would store 0.004 J. Los Alamos National Laboratory classifies this

in Electrical Safety Program P101-13, Rev. 1 as Class 2.1 R&D electronic work, which requires non-energized electrical safety training. The internal lithium-ion battery has a sophisticated integral protection circuit that protects the battery from over voltage when charging (≥ 4.3 V to any cell), over – current when discharging (≥ 8.25 A), and high temperatures when discharging ($>75^{\circ}\text{C}$).

1.3.2. Additional Precautions

The MC-15 weighs 47 pounds and care is required in lifting. Each unit contains 15 Reuter-Stokes ^3He proportional gas tubes 1 inch in diameter by 17 inches long and pressurized to 10 atm. These tubes are classified as a class/division hazard rating of 2.2 by the U.S. Department of Transportation.

2. Hardware Description

2.1. Electrical requirements

The MC-15 can be powered with a wide range of DC voltage. The standard charger that is supplied with each MC-15 gives an output of 24 V @ 3.25 A. The smallest voltage that will charge the MC-15 is 14 V at 5 amps. The maximum voltage allowed for the MC-15 is 24 V.

Minimum Voltage Power Supply	8 V @ 5 amps (no charging)
Maximum Voltage Power Supply	24 V @ 3.25 Amps

2.2. ^3He Tubes and Moderator

Figure 1 shows the fully assembled MC-15 in a canvas bag, and Figure 2 shows the layout of the ^3He tubes inside the assembled MC-15 high density polyethylene (HDPE) moderator. The dimensions of the MC-15 are 38×13×56 cm (15.0 x 5.1 x 22.0 inches), the weight is 21 kg (47 lbs) and the total volume of ^3He at standard temperature and pressure is ~28 liters (1,710 cubic inches). The canvas bag, which can be removed if desired, contains a 0.08 cm (0.030 inch) thick sheet of cadmium on the front side of the detector facing the source to minimize the detection of neutrons that have been thermalized by the environment.

Each ^3He tube is connected to a single preamp. This reduces the effects of dead time on the collected data. The preamps were designed at LANL and are the model KM-100. These preamps have been optimized for the ^3He tubes in the MC-15 in order to decrease dead time and increase the maximum counting rate.

For identification purposes the tubes are numbered 1 through 15 in the configuration outlined below in Figure 2. The tubes are embedded in a high density polyethylene moderator. The density of HDPE is Class 4, which specifies a density > 0.96 g/cc per ASTM D4976-12a.




Figure 2. Layout of the ^3He tubes in the MC-15.

The MC-15 has a removable front panel. This modularity allows for the MC-15 to be reconfigured in a variety of geometries.

The front panel is secured to the main body by four thumb screws as shown in Figure 3.



Figure 3. Locations of screws fastening the MC-15 front panel to the back panel.

When the front panel is removed from the main body, the front panel needs to be connected to the main body with the special cable as shown in Figure 4. The connectors are water-tight and labeled "Front Panel" and "Facing Out" for the back panel.




Figure 4. Special cable connecting the front panel (left) and the back panel (right).

During reassembly the front panel is aligned with the main body by a guide to insure the electrical connection is aligned as shown in Figure 5. Care must be taken to insure the electrical connection on the front body is on the same side as the electrical connection on the main body.




Figure 5. Front panel partially assembled with the main body (left in photos) and the back panel (bottom in photos)

Figure 6 shows a comparison of the MC-15 fully assembled and a configuration with the front panel separated.




Figure 6. MC-15 with front panel attached (left); front panel separated (right).

2.3. Interface

The main operation of the MC-15 is performed through an LCD screen on the top of the main body. All of the control electronics with firmware is contained in the detector head shown in Figure 7. The high resolution 3.5-inch color LCD touch screen has an extended temperature range and is readable in sunlight.




Figure 7. MC-15 detector head.

The electrical ports used for data transfer, interfacing with another detector or a computer, and charging are housed on the left side on top of the MC-15. These connections are protected from the environment by a cover with a gasket. These connections are not water proof, and it is recommended that the cover be closed when taking measurements in a wet environment. The interface connectors are shown in Figure 8.

The external power input (EXT PWR in Figure 7) is used to charge the batteries or to run the unit on AC power. The red LED is on when the battery is charging and goes off when the battery is fully charged. The green light is on whenever external power is connected. The top battery can be removed and charged separately with a charger (see Section 6). A stand-alone computer can be connected to the “LCD PGM” port to modify the display software for different applications. Data can be transferred to a USB device, such as a memory stick, via the “USB” port. The data can be read out in condensed mode (CMX) (See Appendix B), which creates Feynman histograms at a select number of gates, or in list mode

(LMX), which produces a larger file with time tagged data for each event (See Appendix A). The Ethernet port allows a stand-alone computer to be connected via a long cable for remote control and data download. This is essential for data collection in a high radiation environment. The “DETECTOR 2” allows a second MC-15 to be connected. For this configuration the data for all 30 channels, 15 for the primary detector and 15 from the secondary detector, are stored in the primary MC-15. The “IN” connector (pulse on $> \sim 2.0\text{v}$, off $< 2.0\text{v}$, max 6.5 v , minimum width $\sim 50\text{ns}$) allows an external detector to be counted in a 16th channel. Usually this detector is the SNAP, which provides an accurate measurement of the neutron source strength, but future measurement protocols may make this detector unnecessary. The “OUT” connector (pulse on $> 2\text{v} - 3.3\text{v}$, off $< 2\text{v}$, width 60ns) provides an output pulse whenever the MC-15 registers a count in any channel. Using the “IN” and “OUT” connectors allow several MC-15 instruments to be daisy chained together for greater detection efficiency. However, the individual channels are not recorded in list mode when more than two MC-15s are connected together in daisy chained mode. The “GATE” connector (pulse on $> \sim 2.3\text{v}$, off $< 1.5\text{v}$, max 6.5v , minimum width $\sim 50\text{ns}$) allows active interrogation techniques to veto data acquisition to minimize the size of data files.




Figure 8. MC-15 I/O interface connectors.

2.4. Batteries

Power is provided by two lithium ion smart batteries like the one shown in Figure 9. The locations of the two lithium ion batteries in the MC-15 are shown in Figures 10 and 11. The top battery can be hot swapped. A Phillips #1 screwdriver is need to access the bottom battery and is usually not swapped. The charging level on each battery is indicated by bars on the display on the end of each battery as well as on some of the MC-15 screens. A single battery provides seven hours of run time and two batteries provide 14 hours. The charging time is six hours.



Figure 9. Lithium ion smart battery.



Figure 10. The top battery can be hot swapped.




Figure 11. Tools are required to access the bottom battery.

3. Deployment

The front of the MC-15 should be placed toward the source. In order to measure accurate absolute neutron outputs, the MC-15 should be placed as accurately as possible so that the center mark on the front face of the MC-15 is centered vertically and horizontally on the neutron hot spot and 30 cm from the center of the object of interest. Sections 3.1-3.2 give preferred configurations, but if these are not possible, the user should record the dimensions for future analysis. While data is being collected, personnel movement between the MC-15 and the object must be avoided. The internal memory and external USB memory stick should be erased at the start of a deployment exercise. The minimum data acquisition time is 5 minutes

3.1. Passive Measurement Setups

3.1.1. Single MC-15

A simple passive configuration setup with one MC-15 is shown in Figure 12. The source output is assumed to be low enough so that the operator can be present at the MC-15. With the front of the MC-15 toward the source, the operator on the opposite side will be able to read the display. The source is centered in front of the MC-15 and at the same height from the floor as the + on the MC-15. A distance of 30 cm from the front face of the MC-15 to the center of the source is recommended. Typically, a 5 minute count is sufficient for adequate statistics. However, a longer time might be necessary if the count rate is very low.



Figure 12. Setting up the MC-15 for a standard measurement.

The efficiency of the detection system can be increased by removing the front panel of tubes and placing it separately near the source as shown in Figure 13. Earlier in this manual, Figure 3 shows how the front panel is attached to the back panel, Figure 4 shows the special cable connecting the front and rear panels, and Figure 5 shows how the panels must be aligned to reassembly them. Separating the panels increases the efficiency for low energy neutrons but reduces it for high energy neutrons. The front panel must be connected to the rest of the MC-15 with the special cable provided (Figure 4).



Figure 13. Two panel passive setup.

3.1.2. Two MC-15s

A simple passive configuration setup with two MC-15s is shown in Figure 14. The dimensions are the same as for the single MC-15. A distance of 30 cm from the front face of each MC-15 to the center of the source is recommended. The secondary MC-15 must be connected to the Detector 2 connector on the primary MC-15 using a special cable provided. One end of this cable is labeled “Primary” and the other is labeled “Secondary.” Connecting this cable will force the MC-15s to be primary and secondary to correspond to the cable labels. This configuration choice will be shown on the “Analysis Settings” screen (Figure 24 or 25) and the “Active Data Channels” screen (Figure 37).




Figure 14. Setting for a standard measurement with two MC-15s.

Figure 15 shows a measurement in progress with a passive setup using two MC-15s with the front panels removed.



Figure 15. Example of a two MC-15 setup with the front panels removed to enhance the efficiency for low energy neutrons.

3.1.3. One MC-15 with SNAP

A SNAP (Shielded Neutron Analysis Probe), which should provide a more accurate value of the neutron source output than the MC-15, can be added to the setup. The SNAP should be located directly opposite the MC-15 at 100 cm as shown in Figure 16. The SNAP center reference point is indicated by a + on the top of the SNAP. The recommend distance for the MC-15 is still 30 cm. The SNAP should be connected to the “IN” input of the MC-15 with a LEMO cable, and recording of this input is enabled by default but can be confirmed by checking that the Ext. Input box in the Active Data Channels screen (Figure 37) is checked.



Figure 16. Simple passive setup with the SNAP (left) and MC-15 (right).

3.1.4. Two MC-15s with SNAP

The SNAP can also be used with two MC-15s as shown in Figure 17. The recommended distances are 30 cm for the MC-15s and 100 cm for the SNAP. Obviously, the detectors must be arranged so that an MC-15 is not between the source and the SNAP. Also, the reflection of neutrons by the MC-15s to the SNAP should be minimized by placing the MC-15s to the side as shown in Figure 17 instead of behind the test object and opposite the SNAP.



Figure 17. Simple setup with two MC-15s and the SNAP

3.2. Active Measurement Setups

As with the passive measurement setup, many detector configurations are possible. The main difference is that a DT neutron generator is placed close to the object of interest. The neutron generator output point should be at the same height as the heights of the object being interrogated and the MC-15(s). The generator must be placed so as not to attenuate the delayed neutrons traveling toward the detector(s). For radiation safety, the control of the DT generator must be accomplished remotely. The MC-15 could also be controlled remotely, or the operator could start the MC-15 and retreat to a safe distance before starting the DT generator remotely. More details about active measurements, including details about the operation of the generator, are available in reference [1].

3.2.1. Active with One MC-15

Figure 18 shows a simple active configuration setup that includes a DT neutron generator, MC-15, and SNAP. The SNAP detector is connected to the IN port, which is recorded in channel 16. The MC-15 may be controlled by a remote laptop via a cable connected to the Ethernet connector. A signal from the DT generator is connected to the GATE input on the MC-15 to gate off the data collection during the interrogating neutron pulse in order to reduce the size of the data files. The Veto Active, Veto Duration, and Veto True State are selected in the Setup screen (Figure 32).



Figure 18. Simple active configuration setup.

3.2.2. Active with Two MC-15s


Figure 19 shows a simple active configuration setup that includes a DT neutron generator, two MC-15s, and a SNAP. In this setup, the MC-15s are at 50 cm but the SNAP is still at 100 cm.



Figure 19. Active setup with two MC-15s and the SNAP.

4. MC-15 Control

All of the MC-15 operations are controlled either on the instrument itself using the high-resolution color touch screen LCD or remotely by a laptop connected via an Ethernet cable. Figure 17 is a flow chart of the commands. The following sections discuss the various screens.






Figure 20. Flow chart of the MC-15 control screens on two pages.

Starting Operation

When the MC-15 is powered on, the system immediately begins booting up. First the screen displays the boot screen shown in Figure 21 with a progress bar indicating the status. The line below the status

bar reports the sequence of operations: reading configuration, building user interface, setting parameters, building file list, and cleaning up. This screen is only displayed during boot up. Possible warning pop-up messages are “Memory Low (< 5%)” and “Battery Low (<10%)”. When the boot up is finished, the Main screen shown in Figure 22 is displayed. The boot up typically takes 20 seconds.




Figure 21. Boot up screen.

Main Screen

The main screen (Figure 22) shows the date, time, and neutron count rate in the ^{15}He tubes in the MC-15. The secondary rate is the rate in a second MC-15 if a second MC-15 is connected to the primary MC-15. The time and rates are updated once per second. The bars in the upper right corner indicate the charging level of the upper and lower lithium batteries. The battery color is blue if normal, green if charging, and red if low. Three links are provided that navigate to their respective interface elements. This screen will power off if there is no activity in 15 minutes.




Figure 22. Main navigation screen.

Neutron Analysis Mode

If [Collect Data](#) is selected in the Main screen (Figure 22), the Neutron Analysis screen shown in Figure 23 is displayed. This screen gives the time completed and total run time, the name of the run when completed, the total counts, and the average rate. Ym is the reduced Feynman Variance with a bin width of 256 μ s. C-bar and C2-bar are the moments in the Feynman Variance analysis.

Row1/Row2 is the ratio of the count rates in four (2+3+5+6) of the front seven detectors to four (9+10+11+12) of the middle six detectors. Similarly, Row1/Row3 is the ratio of the count rates in the same four of the front seven detectors to the back two detectors. These ratios give some indication of the energy of the neutrons. When this screen is first accessed after the boot up, all of the values in the box will be zero. The count rate at the bottom is updated every second. The battery charging levels are displayed in the upper right corner. This screen provides links to other screens and functions: Start, Clear, Settings, and Save Run.




Figure 23. Neutron Analysis screen.

Analysis Settings with Start Option

If [Start](#) is selected in the Neutron Analysis screen (Figure 23), the Analysis Settings with Start Option screen shown in Figure 24 is displayed. Touching the Run Description panel opens a key pad so that a description can be entered. The pull down menu in the Analysis Settings panel allows the selection of either a single detector or two detectors. Checking the box labeled "Use Ext Input in Ym" enables the use of the channel 16 counts, received via the IN input, in the calculation of Ym. If there are two MC-15s, channel 32 counts are also used. Touching the Run Time: Seconds, Repetitions, Distance to Object [cm], or Distance to Floor [cm] box opens a numeric keypad so that the value can be changed. The maximum number of repetitions is limited by the number of files available in internal storage (200 max). A run can be started by touching [Start](#). The maximum data file size is 2GB and the maximum sample count is 260 million. When the sample count exceeds this value, a message will appear on the LCD

telling the operator that the assay has been cancelled, and the EndOFAssayStatus in the Imx header will show “Cancelled by file size limit.” (See Appendix A)

Back Confirm Settings Start

No Description Provided

Single Detector ☐ Use Ext Input in Ym

1000 <-- Run Time : Seconds

1 <-- Repetitions

0 <-- Object to Detector [cm]

0 <-- Floor to Detector [cm]


Figure 24. Analysis Settings with Start Option screen

Clear

If [Clear](#) is selected in the Neutron Analysis screen (Figure 23), the elapsed time is set to 0, the run label changes to blank, and all of the values in the box are set to 0. The current rate will continue to update every second.

Settings

If [Settings](#) is selected in the Neutron Analysis screen (Figure 23), the Analysis Settings Only screen shown in Figure 25 is displayed. Touching the Run Description panel opens a key pad so that a description can be entered. The pull down menu in the next panel allows the selection of either a single detector or two MC-15 detectors. Checking the box labeled “Use Ext Input in Ym” enables the use of the channel 16 counts, received via the IN input, in the calculation of Ym. If there are two MC-15s, channel 32 counts are also used. Touching the Run Time : Seconds, Repetitions, Distance to Object [cm], or Distance to Floor [cm] box opens a numeric keypad so that the value can be changed. The only link from this screen is [Back](#) to the Neutron Analysis screen (Figure 23).




The 'Analysis Settings' screen features a 'Back' button at the top left. Below it is a text field containing 'No Description Provided' with scroll arrows on the right. A dropdown menu is set to 'Single Detector', followed by an unchecked checkbox labeled 'Use Ext Input in Ym'. Below these are four input fields with labels: '1000' for '<-- Run Time : Seconds', '1' for '<-- Repetitions', '0' for '<-- Object to Detector [cm]', and '0' for '<-- Floor to Detector [cm]'.

Figure 25. Analysis Settings Only screen.

Save Run

If **Save Run** is selected in the Neutron Analysis screen (Figure 23), the Save File Confirmation screen shown in Figure 26 is displayed. The Save Run button is only enabled once an analysis has been collected. Select the file types to proceed, either LMX, CMX or both. Touching **Save** will enable saving to the USB port if a USB storage device is connected to the USB port. The operation can be canceled by touching **Cancel** and the system returns to the Neutron Analysis screen (Figure 23).



The 'File Utilities' screen shows a 'Back' button, '92% Mem', and '74 Files'. A blue banner at the top states 'Estimated transfer time: 0.1 minutes!'. The main heading is 'Select file type(s) to proceed', with checkboxes for 'LMX File' (checked) and 'CMX File' (unchecked). 'Cancel' and 'Save' buttons are at the bottom. A list at the bottom shows file names and sizes: '1970_01_01_000949' (4532136) and '1970_01_01_000935' (4532400).

Figure 26. Save File Confirmation screen.

Cumulative Tube Events

If **Display Counts** is selected in the Main screen (Figure 22), the Cumulative Tube Events screen shown in Figure 27 is displayed. The 15 bars give the total counts per tube or the current count rate. The seven dark blue bars are for the tubes closest to the front surface, the light six light blue bars are for tubes in the next row, and the two aqua bars are for the tubes in the rear. The total count rates in tubes 1-15 and the count rate in channel 16 are displayed at the bottom. Initially the total counts per tube are displayed, but this can be switched to the 1 second count rates as shown in Figure 27 by touching **Mode**.

The displayed values can be reset to zero by touching [Restart](#). If [Toggle](#) is touched, the Count History screen as shown in Figure 28 is displayed.




Figure 27. Cumulative Tube Events screen.

Count History

The Count History screen shown in Figure 28 is selected by touching [Toggle](#) in the Cumulative Tube Events screen (Figure 27). The sum of the count rates in all tubes is displayed. If there is only one MC-15, then only one plot for channels 1-16 is shown. If there is a second MC-15 connected, then a second plot for channels 17-32 is shown. However, if the Combine Detector Results box is checked in the Count History Options screen (Figure 29), then only the sum of 1-16 and 17-32 is plotted. The rate at the bottom in Figure 28 is the sum of channels 1-16 if there is only one MC-15 connected or the sum of channels 1-32 if there are two MC-15s connected. The vertical scale will automatically scale if the Auto Scale box is checked or the scale can be manually changed with the [Up](#) or [Dn](#) buttons. [Pause](#) stops the updating of the display. If [Back](#) is touched, the Main screen is again displayed (Figure 22). If [Toggle](#) is touched, the Cumulative Tube Events screen is again displayed (Figure 27). If [Settings](#) is touched, the Count History Options as shown in Figure 29 is displayed.





Figure 28. Count History screen.

Count History Options

The Count History Options screen shown in Figure 29 is selected by touching [Settings](#) in the Count History screen (Figure 28). Touching the Dwell Interval (s), Time Range (s), or Alarm Threshold (n/s) will bring up a numeric key pad for changing the values. If the Combine Detector Results box is checked, the sum of 1-16 and 17-32 is plotted in the Count History screen (Figure 28) rather than separate plots for each MC-15. Touching [Clear Plot](#) will clear the Count History. If [Back](#) or [OK](#) is touched, the Count History screen (Figure 28) is again displayed.



Back Count History Options

1 <--Dwell Interval (s)

300 <--Time Range (s)

200 <--Alarm Threshold (n/s)

☐ Audible Alarm


☐ Combine Detector Results

Clear Plot OK

Figure 29. Count History Options screen.

Options

If [Options](#) is selected in the Main screen (Figure 22), then the Options screen is displayed as shown in Figure 30. Four links are provided to navigate to their respective screens.



Back Options

File Utilities

Setup

Diagnostics

Advanced Tools

Figure 30. Options screen.

File Utilities

If [File Utilities](#) is selected in the Options screen (Figure 30), then the File Utilities screen shown in Figure 31 is displayed. All of the data files stored in the internal memory are listed. To save a single file to an external USB memory stick, highlight it, then touch [Save](#). The Save File Confirmation screen shown in Figure 26 will be displayed. The internal memory capacity is 128 GB and the maximum number of files is 200. To save all of the files, touch [Save All](#). The screen shown in Figure 26 will be displayed. To erase all the files in the local files system, touch [Erase All](#). Touching [Back](#) returns to the Options screen (Figure 30)




Figure 31. Files Utilities screen.

Setup

If [Setup](#) is selected in the Options screen (Figure 30), then the Setup screen shown in Figure 32 is displayed. The top pull down menu allows the system to be configured for one or two MC-15s. The HV Reading is updated every second. The next three boxes allow the HV, Channel Dead Time, and Veto Duration to be entered. The bottom pull down menu allows the choice of “low” or “high” for the veto gate input. Check the Veto Active box to enable the veto gate. Audible Feedback can be enabled by checking the corresponding box. Audible Feedback controls the beeping at the end of an analysis. To set the date and time touch the [Set Time](#) to display the Instrument Date/Time screen shown in Figure 33. Touching [Back](#) returns to the Options screen (Figure 30).

Back Single Detector

HV Reading = 1692

1680 <-- HV set : 0-2000V

1000 <-- Channel Dead Time (nS)
(In 100 nS increments)

20 <-- Veto Duration (uS)
(In 100 nS increments)

High <-- Veto True State

☐ Veto Active

☒ Audible Feedback

Set Time

Figure 32. Setup screen.

Instrument Date/Time

The Instrument Date/Time screen in Figure 33 is accessed by touching [Set Time](#) in the Setup screen (Figure 32). The four boxes allow the Zulu date and time to be entered. Touching [Back](#) returns to the Options screen (Figure 30).

Back Instrument Date/Time

2018

August

8 / 9 / 2018

09:30:19 Zulu Time

Figure 33. Instrument Date/Time screen.

System Information

If [Diagnostics](#) is selected in the Options screen (Figure 30), then the System Information screen shown in Figure 34 is displayed. Information about the hardware and software/firmware versions is provided. The temperatures and HV Reading are updated every second. Touching [Power](#) opens the Power Information screen shown in Figure 35. Touching [Back](#) returns to the Options screen (Figure 30).




Figure 34. System Information for diagnostics.

Power Information

The Power Information screen (Figure 35) is accessed by touching [Power](#) in the System Information screen (Figure 34). The status of the two batteries is shown. Touching [Back](#) returns to the System Information screen (Figure 34).




Figure 35. Power Information screen.

Advanced Tools

If [Advanced Tools](#) is selected in the Options screen (Figure 30), then the Advanced screen shown in Figure 36 is displayed. Three links are provided to navigate to their respective screens.




Figure 36. Advanced screen.

Advanced Channels

If [Active Channels](#) is selected in the Advanced screen (Figure 36), then the Active Data Channels screen shown in Figure 37 is displayed. The top pull down menu allows the selection of the system configuration: a single detector or two detectors. The second pull down menu allows the user to control the active channels recorded on either the primary or secondary detector. This menu is only enabled if the top menu selection is "2 Detector Primary". If the Ext. Input box is checked when "Primary Detector" is displayed in the second pull-down menu, then counts received at the IN input of the primary MC-15 are recorded in channel 16. Similarly, if the Ext. Input box is checked when "Multiple Detectors" is displayed in the second pull-down menu, then counts received at the IN input of the secondary MC-15 are recorded in channel 32. Whether these counts are used in the Ym calculation is determined by the "Use Ext Input in Ym" box on the Analysis Setup Screens (Figures 24 and 25). Boxes 1 -15 can be used to select which tubes are recorded. Similarly, if "Multiple Detectors" is displayed, boxes 17-31 can be used to select which of these tubes are recorded. Touching [Back](#) returns to the Advanced screen (Figure 36).

Back

Active Data Channels

Single Detector

Ext. Input

Primary Detector

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

Figure 37. Active Data Channels screen.

Current Tube Rates

If [Chan Count Graphic](#) is selected in the Advanced screen (Figure 36), then the Current Tube Rates screen shown in Figure 38 is displayed. The arrangement of the circles presents the location of the tubes in the MC-15. The top seven are for the row of tubes closest to the front surface, the next six are for tubes in the next row, and the bottom two are for the tubes in the rear row. The current tube rate in each tube are color coded with the scale at the bottom. The maximum of the color scale is set to 20000 n/s. Touching [Scaling](#) sets the color maximum to the maximum tube rate as shown in Figure 39. The ratios of the rates in the rows (lower left) and the total rate in all of the tubes (bottom) are shown. All of the values are updated every second. Touching [Reset](#) zeros all of the values. Touching [Mode](#) changes to the Cumulative Tube Events display shown in Figures 40 and 41. Touching [Chans](#) switches to a display of channels 17-31, which are the values for the 15 tubes in a second detector, if present. Touching [Back](#) returns to the Advanced screen (Figure 36).




Figure 38. Current Tube Rates screen.




Figure 39. Color-scaled Current Tube Rates screen.




Figure 40. Cumulative Tube Events.




Figure 41. Color-scaled Cumulative Tube Events.

Miscellaneous Command

If [Debug Items](#) is selected in the Advanced screen (Figure 36), the Miscellaneous Command screen shown in Figure 42 is displayed. If “Count Rates,” “Feynman Data,” “Battery Info,” or “Temp Info” are unchecked, the information is not sent to the user interface. “Audible Rates” controls beeping with a frequency corresponding to the count rate. Touching [Com In](#) brings up the Standard Input (COM1) screen shown in Figure 43. Touching [Com Out](#) brings up the Standard Output (COM1) screen shown in Figure 44. These COM screens show the dialog being transmitted between the user Ethernet interface and the embedded computer within the MC-15. They are generally only used for debug purposes. Touching [Back](#) returns to the Advanced screen (Figure 36).

Back	Com In
<input checked="" type="checkbox"/> Count Rates	Com Out
<input checked="" type="checkbox"/> Feynman Data	Clear Cmds
<input checked="" type="checkbox"/> Battery Info	C code Exit
<input checked="" type="checkbox"/> Temp Info	GUI Exit
<input type="checkbox"/> Audible Rates	
<input type="checkbox"/> Verbose C Code	
Cmd Queue Length = 1	
Response Queue Length = 2	

Figure 42. Miscellaneous Command screen

Standard Input

If [Com In](#) is selected in the Miscellaneous Command screen (Figure 42), the Standard Input (Com1) screen shown in Figure 43 is displayed. When turned on by touching [On](#), it shows the commands being received by the embedded system via the user Ethernet interface. Touching [Clear](#), clears the display box. Touching [Back](#) returns to the Miscellaneous Command screen.




Figure 43. Standard Input (COM1) screen.

Standard Output

If **Com Out** is selected in the Miscellaneous Command screen (Figure 42), the Standard Output (Com1) screen shown in Figure 44 is displayed. When turned on by touching **On**, it shows the commands being sent by the embedded system via the user Ethernet interface. Touching **Clear**, clears the display box. Touching **Back** returns to the Miscellaneous Command screen (Figure 42).




Figure 44. Standard Output (COM1) screen.

5. Shipping Case

The MC-15 and all essential accessories can be packed in one shipping case as shown in Figures 45 and 46.



Figure 45. MC-15 and shipping case.




Figure 46. Essential accessory items in the MC-15 shipping case.

6. Lithium Ion Smart Battery

The NI2040ED Smart Battery shown in Figure 47 includes a rechargeable Lithium Ion battery and a Battery Management Module. The battery consists of (9) Lithium Ion rechargeable cells of 18650 size, assembled in a 3 series / 3 parallel (3S 3P) configuration. Each cell has an average voltage of 3.6V and a typical capacity of 2.9Ah giving a battery pack of 10.8V and 8.7Ah typical. The battery is capable of communicating with the host or the charger through the System Management Bus (SMBus). The battery is fully SMBus and SBDS Revision 1.1 compliant. Protection is provided for over-charge, over discharge and short circuit. For redundancy, passive safety devices have been integrated into the pack to protect against over-current and over-temperature, and secondary over-voltage has been implemented with a logic-fuse and controller.



Figure 47. Lithium ion smart battery.

Table 1. NI2040 Smart Li Ion Battery Specification Summary

Part Number:	Ni2040ED29	2.9Ah cells, Extreme discharge
Chemistry / Cell Array	Nine Lithium Ion 18650 cells in a 3S3P array	
Capacity / Energy	8.7Ah / 94Wh	
Voltage	$V_{max}=12.6V$ $V_{nom}=10.8V$ $V_{cutoff}=7.2V$	
Equivalent Lithium Content	7.83g	
Maximum continuous discharge current	6A	
Maximum continuous discharge power to cutoff voltage	57W	
Weight:	496g, 1lb 2oz	
Communications:	Fully SMBus rev 1.1 Compliant	
Data System:	Fully SBDS rev 1.1 compliant	
Charging:	Inspired Energy "CH" range desktop charger or other SMBus Level 3 compliant charger (Linear Tech 4100 recommended).	
Height:	22.3 mm, 0.9"	
Length:	214.6 mm, 8.5"	
Width:	58.9 mm, 2.4"	
Recommended battery cavity:	23.5mm x 60mm cross section. The battery should be centered within this cavity.	
Engineering Data Sheet	Battery NI2040ED29 spec v1.1.pdf	

7. Specifications for the Reuter-Stokes ³He Proportional Gas Tubes

Model Number: SA-P4-0815-105

Mechanical

Maximum diameter: 1.03 in (26.2 mm)
Maximum overall length: 17.03 in (432.6 mm)
Maximum weight: 200 g

Material

Body material: 1100 aluminum
Internal insulator: alumina ceramic
Primary gas: ³He
Total pressure: 150 psia (1.03 MPa)

Electrical

Resistance: 1.00×10^{12} ohms
Capacitance: 8 pf $\pm 20\%$

Maximum Ratings

Maximum Voltage: 2500 volts
Maximum Temperature: 100°C

Operating Characteristics

Thermal neutron flux range: 9.2×10^{-4} to 7.72×10^2 nv
Thermal neutron sensitivity: 66 cps/nv $\pm 10\%$
Voltage Range: 1200 to 1700 volts
Resolution (FWHM): < 12 %

Temperature range: -25 to 100 °C

8. Reference

[1] C. E. Moss, S. M. Gonzales, J. D. Hutchinson, G. E. McKenzie IV, W. L. Myers, R. B. Rothrock, S. A. Salazar, M. A. Smith-Nelson, E. R. Sorensen, G. M. Sundby, "Field Guide for Active Interrogation (U)," LANL report LA-UR-17-22226, issued 2017-03-16.

Appendix A

*.lmx File format

revised 3-12-13

The List Mode eXtension (LMX) file format consists of two parts; a header section that contains information about the data collection and a binary section that contains the list mode data in binary format.

Header

The header section is in text (ASCII) for easy readability. Multiple spaces may exist at the end of each header line to pad each line to a specific number of characters. Each header line must end with carriage return and line feed (ASCII Codes = 13/10. Hex 0D 0A, format \r\n)

Each line of the header has three parts; a name, the associated value, and units if applicable and takes the form:

keyword : *value* [*units*]

If values have units, they must be present within brackets. If units are not applicable this term is not present.

Not all header lines need to be present and they do not need to be in any specific order (except for “ListModeExtensionFileVersion” must be first line and BinaryDataFollows must be the last line as described below.) Mandatory lines are indicated in Table 1. The file may include unique header lines as long as they do not duplicate the keywords in Table 1.

The header names are case sensitive in order to maintain consistency with the *.cmx file (the condensed file format).

Table 1. List of Header Keywords

Name	Value type	[Units]	Example Values	Man dator y	Comments
ListModeDataFileVersion	string		1.01	Y	This tag must be first in the file
InstrumentType	string		Neutron Multiplicity		
InstrumentModel	string		MC15	Y	
SerialNumber	string		S/N 026		
HardwareVersion	string		006		
FirmwareVersion	string		aa4f717		
CCodeVersion	string		90073114		
LCDVersion	string		1.1.0.184		
MeasurementID	string		2012_05_30_14 3526		
MeasurementDescription	string		Run 10A		Entered on LCD
MeasurementSample	string		1 of 30		Run number is a sequence
MeasurementMode	string		Multiple		MC-15 configuration: single or multiple
FrontPanelConfig	string		Together		Front panel configuration: together, separate, missing
AnalysisChannels	hex		0x00007FFF		Channels used in real time analysis and cmx file creation
DateStart	string	YYYY-MM-DD[Z]	2012-05-30 [Z]	Y	Date of start of assay
TimeStart	string	Hrs:min:sec	14:35:26 [Z]	Y	Time the assay was started
DurationRealTime	positive real	s	600 [s]		Duration of assay
EndOfAssayStatus	string		Normal		Normal, Cancelled by user, Cancelled by file size limit
InternalScaler	real		5000		Counts recorded in firmware during assay

FiFoLostCounts	real		0		Number of counts not read from firmware
AverageCountRate	real	Counts/s	8.33 [counts/s]		Calculated after assay. Internal Scaler/duration
DistanceDetFaceToSource	real	cm	23.5 [cm]		“unknown” if not entered
DistanceDetCenterToFloor	real	cm	0.0 [cm]		“unknown” if not entered
HighVoltageSetPoint	real	V	1600 [V]		
HighVoltageActual	real	V	1598 [V]		
TemperatureCPU	string	C	78 [C]		Temperature of FPGA surface
TemperatureLCD	string	C	78 [C]		Temperature of LCD unit
TemperatureInternal	string	C	78 [C]		Temperature of instrument air
HumidityInternal	real	% RH	38 [%]		% relative Humidity inside instrument
BarometerElevation	real	meters	2146 [m]		Not compensated
FirmwareChannelDeadtime	Real	Nanoseconds	1000 [ns]		Time interval that a second event on a channel cannot be recorded
RowRatio(1/2)	real		1.01		Front row/middle row
RowRatio(1/3)	real		0.988		Front row/back row
RowRatio(2/3)	real		0.938		Middle row/back row
Comment					
BinaryDataEventSizeInBytes	Integer		8		Number of bytes that make up every event
BinaryDataChannelFormat	string		unsigned int32		
BinaryDataClockFormat	string		unsigned int32		
BinaryDataClockTickLength	ns		100 [ns]		Time interval each clock tick represents
BinaryDataActiveChannels	hex		0x00007FFF		Channels input to firmware
BinaryDataFollows				Y	Binary data start immediately after a line feed/carriage return (ASCII codes = 13/10)

Here are the 3 possible “EndOfAssayStatus” states that are written into the lmx header. Only 1 will appear.

EndOfAssayStatus	: Cancelled by file size limit
EndOfAssayStatus	: Cancelled by user
EndOfAssayStatus	: Normal

Cancelled by file size limit happens automatically in software.

Cancelled by user happens when, you guessed it, the user pushes the cancel button.

Normal happens when time expires according to “DurationRealTime” being reached.

If the assay is cancelled, “DurationRealTime” is adjusted to the actual end time. It may have started out at 1000s or something, but is written as the actual end time.

Here is an example of the latest version of the lmx header (when file size limit was reached). (Required lines are in Red)

ListModeDataFileVersion	: 2.01
InstrumentType	: Neutron Multiplicity
InstrumentModel	: MC-15
SerialNumber	: S/N 007
HardwareVersion	: 6a
FirmwareVersion	: 1bdf0417
CCodeVersion	: 10170417
LCDVersion	: 1.0.0.182
MeasurementID	: 2017_04_27_110748
MeasurementDescription	: ULTRA HIGH C-252 TEST
MeasurementSample	: 1 of 100
MeasurementMode	: Multiple
FrontPanelConfig	: Together
AnalysisChannels	: 0x7FFF7FFF
DateStart	: 2017-04-27 [Z]
TimeStart	: 11:07:48 [Z]
DurationRealTime	: 174 [s]
EndOfAssayStatus	: Cancelled by file size limit
InternalScaler	: 271198097

FifoLostCounts	: 0
AverageCountRate	: 1558609.88 [counts/s]
DistanceDetFaceToSource	: 0.0 [cm]
DistanceDetCenterToFloor	: 0.0 [cm]
HighVoltageSetPoint	: 1680 [V]
HighVoltageActual	: 1682 [V]
TemperatureCPU	: 39 [C]
TemperatureLCD	: 0 [C]
TemperatureInternal	: 33 [C]
HumidityInternal	: 16 [%]
BarometerElevation	: 2289 [m]
FirmwareChannelDeadtime	: 1000 [ns]
RowRatio(1/2)	: 0.909
RowRatio(1/3)	: 1.353
RowRatio(2/3)	: 1.489
Comment	: If all channels are zero, this indicates channels
Comment	: in next event is a flag (not a real event)
Comment	: Flag = 0x00000001 (Clock rollover occurred)
Comment	: Flag = 0x00000002 (Gate input started)
Comment	: Flag = 0x00000003 (Gate input ended)
Comment	: Flag = 0xFFFFFFFF (End of binary data)
BinaryDataEventSizeInBytes	: 8
BinaryDataChannelFormat	: unsigned int32
BinaryDataClockFormat	: unsigned int32
BinaryDataClockTickLength	: 100 [ns]
BinaryDataActiveChannels	: 0xFFFFFFFF
BinaryDataFollows	:

Binary data

The binary section of data is started with the tag

BinaryDataFollows : n number of spaces \r\n Binary data start immediately after a carriage return/line feed (ASCII codes = 13/10)

Format of Binary Data: [Channel Data] [Time]

Channel Data and Time are written in little endian, i.e., least significant byte first.

Insertion of specific flags/codes into Binary Data:

A flag is captured in 2 event records. The event size in bytes is defined in the header section (BinaryDataEventSizeInBytes :). All zeros in every channel indicate that the next event is a flag (not a neutron event). The all-zero channel data is paired with the time tag associated with the flag occurrence. The next event describes the flag type and is paired with the same time tag.

Examples: Flag = 1: Clock rollover occurred
 Flag = 2: Hardware veto/gate started
 Flag = 3: Hardware veto/gate ended
 Flag = FFFFFFFF : End of File indicator

Clock Rollover: 00000000FFFFFFFF (all zeros in channels along with maximum clock tic before rollover)
 00000001FFFFFFFF (flag type 1 in channels along with maximum clock tic before rollover)

End of File: 000000001B037AD1 (all zeros in channels along with clock tic at end of assay)
 FFFFFFFF1B037AD1 (flag type FFFFFFFF in channels along with clock tic at end of assay)

Appendix B

Example of CMX Data File

```
<CondensedNeutronDataFile version="1.07">
  <InstrumentInformation>
    <InstrumentType>Neutron Multiplicity</InstrumentType>
    <InstrumentModel>MC-15</InstrumentModel>
  </InstrumentInformation>
  <MeasurementData type="multiplicity" sample="1 of 1">
    <MeasurementID>2015_09_28_080322</MeasurementID>
    <MeasurementDescription>No Description Provided</MeasurementDescription>
    <StartTime>2015-09-28T08:03:22 Z</StartTime>
    <RealTime units="s">300</RealTime>
    <FIFOLostCounts>0</FIFOLostCounts>
    <DetectorInformation>
      <SerialNumber>S/N 007</SerialNumber>
      <HardwareVersion>005</HardwareVersion>
      <FirmwareVersion>aa9f0703</FirmwareVersion>
      <CCodeVersion>29000723</CCodeVersion>
      <LCDVersion>167e</LCDVersion>
      <MeasurementMode>Multiple</MeasurementMode>
      <FrontPanelConfig>Together</FrontPanelConfig>
      <ActiveChannels>0xFFFFFFFF</ActiveChannels>
      <AnalysisChannels>0x7FFF7FFF</AnalysisChannels>
      <DistanceDetectorFaceToSourceCenter units="cm">30.0</DistanceDetectorFaceToSourceCenter>
      <DistanceDetectorCenterToFloor units="cm">22.0</DistanceDetectorCenterToFloor>
      <TemperatureCPU units="C">61</TemperatureCPU>
      <TemperatureLCD units="C">36</TemperatureLCD>
      <TemperatureInternal units="C">0</TemperatureInternal>
      <HumidityInternal units="percent">42</HumidityInternal>
      <HighVoltageSetPoint units="V">1680</HighVoltageSetPoint>
      <HighVoltageActual units="V">1680</HighVoltageActual>
      <BarometerElevation units="m">0</BarometerElevation>
    </DetectorInformation>
  </MeasurementData>
</CondensedNeutronDataFile>
```

```

    <FirmwareChannelDeadtime units="ns">100</FirmwareChannelDeadtime>
  </DetectorInformation>
  <ProcessedNeutronData>
    <CountHistory channelmask="0xFFFFFFFF">
      <Interval units="s">1</Interval>
      <CountData>9006 9148 9084 9183 9258 9299 9124 9087 9097 9179 9230 9159 9043 9221 9027 9094 9121 9022 9082 8971 8776 8797
8946 8827 9352 9079 9259 9175 9111 9003 9050 9209 9101 9312 9149 9179 8972 9179 9080 9357 9125 8881 9136 9079 9146 9038 9144 9147
9089 8959 8884 9217 9010 9087 9219 9056 9185 9154 9059 9183 8998 9116 9314 9154 9218 9272 9143 9067 9201 9076 9181 9149 9063 9227
9051 9088 9229 9168 9184 9238 9246 9099 9066 9318 8969 8962 9195 9059 9170 9020 9190 9141 9174 8986 9079 9075 9210 8952 9039 9173
8994 9161 9232 8991 9155 9163 9188 9046 9280 8886 8916 9177 9104 9118 9219 9070 9115 8995 9064 9241 9124 9073 9084 9015 9013 9177
9031 9104 8792 8642 8878 8623 8785 8764 8914 9137 9068 9202 9276 9106 8963 9224 9016 8933 9116 9034 9070 9226 8942 9084 9193 9155
9095 9336 9024 9226 9167 9085 8925 9221 9101 9109 9062 9081 9038 9040 9060 9235 9117 9196 9161 9064 9261 9020 9188 9145 8917 9342
9214 9142 8886 9017 9192 9181 9331 9064 9180 9024 9098 9166 9049 9147 9075 9033 9021 9147 9111 9152 9056 9067 9254 8859 9283 9173
9209 9172 9071 9096 9250 9105 9151 9115 9104 8941 8890 8946 8979 9048 9115 9122 9249 9113 9061 9248 9279 9035 9072 9197 9141 9083
9152 9118 8998 9306 8951 9226 8949 9014 9116 9158 9033 9167 9064 8974 8990 9266 9127 9278 9060 9263 9035 9105 9198 9163 9181 9178
8860 9058 9033 8961 8967 9276 9077 9080 9165 9151 9072 9183 9194 9146 9031 9074 9083 8885 8877 9181 9079 9146 9022 9357 9274 8876
8997 9005 9237 8998 9032 9092 9278 9097 9131 8877 9084 9191 9104 9068 9092 9171 9152 9124 </CountData>
    </CountHistory>
    <RowRatio>
      <ChannelsUsedInRowOne>2 3 5 6</ChannelsUsedInRowOne>
      <ChannelsUsedInRowTwo>9 10 11 12</ChannelsUsedInRowTwo>
      <ChannelsUsedInRowThree>14 15</ChannelsUsedInRowThree>
      <RowOneOverTwoRatio>0.678</RowOneOverTwoRatio>
      <RowOneOverThreeRatio>1.053</RowOneOverThreeRatio>
      <RowTwoOverThreeRatio>1.554</RowTwoOverThreeRatio>
    </RowRatio>
    <GrossCounts channelmask="0xFFFFFFFF">
      <TotalCounts>2730123</TotalCounts>
      <RealTime units="s">300</RealTime>
      <AverageCountRate units="counts/s">9100.41</AverageCountRate>
      <ChannelData>66916 77492 83384 85495 83299 77071 68065 105625 114600 122520 123537 113442 105164 76002 76539 0 66416
75899 81004 83725 81209 75845 66335 104354 112643 120796 120043 108939 103705 74976 75083 0</ChannelData>
      <ExternalInput channel="16">0</ExternalInput>
      <ExternalInput channel="32">0</ExternalInput>

```

```

</GrossCounts>
<FeynmanHistograms gateunits="ns" method="sequential" channelmask="0x7FFF7FFF">
  <FeynmanData gatewidth=" 20096">12475571 2197486 234740 19114 1354 75 4 0</FeynmanData>
  <FeynmanData gatewidth=" 120448">880883 867296 473688 187340 60295 16240 3894 872 162 25 4 2 0</FeynmanData>
  <FeynmanData gatewidth=" 185216">332234 493730 398178 229854 105648 40469 13873 4161 1190 303 69 18 2 1
0</FeynmanData>
  <FeynmanData gatewidth=" 249472">143463 285150 302264 226859 133502 66265 28579 11003 3737 1217 355 108 30 6 1
0</FeynmanData>
  <FeynmanData gatewidth=" 317056">63903 160558 212925 197445 143593 87074 45336 21176 8940 3477 1221 370 145 32 9 1
0</FeynmanData>
  <FeynmanData gatewidth=" 392576">27498 84538 137106 155099 135346 98652 61812 33973 17106 7781 3213 1337 472 160 66 18
6 0</FeynmanData>
  <FeynmanData gatewidth=" 478592">11049 41043 79730 108019 113491 98147 72252 47499 27864 14854 7322 3227 1441 539 233
82 38 3 2 3 0</FeynmanData>
  <FeynmanData gatewidth=" 578944">3818 17739 41049 66514 82675 84384 73722 56787 39030 24392 14048 7346 3688 1750 736
304 125 50 17 5 3 0 1 0 0 1 0</FeynmanData>
  <FeynmanData gatewidth=" 711296">1042 5732 16350 32260 48442 59723 62180 57377 47017 34808 23600 14840 8732 4866 2581
1201 585 250 95 49 26 5 2 1 0 0 1 0</FeynmanData>
  <FeynmanData gatewidth=" 856448">255 1724 5749 13692 24314 35308 43343 46940 45161 39436 31853 23078 15848 10141 6182
3450 1904 1004 492 229 93 53 19 12 2 0 1 0</FeynmanData>
  <FeynmanData gatewidth=" 1030784">43 394 1640 4550 9439 16345 23965 30486 34792 35369 33523 28518 22842 17225 12114
8027 5036 3025 1784 958 498 241 114 65 29 10 3 0 3 1 0 0 1 0</FeynmanData>
  <FeynmanData gatewidth=" 1240192">2 66 332 1074 2796 5641 9779 15196 20240 24265 26492 26534 24850 22002 18426 14129
10365 7318 4878 3087 1910 1194 602 358 172 101 47 21 13 2 1 2 1 1 0 0 1 0</FeynmanData>
  <FeynmanData gatewidth=" 1492608">2 10 56 190 595 1383 2824 5079 8138 11665 15171 18119 19738 20239 19662 17833 15538
12462 9857 7377 5277 3620 2356 1559 971 555 349 170 97 48 25 11 7 3 3 1 0</FeynmanData>
  <FeynmanData gatewidth=" 1796224">0 1 4 25 80 217 505 1196 2147 3692 5648 8061 10470 12435 14162 15148 15352 14738 13527
11848 9752 7790 6019 4489 3325 2249 1549 1033 649 380 222 133 78 48 17 15 9 3 0 0 1 0</FeynmanData>
  <FeynmanData gatewidth=" 2163584">0 0 0 0 9 23 54 182 322 700 1251 2068 3220 4621 6292 7852 9472 10510 11376 11696 11365
10794 9719 8439 7019 5819 4485 3492 2488 1806 1227 840 596 355 225 137 88 47 30 20 11 5 3 0</FeynmanData>
  <FeynmanData gatewidth=" 2605184">0 0 0 0 0 1 4 13 31 62 168 294 550 901 1503 2210 3179 4171 5446 6260 7502 8125 8575 8742
8638 8174 7513 6727 5748 4977 3854 3135 2491 1852 1275 974 695 468 364 225 131 75 34 34 15 10 6 1 1 0 1 0</FeynmanData>
</FeynmanHistograms>
<TimeIntervalAnalysis channelmask="0x7FFF7FFF">

```

<Interval units="ns">128</Interval>

<NumberChannels>1024</NumberChannels>

<ChannelData>1888 3966 3873 3942 3850 3833 3843 3867 3837 4056 3999 3951 4006 3874 3870 3910 3821 3940 3917 3903 3838

3927 3956 3846 3671 3843 3755 3651 3773 3708 3788 3719 3624 3753 3705 3693 3696 3616 3692 3595 3586 3566 3573 3572 3630 3620 3502
3531 3458 3627 3567 3518 3531 3475 3501 3597 3419 3351 3490 3419 3308 3365 3438 3447 3415 3527 3416 3419 3407 3329 3359 3367 3251
3350 3331 3297 3264 3332 3324 3266 3247 3286 3288 3260 3270 3257 3188 3220 3200 3176 3179 3161 3303 3242 3265 3229 3162 3104 3097
3112 3201 3107 3178 3137 3040 3107 3144 3116 3115 3091 3175 3154 3060 3109 3117 2999 2993 3044 3098 2984 2952 3068 3038 3008 2909
2954 2998 2835 2771 3017 2861 2978 2943 2961 2829 2912 2892 2964 2926 2921 2834 2821 2963 2816 2802 2896 2860 2865 2958 2838 2839
2836 2918 2768 2866 2779 2798 2779 2702 2855 2827 2821 2733 2729 2833 2754 2760 2760 2572 2685 2752 2709 2715 2775 2780 2697 2623
2711 2746 2610 2684 2710 2669 2663 2743 2755 2675 2719 2658 2566 2605 2687 2676 2566 2622 2635 2574 2754 2552 2504 2629 2666 2650
2626 2616 2561 2559 2562 2646 2554 2470 2501 2533 2458 2600 2598 2507 2537 2522 2447 2551 2528 2595 2484 2547 2476 2540 2462 2478
2453 2450 2448 2370 2454 2536 2413 2408 2461 2435 2396 2381 2456 2446 2318 2359 2470 2318 2401 2416 2385 2414 2373 2336 2449 2408
2354 2410 2293 2354 2400 2372 2317 2352 2399 2261 2326 2207 2345 2304 2267 2394 2188 2252 2327 2252 2360 2275 2254 2308 2282 2267
2301 2189 2307 2278 2265 2319 2283 2156 2217 2258 2305 2174 2242 2214 2194 2259 2253 2200 2221 2225 2254 2202 2130 2208 2213 2138
2286 2254 2160 2170 2156 2138 2197 2188 2126 2187 2256 2121 2086 2114 2145 2185 2111 2099 2133 2145 2121 2113 2137 2167 2154 2161
2176 2008 2130 2131 2103 2111 2080 2053 2175 2030 2011 2097 2087 2058 2100 2120 2060 2040 1990 2057 2045 2039 1995 2133 2053 1947
2026 1943 2076 1931 2047 2047 2036 1969 1982 1986 2029 2021 2024 1993 1996 2018 1982 1935 1914 1948 2007 1890 2025 2079 1964 1982
2032 1951 1997 1837 2022 1881 1830 1943 1950 1862 1916 2003 1901 1869 1874 1920 1906 1896 1911 1933 1890 1810 1927 1858 1816 1905
1891 1825 1842 1934 1870 1945 1851 1904 1859 1905 1833 1846 1854 1875 1893 1797 1802 1857 1836 1798 1828 1874 1839 1827 1821 1818
1842 1825 1853 1761 1765 1841 1751 1899 1775 1800 1838 1855 1748 1793 1759 1780 1771 1730 1740 1876 1745 1751 1772 1820 1760 1755
1750 1782 1756 1729 1774 1741 1755 1803 1748 1749 1776 1716 1699 1771 1752 1748 1733 1787 1704 1718 1720 1690 1727 1628 1693 1665
1715 1681 1714 1753 1735 1658 1692 1701 1690 1711 1705 1639 1682 1665 1706 1784 1632 1706 1694 1682 1623 1625 1559 1571 1620 1639
1629 1685 1610 1614 1609 1604 1591 1660 1662 1628 1628 1660 1573 1643 1583 1543 1569 1563 1644 1551 1555 1624 1574 1561 1600 1570
1633 1495 1614 1571 1527 1585 1541 1566 1508 1579 1609 1510 1558 1608 1546 1594 1550 1575 1466 1514 1490 1558 1573 1507 1514 1562
1556 1556 1495 1552 1499 1509 1562 1492 1466 1526 1461 1464 1421 1486 1504 1467 1485 1544 1535 1567 1482 1442 1546 1487 1537 1494
1476 1501 1526 1545 1442 1485 1478 1500 1479 1445 1435 1464 1402 1419 1423 1446 1444 1437 1469 1457 1458 1453 1467 1465 1404 1438
1528 1400 1439 1481 1491 1481 1430 1515 1366 1397 1420 1446 1368 1448 1429 1470 1429 1400 1370 1392 1408 1463 1393 1382 1356 1448
1312 1378 1352 1368 1445 1386 1410 1331 1355 1402 1420 1331 1321 1325 1353 1360 1442 1374 1317 1393 1370 1323 1315 1406 1297 1352
1348 1334 1339 1339 1393 1390 1312 1365 1385 1298 1386 1349 1269 1353 1352 1305 1328 1335 1300 1283 1346 1348 1344 1325 1308 1271
1336 1346 1319 1332 1305 1352 1350 1276 1267 1298 1317 1272 1300 1281 1277 1316 1283 1289 1242 1270 1303 1302 1279 1289 1311 1232
1271 1262 1262 1276 1298 1313 1302 1284 1303 1278 1223 1252 1198 1259 1219 1289 1235 1209 1208 1196 1211 1282 1163 1198 1209 1173
1274 1221 1222 1241 1258 1230 1189 1216 1249 1208 1223 1278 1211 1231 1212 1237 1225 1205 1164 1215 1244 1138 1182 1212 1162 1164
1142 1243 1218 1244 1186 1146 1186 1173 1168 1153 1156 1150 1169 1189 1119 1142 1130 1147 1166 1179 1156 1097 1162 1118 1148 1138
1190 1083 1144 1142 1148 1154 1162 1128 1144 1130 1122 1179 1131 1101 1126 1158 1080 1146 1139 1123 1095 1163 1154 1118 1057 1163

1165 1140 1158 1143 1118 1056 1085 1095 1131 1061 1048 1037 1117 1176 1084 1133 1064 1137 1159 1047 1114 1132 1064 1071 1158 1090
1098 1122 1084 1084 1063 1044 1075 1014 1096 1048 1142 1081 1063 1034 1103 1079 1083 1040 1036 1107 1011 1101 1051 1102 1056 1066
1063 1041 1074 1048 1090 1028 1017 1046 1042 1027 1066 1072 1043 1035 1011 1019 1018 1054 1051 1019 1027 1003 1002 1017 1030 1089
1021 983 1038 982 1060 1028 979 1021 1029 1005 1034 1058 1020 999 991 1010 1077 969 1062 1000 1002 1048 1023 961 1006 985 1037 1003
1000 1036 1030 962 945 1013 973 945 985 978 981 1051 955 951 957 942 960 1003 988 968 956 1003 976 966 952 944 954 962 1003 1043 945
1023 980 964 896 966 951 945 946 970 971 968 915 953 953 1032 910 955 882 896 992 971 955 980 925 973 941 950 937 957 925 938 959 932
910 977 859 896 950 890 920 940 931 962 946 920 963 880 872 878 957 911 906 901 938 948 894 900 916 910 911</ChannelData>
 </TimeIntervalAnalysis>
 </ProcessedNeutronData>
</MeasurementData>
</CondensedNeutronDataFile>