LA-UR- 04-4191

Approved for public release; distribution is unlimited.

Title:

STUDY OF THE IDGS TECHNIQUE FOR MIXED PLUTONIUM-URANIUM (MOX) SAMPLES

Author(s):

T. K. Li, D. T. Vo, M. Sumi, T. Suzuki, Kobayashi, and Ohnishi

Submitted to:

45th Annual INMM Meeting Orlando, FL USA July 18-22, 2004 (ABSTRACT)





Study of the IDGS Technique for Mixed Plutonium-Uranium (MOX) Samples

M. Sumi, T. Suzuki, Kobayashi and Mr. Ohnishi Plutonium Fuel Center, Tokai Works Japan Nuclear Cycle Development Institute (JNC) Tokai-mura, Ibaraki-Ken, Japan

> T. K. Li and D. Vo Los Alamos National Laboratory Group N-1, MS E-540 Los Alamos, NM 87545 USA

ABSTRACT

The isotope dilution gamma-ray spectrometry (IDGS) technique has been demonstrated for simultaneously measuring concentrations and isotopic compositions of plutonium in spent-fuel input dissolver solutions. For timely analyzing nuclear materials on the purpose of material accountancy and quality control/assurance, we have performed a feasibility study to implement the IDGS for measuring mixed plutonium-uranium oxide (MOX) samples at the Plutonium Fuel Center (PFC) of Japan Nuclear Cycle Development Institute (JNC). Proof-of-principle experiments and analysis have been conducted for developing simultaneous plutonium and uranium measurements in MOX samples with wide variation of Pu/U ratios including powder, pellets and process scraps from the MOX fuel fabrication plant at PFC.

1. INTRODUCTION

In the IDGS method, the concentration of plutonium in the unknown solution is determined by calculating the differences among the isotopic ²⁴⁰Pu/²³⁹Pu ratios of the spike, the spiked solution, and the unknown solution.

$$C(Pu) = \frac{M_s}{V_u} \bullet \frac{W_s^9}{W_u^9} \bullet \frac{R_m - R_s}{R_u - R_m}, \qquad (1)$$

where

 $M_{\rm s}$ = mass of plutonium in the spike

= volume of dissolver solution taken

= weight fraction of ²³⁹Pu in the spike

= weight fraction of 239 Pu in the unknown solution = the 240 Pu/ 239 Pu ratio in the mixed or spiked solution = the 240 Pu/ 239 Pu ratio in the spike

= the ²⁴⁰Pu/²³⁹Pu ratio in the unknown solution

In this equation, the values of M_s , V_u , W_s^9 , and R_s are known. Only the values of R_u and W_u^9 in the unknown solution and R_m in the spiked solution need to be measured by gamma-ray spectrometry [1,2].




Figure 1. Deviation of the FRAM analytical results from Mass spectrometry results.

From figure 1, we can see that the isotopic results behave expected: they scatter wildly from the expected values when the acquisition times are small (weak statistics) and they all converge to the expected values as more data are obtained, including ²⁴²Pu. Note that ²⁴²Pu is determined from correlation and is not expected to be very accurate. We did not modify the ²⁴²Pu correlation for these data and it is just coincidence ²⁴²Pu results agree that the excellently with the MS results. For typical measurements, we may expect the systematic error for ²⁴²Pu correlation to be about 5-10%.

Table 2 shows the ratio of the gamma rays spectrometry results to that of the MS results of the 70-h spectra, including the FRAM's reported uncertainties.

We see that the results of ²³⁸Pu, ²³⁹Pu, and ²⁴⁰Pu agree well with the

MS results, well within the uncertainties. The errors for ²⁴²Pu are from statistics only and do not include the systematical errors from correlations. From other FRAM analyses, we know that the FRAM reported uncertainty for ²⁴¹Pu is about a factor of two too small. The true errors for ²⁴¹Pu for the results in Table 2 should be about 0.5% instead of 0.20% or 0.25% as shown. With 0.5% uncertainty then the ²⁴¹Pu results would agree with that of the MS within 2 sigmas.

Table 2. Ratios of FRAM's isotopic results to that of the mass spectrometry for the 70-h spectra. The errors are the 1 sigma uncertainties reported by FRAM.

1	Pu238		Pu239		Pu240		Pu241		Pu242		U235	
Samples	FRAM	%										
	/MS	error	/MS	error	/ MS	error	/MS	error	/MS	error	/ MS	error
#1	0.9971	0.46	1.0003	0.10	0.9998	0.27	0.9914	0.20	1.0060	0.47	0.9922	0.48
#2	1.0002	0.58	1.0006	0.12	0.9996	0.34	0.9912	0.25	0.9982	0.59	0.9818	0.60

As for the ²³⁵U results, where the results are several sigmas away from the MS results, it could very well be that the 2 sets of branching ratios for plutonium and uranium used in FRAM are

materials from different sources. That is, each plutonium material from a type of reactor (BWR, PWR, CANDU, etc.) or reactor burning different fresh fuel (uranium or MOX) will have the ²⁴²Pu correlations optimized for that material instead of the current one-size-fit-all ²⁴²Pu correlation formula.

CONCLUSION

We have shown that FRAM can be used with the IDGS technique to simultaneously determine plutonium and uranium isotopic compositions and concentrations in MOX samples at PFC, JNC. The uncertainties of the results are somewhat large due to weak statistics. If better statistics are obtained by either using more plutonium in the measurements, acquire the data for longer time, or using higher efficiency detector then the results can be better. The accuracy of the results can also be improved by a factor of 2-3 by using the generalized IDGS technique [7] instead of this traditional IDGS.

REFERENCES

- 1. T. K. Li, Y. Kuno, S. Sato, K. Nakatsuka, and T. Akiyama, "Determination of Plutonium Concentration and Isotopic Compositions by Isotope Dilution Gamma-Ray Spectroscopy on Resin Beads," *Nulc. Mater. Manage.* XIX (Proc. Issue), 390–396 (1990).
- T. K. Li, J. L. Parker, Y. Kuno, S. Sato, A. Kurosawa, and T. Akiyama, "Development of Isotope Dilution Gamma-Ray Spectrometry for Plutonium Analysis," Proc. 13th ESARDA Symposium on Safeguards and Nuclear Materials Management (ESARDA, Avignon, France, 1991) 24, 175–185.
- 3. T. K. Li, "A Novel Technique to Overcome Difficult Problems in Simultaneously Measuring Plutonium and Uranium in Spent-Fuel Dissolver Solutions," Los Alamos National Laboratory report LA-UR-99-1988 (1999), presented at ESARDA 21st Annual Meeting Symposium on Safeguards and Nuclear Materials Management, Palacio de Exposiciones y Congresos, Sevilla, Spain, 1999.
- 4. T. K. Li, T. Kuno, O. Kitagawa, S. Sato, A. Kurosawa, and Y. Kuno, "Simultaneous Measurements of Plutonium and Uranium in Spent-Fuel Dissolver Solutions," Los Alamos National Laboratory document LA-UR-97-2715 (1997).
- 5. T.E. Sampson, G.W. Nelson, and T.A. Kelley, "FRAM: A Versatile Code for Analyzing the Isotopic Composition of Plutonium from Gamma-Ray Pulse Height Spectra," Los Alamos National Laboratory report LA-11720-MS (1989).
- 6. T.E. Sampson, T.A. Kelley, and D.T. Vo, "Application Guide to Gamma-Ray Isotopic Analysis Using the FRAM Software," Los Alamos National Laboratory report LA-14018 (2003).
- 7. D.T. Vo and T.K. Li, "Generalization of the IDGS technique," Los Alamos National Laboratory report LA-UR-04-???? (2004), to be presented at the INMM 45th Annual Meeting, Orlando, Florida, USA, 2004.