

# LEGIBILITY NOTICE

A major purpose of the Technical Information Center is to provide the broadest dissemination possible of information contained in DOE's Research and Development Reports to business, industry, the academic community, and federal, state and local governments.

Although a small portion of this report is not reproducible, it is being made available to expedite the availability of information on the research discussed herein.

Los Alamos National Laboratory is operated by the University of California for the United States Department of Energy under contract W-7405-ENG-36

LA-UR--87-2837

DE87 014761

**TITLE** DESIGN OF A TRANSURANIC VUV SPECTROMETER

**AUTHOR(S)** A. J. Arko

**SUBMITTED TO** Proceedings of the Conference on Synchrotron Radiation  
Instrumentation, June 22-25, 1987, Madison, WI

**DISCLAIMER**

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

By acceptance of this article the publisher recognizes that the U S Government retains a nonexclusive royalty-free license to publish or reproduce the published form of this contribution or to allow others to do so for U S Government purposes

The Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U S Department of Energy

**MASTER**

**Los Alamos** Los Alamos National Laboratory  
Los Alamos, New Mexico 87545

DISTRIBUTION OF THIS DOCUMENT IS UNLIMITED  
BY THE U S GOVERNMENT

## DESIGN OF A TRANSURANIC VUV SPECTROMETER

A. J. Arko

Los Alamos National Laboratory, Los Alamos, NM 87545

The physics of 5f materials includes such interesting solid state phenomena<sup>1</sup> as itinerant and localized magnetism, heavy fermion superconductivity, strongly interacting Kondo systems, and similar phenomena that can be associated with narrow 5f bands. The radial extent of the 5f wave function in the light actinides is sufficient to form narrow 5f bands (in contrast to 4f systems), while in the heavy actinides the 5f levels are localized. An actinide contraction similar to that of the 4f series occurs in the latter.

The crossover point between localized and delocalized behavior in pure metals occurs somewhere around neptunium or plutonium. The many allotropic phases<sup>2</sup> in these metals suggest participation of the nearly localized 5f electrons in the bonding, which surely must be of very short range. Thus it is anticipated that the most interesting narrow band phenomena would occur in the transuranic materials Np and Pu and their compounds. In spite of the current emphasis on high- $T_c$  superconductors there is likely to be a substantial interest in transuranic research in the future.

Most of our present knowledge about the physics of 5f systems derives from studies of uranium and its compounds. This is particularly true in the case of photoemission measurements where the intense radioactivity has prevented studies at synchrotron sources. The development of a transuranic VUV spectrometer capable of safe operation at a synchrotron source would represent a giant step

in 5f research. We present below a conceptual design of such a spectrometer.

The major hazard of actinide materials is intense  $\alpha$ -radioactivity which is easily stopped by any structural material (even a layer of air). Actinides are harmful only if ingested. The whole safety concept of a transuranic system then is containment of samples in order to prevent ingestion of even minute quantities. Any transuranic apparatus of necessity consists of two parts. The first is mechanical and involves the various safety features which are designed to fail safe in case of an accident. The second, however, is a set of strict operating procedures which must be adhered to for safe operation in order to prevent accidents in the first place. Both are important. We first consider the hardware.

Figure 1 shows a schematic layout of a VUV spectrometer and sample preparation system capable of handling transuranics. The entire system is envisioned as a single portable unit which detaches from the beamline when not in use. The spectrometer chamber is kept as a clean system where both "hot" (radioactive) and "cold" (non-radioactive) samples can be measured. It is connected to the beamline via a small slit which serves the dual purpose of preventing migration of small particles up the beamline, and as a delay line (allowing time to close fast valves) in case of a catastrophic loss of vacuum. The "hot" sample preparation chamber is entirely enclosed in a glovebox where all surface preparation of "hot" samples is done. A series of right angle transfers are required to introduce the "hot" samples into the spectrometer chamber. This again serves a dual purpose. First the right angles prevent a direct line-of-sight into the spectrometer chamber, and second, the three transfers vastly improve the chances of keeping the spectrometer chamber

contamination-free since each subsequent transfer line is "colder" than the previous transfer line. Only the sample itself and the sample holder are "hot" in the spectrometer.

Pre-cleaned "hot" samples are introduced into the preparation chamber through an interlock system, as shown in Fig. 1. They are brought to the synchrotron from the home lab in a sealed UHV container and transferred into the glovebox via a magnetically-coupled sample rod so that at no time are they exposed to the atmosphere. All roughing pump lines are exhausted through a series of filters to the roof. Key valves in the transfer lines are interlocked so that one cannot open to the beamline when the prep chamber is open to the spectrometer.

Other safety features not shown in the figure include: a) Fast-acting valves in the beamline that are interlocked to the spectrometer vacuum; b) Radiation monitoring devices strategically located in the spectrometer and in the beamline; c) Safety covers on all optical ports; d) Inert atmosphere at a slight under-pressure in the glovebox; e) External dry active waste disposal container for all suspect materials; f) External radiation counters; g) Water cooling of glovebox to dissipate heat during bakeout of sample prep chamber; h) Ultrasonic vibrator in prep chamber to remove loose dust from scraped samples prior to transferring. While the system described above should be able to withstand safely a catastrophic vacuum failure; it is nonetheless essential to have a set of Strict Mandatory Operation Procedures as an integral part of any transuranic system. A complete operating manual should be developed. It should contain, at a minimum, the basic procedures outlined here: 1) The entire system must be limited to sub-gram quantities of material so that even in a "worst case

scenario" we are dealing only with minute quantities of material; 2) Samples are brought to the ring just prior to measurement and removed immediately upon completion of experiment; i.e., no long-term storage of samples; 3) Samples are stored in a double containment mode, except in the spectrometer chamber during actual measurement when a qualified operator is present; 4) System is confined to stable materials only - single crystals when possible. This excludes gas-phase work; 5) All surface preparation is done in the glovebox-enclosed sample preparation chamber; 6) Samples are pre-cleaned in the home lab; 7) Strict inventory of all materials is to be maintained; 8) Frequent (possibly daily) external and internal radiation monitoring to detect problems before they develop; 9) Monthly check of exhaust filters; 10) Cleaning of spectrometer chamber in the home lab whenever background  $\alpha$ -radiation exceeds  $\approx 10^4$  dpm (most likely once a year).

Similar systems have operated successfully at Los Alamos and at Karlsruhe.<sup>3</sup> Substantial contamination was accumulated in the Karlsruhe experimental chamber because of the single transfer rod which served the dual purpose of transfer rod and sample manipulator. Our multiple transfer system overcomes the problems encountered there. Indeed, with the implementation of all the precautions noted above, there is no reason why safe operation at a synchrotron source would not be possible.

## References

1. See for example "The Handbook on the Physics and Chemistry of the Actinides," A. J. Freeman and G. H. Lander eds (North Holland, Amsterdam).
2. J. L. Smith, Z. Fisk, and S. S. Hecker, in Proceedings of the Conference on Electronic Structure and Properties of Rare Earth and Actinide Intermetallics, St. Polten, Austria, 1984, pp. 151-158.
3. J. R. Naegele, J. Ghijsen and L. Manes, in Structure and Bonding 59/60 (Springer-Verlag, Berlin, Heidelberg) 1985, pp. 198-262.

### Figure Captions

Fig. 1 Schematic layout of a transuranic VUV spectrometer, emphasizing the essential features of such a system.

1  
2  
3



