

LA-UR-12-26024

Approved for public release; distribution is unlimited.

Title:	Micro X-ray Fluorescence for Spatially Resolved Elemental Materials Characterization
Author(s):	Havrilla, George J
Intended for:	Invited seminar at NIST Boulder



Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by the Los Alamos National Security, LLC for the National Nuclear Security Administration of the U.S. Department of Energy under contract DE-AC52-06NA25396. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

Micro X-ray Fluorescence for Spatially Resolved Elemental Materials Characterization

George J. Havrilla
NIST Boulder
November 8, 2012



Operated by Los Alamos National Security, LLC for NNSA

UNCLASSIFIED

LAUR 12-XXXXX



Overview



- Introduction – XRF and elemental imaging
- Experimental - optics
- Instrumentation
 - Polycapillary – 2D EDAX Eagle III
 - DCC – 2D hiRX R&D instrument
 - Polycapillary – 3D R&D instrument
- Matrices – simple to complex
- Summary
- Future Directions



UNCLASSIFIED



XRF Process



MXRF

- Spatially resolved elemental analysis
- Spatial resolution down to 10 micrometers with commercial X-ray optics
- Limited sensitivity due to energy filtering by polycapillary optic
- Multi-element mapping
- Commercial laboratory based instrumentation









UNCLASSIFIED



MXRF Capabilities

- Single point analysis, spatially resolved, quantitative
- Line scans - concentration gradients, qualitative
- Maps or elemental imaging - large area heterogeneity, PIWTA, qual and quant



UNCLASSIFIED

Elemental Imaging


- Spatially resolved elemental characterization
- Single point spectra, line scans and elemental maps/elemental imaging
- X-ray optics need to be selected for the application
 - Spatial resolution
 - Sensitivity
 - Spectral range/resolution

Optics

- **Polycapillary 2D** – EDXRF EDAX Eagle III MXRF instrument, 40 kV, 1 mA, Rh, ~50 micrometer spot size
- **DCC (doubly curved crystal) 2D** – research instrument, 50 kV, 1 mA, Rh, 200 μm
- **Polycapillary 3D** – research instrument, 50 kV, 0.5 mA, Ag, 30 x 30 x 60 μm



Polycapillary optics

- Shaped bundle of glass capillaries
- Polychromatic
- High X-ray intensity
- High sensitivity
- Compromised spatial selectivity
- Rapid data collection



Monolithic Polycapillary X-ray Optics

Polycapillary




EDAX Eagle III MXRF



Pu contaminated soil particles on filter


white light image



P1 - magnified white light image




autoradiograph




800 x 600 x 400 micrometers

EDXRF Elemental Maps of Soil Particles on Filter




EDXRF Pu map of P1 particle



50 micrometer
spot size
200 ms dwell
per point

EDXRF spectra of P1 soil particle



187 ng Pu
in particle

MXRF

- Commercially available lab-based instrumentation
- Low power X-ray tube sources
- Vendor software for elemental mapping
- Good sensitivity and spatial resolution


hiRX

- high Resolution X-ray – based on MWDXRF technology developed for direct compositional characterization of spent nuclear fuel
- MWDXRF – (monochromatic wavelength dispersive XRF) utilizes doubly curved crystals (DCC) to focus and monochromatize the X-rays passing through the optic both for excitation and detection

hiRX


- DCC optics – focus and monochromatize X-rays
Based on Bragg diffraction
 - Select both excitation and detection energies
- MWDXRF – monochromatic wavelength dispersive XRF
 - Monochromatic excitation using Rh Ka line 20.2 keV
 - Monochromatic detection – selected analyte target Pu at 14.28 keV

DCC X-ray optics




- Monochromatic optics
- Si, Ge, quartz, graphite or mica crystal material
- 10-20% reflectivity
- Useful for X-rays of energy 1.5 keV-35 keV

hiRX Schematic Diagram



hiRX Prototype Instrument




hiRX Technology

- Doubly curved crystal (DCC) optics for excitation and collection
- Small spot excitation – several hundred micrometers
- Small sample requirements – 1 microliter or less
- Collection optic to reject background and collect only analyte signal


Simple matrix - Aqueous Pu Drop

270 μg Pu deposit


white
light
image of
dried
residue



EDAX
PuLa
image
with 50
 μm spot
size




EDXRF spectrum Pu Drop #1




hiRX spectrum of Pu Drop #1

hot spot



Moderate matrix - Soil particles on filter


white light image



P1 - magnified white
light image




autoradiograph



800 x 600 x 400 micrometers


EDXRF and hiRX Elemental Maps of P1

EDXRF




50 micrometer spot size
200 ms dwell per point

hiRX




200 micrometer spot size
3 s dwell per point

hiRX-EDXRF spectra of P1 soil particle




Moderate matrix - Soil particles on filter P2


white light image




autoradiograph



hiRX elemental map of P2 particle



hiRX spectrum P2 particle



No EDXRF
spectrum


Complex matrix – Spiked Synthetic Spent Fuel

Element	X-Ray Energies	Element	X-Ray Energies
AgL	2.98	NiK	7.47
AsK	10.54	OsL	8.91
AuM	2.12	PdL	2.83
BaL	4.46	PrL	5.03
BiM	2.42	PtM	2.05
CdL	3.13	RbK	13.39
CeL	4.84	ReL	8.65
CsL	4.28	RhL	2.69
CuK	8.04	RuL	2.55
DyL	6.49	SbL	3.6
ErL	6.94	SeK	11.22
EuL	5.84	SmL	5.63
GaK	9.25	SnL	3.44
GdL	6.05	SrK	14.16
GeL	1.18	TaL	8.14
HfL	7.89	TaM	1.71
HoL	6.72	TbL	6.25
I L	3.937	TeL	3.76
InL	3.28	TmL	7.18
IrL	9.17	W L	8.39
LaL	4.65	W M	1.775
LuL	7.65	Y K	14.95
MnK	5.89	YbL	7.41
MoL	2.29	ZnK	8.6
NdL	5.23	ZrL	2.04

nominal
concentration of
each element is
100 ng


U is spiked at 10 μg
Pu at 30 ng

Spiked Synthetic Spent Fuel




Spike deposit is
~1800
micrometers in
diameter

hiRX elemental map of 30 ng Pu spike in synthetic spent fuel



Spectra overlay hiRX and EDXRF



30 ng Pu spike



DCC optics

- Monochromatic excitation – low background, high intensity
- Monochromatic collection/detection – high sensitivity, high selectivity
- Single element mapping

Confocal Geometry



Confocal MXRF


Beam dimensions at 30 kV, 0.3 mA

z: $\sim 65 \mu\text{m}$


y: $\sim 30 \mu\text{m}$

x: $\sim 30 \mu\text{m}$


* Based on the first derivative *



beam is ellipsoidal



3D Image of Pu P1 Particle




orange - Fe
red - Pu




Confocal 3D images P1 particle

Top view



Side view



Confocal

- 3D nondestructive elemental imaging
- Depth penetration from 10's to 1000's of micrometers depending upon analyte and matrix
- Multi-element depth dependent imaging
- Unique elemental characterization technique

Optic Comparison

	Polycapillary	DCC	Monocapillary	Aperture
Smallest Spatial Resolution			✓	
Highest Flux (cps)		✓		
Highest Flux Density (cps/area)	✓			
Fastest Data Acquisition	✓			
Lowest Detection Limits		✓		

Applications

- Catalysts before and after use
- Fossil visualization
- Bone and tissue mapping
- Meteorite composition
- Minerals
- Photovoltaics, batteries, fuel cells
- Integrated spectrometry with Raman and FTIR

Summary

- X-ray optics provide unique capabilities for spatially resolved nondestructive elemental analyses
- No one optic provides complete capabilities of spatial resolution, energy transmission and sensitivity
- Focus on laboratory based instrumentation, for MXRF

Summary

- Select appropriate X-ray optics to optimize data acquisition
- Polycapillary offers good spatial resolution, limited by energy filtering
- hiRX DCC offers highest sensitivity, limited to single element detection

Acknowledgements

- The support of the Next Generation Safeguards Initiative (NGSI), Office of Nonproliferation and International Security (NIS), National Nuclear Security Administration (NNSA).
- Velma Montoya
- Christopher Galvin
- Hakim Boukhalfa
- Michael Collins

