

## LA-UR-12-26909

Approved for public release; distribution is unlimited.

Title: Applications of Neutron Diffraction Protein Crystallography

Author(s): Gonzalez, Javier M.


Intended for: Invited speaker at Institute of Molecular and Cell Biology of Rosario, Argentina.



### Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by the Los Alamos National Security, LLC for the National Nuclear Security Administration of the U.S. Department of Energy under contract DE-AC52-06NA25396. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

# Applications of Neutron Diffraction Protein Crystallography



**IBR-CONICET**


**Rosario, Argentina. 12-17-2012**

Javier M. González  
Bioscience Division, PCS, LANL.

LA-UR-XX-XXXXX

UNCLASSIFIED

# Neutron vs. X-ray scattering



Atom	Coherent Scattering	Incoherent Scattering
<b>H</b>	<b>-3.7</b>	<b>80.3</b>
<b>D</b>	<b>6.7</b>	<b>2.1</b>
<b>C</b>	<b>6.6</b>	<b>0.0</b>
<b>O</b>	<b>5.8</b>	<b>0.0</b>
<b>N</b>	<b>9.4</b>	<b>0.5</b>
<b>S</b>	<b>2.8</b>	<b>0.0</b>

# The Protein Crystallography Station at LANSCE

**Team:**


**PI: Cliff Unkefer**

**Instrument Scientist: Zoe Fisher**

**Technologist: Mary Jo Waltman**

**Postdocs: Oksana Gerlits, John-Paul Bacik, Javier Gonzalez**

**Senior Advisor: Benno Schoenborn**



# Neutron Xtallography: Pros & Cons




- No radiation damage
- Elucidation of *all* H atom positions
- Protonation states, H-bonds and solvent structures ( $\text{OH}^-$ ,  $\text{H}_2\text{O}$ ,  $\text{H}_3\text{O}^+$ ,  $\text{H}_5\text{O}_2^+$  ...)
- Information obtained from deuterium exchange (solvent accessibility, dynamics)
- Still in development





- Needs neutrons (\$\$\$)
- Usually at room temperature
- Sample (per-) deuteration
- Crystal size of  $\text{mm}^3$
- Available neutron sources are few and weak (long data collection, poor statistics)
- Diffraction quality and unit cell limitations
- Still in development






$\sim 1 \text{ mm}^3$ , 32 hrs exposure ( $2.0 \text{ \AA}$  resolution)

# Bigger is better





$\sim 30 \text{ mm}^3$ , 12 hrs exposure ( $1.8 \text{ \AA}$  resolution)


# Typical nuclear density maps





His




Arg



Lys




Gln



D<sub>2</sub>O

# Protein Crystallography Station (PCS)

- PCS is a high performance, state-of-the-art neutron crystallography beamline (DOE-OBBER funded)
- Unique - 1st detector of its kind built at a spallation source
- User program and instrument commissioned in '02
- Accommodate 25-50 users a year, 15-25 proposals a year



UNCLASSIFIED




# Protein Crystallography Station (PCS)


- $^3\text{He}$  position-sensitive neutron detector (BNL)
- $120^\circ$  coverage
- K-circle goniometer
- Hüber motors
- Usable wavelength range  $0.6 - 7.0 \text{ \AA}$  (TOF, Laue)
- $> 80\%$  counting efficiency (very sensitive)
- Relatively insensitive to H background (can collect H/D samples with great success)
- Smallest done here so far  $0.3 \text{ mm}^3$ , most are in the  $1 - 7 \text{ mm}^3$  range



# Support for Users




BDL (Bio-Deuteration Laboratory)




PPCL (Protein Purification and Crystallization Lab)


Grow large crystals of target proteins



PCS - Collect and process neutron and X-ray diffraction data




Support for using joining neutron and X-ray software




# PCS Support Labs

- Protein purification facilities
- Protein crystallization facilities
- Liquid handling robot
- In-house Rigaku HighFlux X-ray home source
- Biodeuteration lab
- Algae farm to make Altone
- **Data processing software (*n*CNS, phenix.refine)**



# PCS Support Labs



**Figure 3**

A schematic representation of the process used for perdeuterating proteins at the BDL at Los Alamos National Laboratory (LANL). A bacterial culture medium, which we designate Altone, is made from the hydrolyzate of algae (e.g. *Scenedesmus obliquus*) grown in D<sub>2</sub>O and used for protein expression in *E. coli*. Steps represented by framed boxes can be performed robotically at LANL. Steps represented by shaded boxes can be performed either at LANL or the user's home laboratory.

Paul Langan et al. *J. Synchrotron Rad.* (2008). **15**, 215–218



# Neutron structures determined at the PCS



Macromolecule	Space group	Resolution, Å	Completeness, %	Crystal size, mm <sup>3</sup>	Data Collection days	PDB ID	Citation
Crambin	$P2_1$	1.1	78.8	~5	21	---	Chen <i>et al.</i> , 2012
Rubredoxin, W3Y mutant	$P2_12_12_1$	2.1	72	~2	5	---	Li <i>et al.</i> , 2004
<i>E. coli</i> DHFR/Methotrexate	$P6_1$	2.2	80	0.3	23	2INQ	Bennett <i>et al.</i> , 2006
XI-2Co <sup>2+</sup>	$I222$	1.8	78	8	19	2GVE	Katz <i>et al.</i> , 2006
Z-DNA, d(CGCGCG)	$P2_12_12_1$	1.6	62	0.7	18	---	Langan <i>et al.</i> , 2006
PYP	$P6_3$	2.5	89	0.8	14	2QWS	Fisher <i>et al.</i> , 2007
Endothiapepsin/ <i>gem</i> -diol inhibitor	$P2_1$	1.8	80	2.7	37	2VS2	Tuan <i>et al.</i> , 2007; Coates <i>et al.</i> , 2008
Deoxy human hemoglobin (HbA)	$P2_1$	2.0	87	20	18	3KMF	Kovalevsky <i>et al.</i> , 2008 & 2010
XI-2Mg <sup>2+</sup> -D-xylulose	$I222$	2.2	89	4	37	3CWH	Kovalevsky <i>et al.</i> , 2008
Carbonic anhydrase (pH=10)	$P2_1$	2.0	85	1.2	55	3KKX	Fisher <i>et al.</i> , 2009; Fisher <i>et al.</i> , 2010
DFPase	$P2_12_12_1$	2.2	82	0.43	37	3BYC	Blum <i>et al.</i> , 2009
Amicyanin	$P2_1$	1.8	68	2.6	21	3L45	Sukumar <i>et al.</i> , 2005; Sukumar <i>et al.</i> , 2010
XI-2Ni <sup>2+</sup> -D-glucose	$I222$	1.8	84	15	20	3KCO	Kovalevsky <i>et al.</i> , 2010
Equine CNmet hemoglobin	$C2$	2.0	80	10	28	---	Kovalevsky <i>et al.</i> , 2010
Human glycosyltransferase	$C222_1$	2.5	84	0.4	14	---	Schuman <i>et al.</i> , 2011
Bovine $\gamma$ -chymotrypsin	$P4_22_2$	2.0	85	1.5	12	---	Lazar <i>et al.</i> , 2011
<i>apo</i> -XI (pH=5.9)	$I222$	2.0	87	9	10	3QZA	Kovalevsky <i>et al.</i> , 2011
Xylanase	$P2_12_12_1$	1.8	88	7	15	---	Kovalevsky <i>et al.</i> , 2011
Z-DNA	$P2_12_12_1$	1.4	-	-	-	3QBA	Fenn <i>et al.</i> , 2011
HCA II + acetazolamide	$P2_1$	2.0	85	2.0	22	---	Aggarwal <i>et al.</i> 2012 in preparation
HCA II (pH=7.8)	$P2_1$	2.0	85	1.7	20	3TMJ	Fisher <i>et al.</i> , 2011

# Publications

Structure 18, 688–699, June 9, 2010


## Metal Ion Roles and the Movement of Hydrogen during Reaction Catalyzed by D-Xylose Isomerase: A Joint X-Ray and Neutron Diffraction Study

Andrey Y. Kovalevsky,<sup>1,\*</sup> Leif Hanson,<sup>2</sup> S. Zoe Fisher,<sup>1</sup> Marat Mustyakimov,<sup>1</sup> Sax A. Mason,<sup>3</sup> V. Trevor Forsyth,<sup>3,4</sup> Matthew P. Blakeley,<sup>3</sup> David. A. Keen,<sup>5</sup> Trixie Wagner,<sup>1,7</sup> H.L. Carrell,<sup>6</sup> Amy K. Katz,<sup>6</sup> Jenny P. Glusker,<sup>6</sup> and Paul Langan<sup>1,2,\*</sup>



UNCLASSIFIED

# Publications




Angew. Chem. Int. Ed. 2011, 50, 7520–7523

DOI: 10.1002/anie.201101753


## X-ray/Neutron Crystallography

### Identification of the Elusive Hydronium Ion Exchanging Roles with a Proton in an Enzyme at Lower pH Values\*\*

Andrey Y. Kovalevsky,\* B. L. Hanson, S. A. Mason, T. Yoshida, S. Z. Fisher,  
M. Mustyakimov, V. T. Forsyth, M. P. Blakeley, D. A. Keen, and Paul Langan\*



pH < 6  
pH ≈ 8



# Publications

*J. Am. Chem. Soc.* 2012, 134, 14726–14729



## Neutron Diffraction of Acetazolamide-Bound Human Carbonic Anhydrase II Reveals Atomic Details of Drug Binding

S. Zoë Fisher,<sup>†</sup> Mayank Aggarwal,<sup>‡</sup> Andrey Y. Kovalevsky,<sup>†</sup> David N. Silverman,<sup>§</sup> and Robert McKenna<sup>\*,‡</sup>

<sup>†</sup>Bioscience Division, Los Alamos National Laboratory, Los Alamos, New Mexico 87545, United States

<sup>‡</sup>Department of Biochemistry and Molecular Biology, University of Florida, P.O. Box 100245, Gainesville, Florida 32610, United States

<sup>§</sup>Department of Pharmacology and Therapeutics, University of Florida, P.O. Box 100247, Gainesville, Florida 32610, United States



**Figure 2.** Stick representation of AZM bound to HCA II. Zinc is shown as a magenta sphere, D atoms are in cyan (indicated by arrows), and H atoms in white (positions calculated). The nuclear density map is shown in yellow and is contoured at 1.5  $\sigma$ , and the electron density map is shown in blue and is contoured at 2.0  $\sigma$ .




## Direct observation of hydrogen atom dynamics and interactions by ultrahigh resolution neutron protein crystallography

Julian C.-H. Chen<sup>a,1</sup>, B. Leif Hanson<sup>a</sup>, S. Zoë Fisher<sup>b</sup>, Paul Langan<sup>a,c</sup>, and Andrey Y. Kovalevsky<sup>b</sup>

<sup>a</sup>Department of Chemistry, University of Toledo, Toledo, OH 43606; <sup>b</sup>Bioscience Division, Los Alamos National Laboratory, Los Alamos, NM 87544; and

<sup>c</sup>Biology and Soft Matter Division, Oak Ridge National Laboratory, Oak Ridge, TN 37831



**Fig. 5.** Resolving a hydrogen bond network. (A) 0.85 Å X-ray electron density map showing potential H bond partners for W1022. (B) 1.1 Å nuclear density map showing that W1022 interacts with W1023, but not W1007, W1024, or W1026. (C) Hydration environment around Y44. A potential O—D... $\pi$  interaction is highlighted involving W1023. Distances are indicated in Å.

# Take-home messages

- Neutron crystallography has the potential to describe molecular details untraceable through ultra-high resolution X-ray crystallography
- Different isotopes scatter neutrons differently (e.g.  $^1\text{H}/^2\text{H}$  exchange), increasing the experimental possibilities
- Accurate description of protein-ligand interactions is particularly useful for Structural Enzymology
- The User Program at PCS offers a full support for protein neutron crystallography studies