

LA-UR-13-25116

Approved for public release; distribution is unlimited.

Title: Pressure Study of Savillex Vessels

Author(s): Colletti, Lisa M.
Dickson, Frank P. III

Intended for: Group meeting Program reviews

Issued: 2013-07-10



Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by the Los Alamos National Security, LLC for the National Nuclear Security Administration of the U.S. Department of Energy under contract DE-AC52-06NA25396. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

70 YEARS OF CREATING TOMORROW



Los Alamos
NATIONAL LABORATORY

Pressure Study of Savillex Vessels

Lisa Colletti, Frank
Dickson

UNCLASSIFIED



Why did we decide to test our vessels


- Previous glass vessels under went a violent disintegration in September of 2011, despite the existing process documentation (1972 paper) that indicated only 80 psi was generated. Investigation revealed changes occurring to process over 20 years prior that unknowingly changed safety envelope.
- Wanted a through understanding of our process conditions to prevent conditions that might cause vessel disintegration again (once bitten twice shy).
- A recent conversation with Dave Gallimore indicated that PS had gone with this system solely on the manufacture's specifications of vessels and had not done any independent testing of vessels.
- Manufacture states that the vessels and caps have a 100 psi rating only up to 38 °C which then drops 1 psi for every degree over that. Dissolution methods are at 120-150 °C which theoretically drops the pressure rating to 18 psi to less than 0 psi.

UNCLASSIFIED





Manufacturer's Specification for Vessel Pressure and Rating vs Temp



UNCLASSIFIED





Room temperature testing


Pneumatic Testing

- Vessels pressurized with Argon gas at room temperature leaked at 20 psi and less
- All 3 valves tested leaked

Hydrostatic Testing


- Vessels filled and pressurized with water
- Bubbles seen at 20 psi and less, as low as 8 psi
- All 6 valves tested leaked

BUBBLES FROM
RELIEF PORT AT 20 PSI
(ROOM TEMP)



UNCLASSIFIED

Pressure vessel set up



Without PRV lid

With PRV lid





PRV lid cutaway

PFA DIAPHRAGM
CONTINUOUS WITH VALVE
STEM AND BODY

PLASTIC VALVE SPRING

VALVE STEM
SEALS WITH ORIFICE IN
LID


VALVE STOP
LIMITS SPRING TRAVEL




UNCLASSIFIED



Empty Vessel Heating Characteristics



THERMAL PROPERTIES	Conductivity	Specific Heat
	W/m*K	KJ/Kg*K
Air (20 C)	0.026	1.005
Graphite	400 - 1700	0.71
PFA	0.19	1.17




UNCLASSIFIED



How does the use of a hotplate affect temperature profile?


60 mL vessel, Hot plate heating



UNCLASSIFIED



Max Lid temperatures with acids in vessels




UNCLASSIFIED



Vessel Pressure vs Temp

5 and 10 mL 33% HCl Acid in 60 and 120 mL Vessel




UNCLASSIFIED



Vessel Pressure vs Temp



5 and 10 mL OPTIMA (69%) HNO₃ in 60 and 120 mL Sealed Vessel



UNCLASSIFIED



Comparison of sealed vs PRC




UNCLASSIFIED



What a full run looked like

10 mL Concentrated HCl in 120 mL Vessel with Pressure Relief Closure
(Volume Loss 1.3 mL)
(Initial Molarity 10 M, Final Molarity 8.8 M)



UNCLASSIFIED



Acid loss summary

Acid	Acid	Vessel	Filler	Time	Temp	Diff. in Acid Weight	Concentration, M		
	mL	mL	Block	hr	C	g	Initial	Final	Diff
Conc. HCl	5	60	Yes	5-Apr	150	0.75	10.3	8.7	1.6
Conc. HCl	5	60	No	5-Apr	150	0.49	10.3	9.5	0.8
Conc. HCl	5	60	Yes	5-Apr	100	0.46	10.3	10.3	0
Conc. HCl	5	60	No	5-Apr	100	0.33	10.3	10.3	0
Optima HNO3	5	60	No	5-Apr	150	0.04			
Optima HNO3	5	60	Yes	5-Apr	150	0.14	15.3	15.3	0
Optima HNO3	5	120	NA	5-Apr	150	0.26	15.3	15.3	0
Optima HNO3	10	120	NA	5-Apr	150	0.08	15.3	15.3	0
Optima HNO3	5	60	No	5-Apr	100	0.04			


UNCLASSIFIED





HNO₃, HCL, HF are all negative azeotropes

- Azeotropic point: The solution composition exactly matches the composition of vapor phase. You cannot separate the two components by basic distillation processes.
- for a system pressure of 1 bar
 - HNO₃: BP = 120.5°C, and a HNO₃ concentration of 68 wt%, (15.2M)
 - HCl: BP = 108.6°C, a HCl concentration of 20.2 wt% (6.15M).
 - HF BP = 120 °C , HF concentration of 35.35% (20M).



- By starting at 12 M HNO₃, PS is potentially distilling out water shifting concentration towards azeotropic point.
- By starting with 15.3 M HNO₃ in PA, concentration never changes even with solution loss.
- By starting with 12M HCl in PA, we distill out HCl lowering the HCl concentration in solution towards azeotropic point.

UNCLASSIFIED



Is solution loss a problem... YES!

- Some materials are volatile (SiF_4 , GaCl_3 , etc.) which would affect analyses of trace elements.
- Solution chemistry is critical during the dissolution of PuO_2
 - Use of HCl allows higher HF concentrations, a critical factor in dissolution of high fired materials or materials with large amounts of refractory materials.
 - Loss of HCl affects the solution chemistry by lowering acidity values and increasing free F^- concentrations.
 - Loss of either HCl/ HNO_3 will lower the volume and increase HF and F^- concentrations
 - Both of previous effects can drive formation of PuF_4 , an insoluble precipitate
 - Prevents true understanding of the dissolution process as the dissolution conditions change during the experiments.
- Loss of acid changes the matrix of the material being analyzed which is critical in
 - Separations
 - Matrix matching for spectrophotometric methods

UNCLASSIFIED





Can we just use HNO_3 ? Not desired...

- Though we can just use HNO_3 , this would cause massive disruptions to existing processes.
- The required changes to PA methods would cause significant increase in sample prep time and increased uncertainties.
 - We've seen this in our 238 spectroscopy system. A rough doubling of uncertainty has been seen since switching to HNO_3 for dissolution. On going work to reduce back to past levels.


UNCLASSIFIED





Options Going Forward

- Could potentially change to hotplate to have lower overall temperature on lid.
- Could find alternate vessels.
 - Q-tube disposable glass vessel
 - Pressure release at 200 psi, tube rated to 500 psi.
- Carry on with current method, but..
 - Some elements may be lost.
 - May require dissolution spikes.
 - May have reruns due to precipitation of PuF_4 .



UNCLASSIFIED