

## LA-UR-14-22067

Approved for public release; distribution is unlimited.

Title: Biomedical Applications of Microfluidic Technology

Author(s): Gao, Jun

Intended for: Presentation at New Mexico State University

Issued: 2014-03-28



### Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by the Los Alamos National Security, LLC for the National Nuclear Security Administration of the U.S. Department of Energy under contract DE-AC52-06NA25396. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

## Biomedical Applications of Microfluidic Technology


Jun Gao

B-10

Bioscience, LANL

### What are the microfluidic lab-on-chips?


- **Lab-on-chips** : A class of submillimetre size bioanalytical devices.
- **Perform**: fluidic processes, sensing, analysis and separation of biochemical samples.
- **Integrate**: fluidics, electronics, optics and biosensors.
- **Analyse**: metabolites, molecules, proteins, nucleic acids, cells and viruses.



### Advantages and Applications

- **Low fluid consumption** (less waste, lower reagents costs and less required sample volumes for diagnostics)
- **Faster analysis and response times** due to short diffusion distances, fast heating, high surface to volume ratios, small heat capacities.
- **Better process control** because of a faster response of the system (e.g. thermal control for exothermic chemical reactions)
- **Small footprint**
- **Massive parallelization** allows high-throughput analysis
- **Lower fabrication costs** allow mass production
- **Safer platform for chemical, radioactive or biological studies** because of integration of functionality, smaller fluid volumes and stored energies
- PCR
- Immunoassays
- drug screening
- cell biology
- Genomic & proteomics
- single cell manipulation
- Water & food quality
- Environmental monitoring
- Improvements in synthetic yields
- Novel reactions by control surface chemistry, local heat and mass transfer
- Controlling reagent concentrations


### Microfluidic Device Fabrication Using Polydimethylsiloxane (PDMS)




### Microfluidic Applications at Bioscience Current status

- Tissue engineering
  - Liver tissue/vascular tissue engineering
  - Control neuron differentiation
- Droplet-based microfluidic Technology
  - Cell and bacteria encapsulation in W/O droplet
  - Algae encapsulation in agar droplet
  - Double emulsion (W/O/W droplets)
- Applications
  - Digital PCR
  - Protein crystallization
- Future directions



### Microfluidic studies Develop Liver Microvascular Tissue



### Neuron Microfluidic Device




### Droplet-based Microfluidic Technology




The combination of aqueous phase is flexible. It can be chemical reaction solution, cell culture media with cells, agarose/beads and other buffers.

### Microfluidic droplet generator (W/O)




### Real-time Changing the Size of Droplet (50-480 $\mu\text{l/h}$ )




Real-time changing the droplet size: The droplet will contain different amount of reaction solution for reaction efficiency evaluation.


### Stability Evaluation




### Single Bacteria Encapsulation BL21 (DE3) Gold (pETck3)




### BL21 (DE3) Gold (pETck3) (Overnight incubation @ 37°C)




### *E. Coli* Proliferation in W/O Droplets




### Cell and Bacteria Coencapsulation



### Encapsulated Human Jurkat T Cells




### Phagocytosis (THP1 and *E. Coli*)




### Encapsulation Efficiency

(Poisson Distribution)




### Move from Aqueous Droplet to Hydrogel Droplet




### Nannochloropsis in Agar Droplets

#### In F2 media




#### In Oil




### Tetraselmis in Agar Droplet


#### In F2 media




#### In Oil




### Double Emulsion (W/O/W)



### Microfluidic PCR




### Microfluidic PCR






### mAb PCR Gel



1. Regular PCR
2. Regular PCR using droplets
3. Regular PCR w/o template
4. Regular PCR w/o template using droplets
5. Microfluidic PCR w/flow rate 20 $\mu$ l/hr
6. Microfluidic PCR w/o template 20 $\mu$ l/hr
7. Microfluidic PCR w/flow rate 200 $\mu$ l/hr
8. Microfluidic PCR w/flow rate 200 $\mu$ l/hr (leftover)


Sample loading in gel: 3 $\mu$ l + 3 $\mu$ l dye  
Reaction: 2.5 $\mu$ l template was used in 50 $\mu$ l (total volume)

### Digital RT-PCR Taqman assays in microfluidic droplets




PLOS, Vol 8 (4) 2961

### Principle of cell sorting




Nature Protocols 8, 870–891 (2013)

### Sorting and Detection



Nature Protocols 8, 870–891 (2013)

### Lysozyme crystallization in droplet (20hrs)



### Future

- Multiplexing
- Not just pathogen detection, but detection of virulence and antibiotic resistance markers
- Organism sensing and detections of biomarkers for infection together
- Point of care application in resource limited setting

### Acknowledgment

- Dr. Andrew Bradbury and his team
- Dr. Elizabeth-Hong Geller
- Dr. Scott Twary
- Dr. Jeri Sullivan