

LA-UR-18-28459

Approved for public release; distribution is unlimited.

Title: Characterization of EJ-301 liquid scintillators

Author(s): Redpath, Thomas Hundley

Meierbachtol, Krista Cruse

Baugher, Travis Ray

Smith, Karl

Intended for: Summer student presentation at home institution

Issued: 2018-09-05



Characterization of EJ-301 liquid scintillators

Thomas Redpath^{1,2}, Krista Meierbachtol¹, Travis Baugher¹, Karl Smith³



¹Los Alamos National Laboratory
Nuclear Engineering and Nonproliferation
NEN-2

²Michigan State University
Department of Physics and Astonomy

³Los Alamos National Laboratory
Space Science and Applications
ISR-1

Thomas Redpath (NEN-2)

- Educational Background
 - O James Madison University, 2011
 - O Central Michigan University, 2014
 - O Michigan State University, present







Division: NEN

O Group: NEN-2

O Krista Meierbachtol









- Research
 - O Characterization of liquid scintillators for nγ detection
 - O Lifetime measurement of neutron-unbound ²⁶O

Research Overview and Motivation

- Compact and portable detectors for neutron spectrometry
- ID neutron sources (e.g. ²⁴⁰Pu, ²⁵²Cf)
- Streamline data acquisition and analysis
- Calibration with gamma sources
- Translate detector response into an energy spectrum unfolding
- Accurate unfolding increases the sensitivity of assays performed on nuclear material





J.A. Mullens, J.D. Edwards, S.A. Pozzi, Analysis of scintillator pulse-height for nuclear material identification, in: Institute of Nuclear Materials Management 45th Annual Meeting, July 18–22, Orlando, FL, 2004.

Portable Detectors and Electronics

 Four Eljen EJ-301 liquid scintillators, bubble-free

PROPERTIES	EJ-301
Light Output (% Anthracene)	78
Scintillation Efficiency (photons/1 MeV e-)	12,000
Wavelength of Maximum Emission (nm)	425
Decay Time, Short Component (ns)	3.2
Mean Decay Times of First 3 Components (ns)	3.16 32.3 270
Bulk Light Attenuation Length (m)	2.5 - 3
Specific Gravity	0.874
Refractive Index	1.505
Flash Point (°C)	26
Boiling Point (°C at 1 atm)	141
Vapor Pressure (mm Hg, at 20°C)	6
H Atoms per cm ³ (×10 ²²)	4.82
C Atoms per cm³ (×10 ²²)	3.98
Electrons per cm³ (×10 ²³)	2.27





Caen DT5730 desktop digitizer




Caen DT5533 High Voltage Supply





Neutron-gamma Discrimination

 Caen Pulse Shape Discrimination firmware (DPP-PSD)



B. Wan *et al.*, Digital pulse shape discrimination methods for n-y separation in an EJ-301 liquid scintillation detector, Chinese Physics C, Vol 39, Issue 11, 2015






Calibration of light output – gamma sources

Compton Edges

	E _y [keV]	EC [keV]
Cs-137	661.6	477.3
Co-60	1173.2	963.4
Co-60	1332.5	1118.1
Co-60 (avg)	1252.9	1040.6
K-40	1460.8	1243.3
TI-208	2614.6	2381.8









Figure 1. Ideal linear response function: r(E), and the convolved response function: R(E).



M.J. Safari *et al.*, Differentiation method for localization of Compton edge in organic scintillation detectors, https://arxiv.org/abs/1610.09185


Neutron light deposited – ²⁵²Cf





J. W. Boldeman, D. Culley, and R. Cawley, "The Fission Neutron Spectrum from the Spontaneous Fission of ²⁵²Cf," *Transactions of the American Nuclear Society* 32, 733 (1979).

 $N(E) = \sqrt{E}e^{-E/T}; T = 1.43 \text{MeV}$




Michael T. Febbraro, University of Michigan dissertation, 2014

W. P. Poenitz and T. Tamura, "Investigation of the Prompt Neutron Spectrum for Spontaneously-Fissioning ²⁵²Cf," Proc. Int. Conf. Nucl. Data Sci. Techno., Antwerp, Belgium, Sept.


Research Approach

Four Eljen EJ-301 liquid scintillators

PROPERTIES	EJ-301
Light Output (% Anthracene)	78
Scintillation Efficiency (photons/1 MeV e-)	12,000
Wavelength of Maximum Emission (nm)	425
Decay Time, Short Component (ns)	3.2
Mean Decay Times of First 3 Components (ns)	3.16 32.3 270
Bulk Light Attenuation Length (m)	2.5 - 3
Specific Gravity	0.874
Refractive Index	1.505
Flash Point (°C)	26
Boiling Point (°C at 1 atm)	141
Vapor Pressure (mm Hg, at 20°C)	6
H Atoms per cm ³ (×10 ²²)	4.82
C Atoms per cm³ (×10 ²²)	3.98
Electrons per cm³ (×10 ²³)	2.27



 Caen Pulse Shape Discrimination firmware (DPP-PSD)





Caen DT-5730 desktop digitizer

