

LA-UR-18-29524

Approved for public release; distribution is unlimited.

Title: Thermal Studies of HMX Decomposition

Author(s): Salinas, Daniela Evodia

Intended for: Senior Seminar Class

Issued: 2018-10-05




THERMAL STUDIES OF HMX DECOMPOSITION

W-10, WEAPONS RESPONSE CPCS, PHYSICAL CHEMISTRY AND APPLIED SPECTROSCOPY

DANIELA SALINAS





HMX

- $C_4H_8N_8O_8$
- Powerful insensitive nitroamine high explosive
- Detonator in nuclear weapons
- HMX most effective chemical explosive
- High melting point and molecular weight

OBJECTIVE

- Understand the kinetics of the HMX decomposition process
- Isothermal tests at a new temperature range 205°C-280°C
- Contribute to the global chemical kinetic model of HMX high explosive (Henson)

TG/DSC METHODS

- Thermogravimetric/ Differential Scanning Calorimetry
- STA 409 PC Instrument
- Simultaneously measures weight change and heat flow



TG/DSC METHODS

Run a baseline measurement

Measurement Type: 'Correction'



Calibration measurement with sapphire

Measurement Type: 'Sample+Correction'



Sample measurement

Measurement Type: 'Sample+Correction'

OPEN AND CLOSED PANS

- Open
 - Decomposition reaction characteristics
 - Deflagration and detonation processes


- Closed
 - Confined spaces
 - Secondary reactions
 - Accelerated decomposition reactions

TEMPERATURE PROGRAM


- Initial (25°C)
 - Dynamic (10°C/min)
- Isothermal (205°C-285°C)
- Final (+20°C)
- Final Standby/ Emergency (+20°C)

- NETZSCH Proteus Analysis
 - Temperature Calibration
 - Sensitvity Calibration
 - Sensitivity Calibration via Cp


INDIUM, TIN, ZINC



Area under the curve specific enthalpy of fusion $(\mu Vs/mg)$




Endothermic and Exothermic Peaks Melting and Recrystallization




Hysteresis can be observed

CALIBRATIONS




Temperature Calibration: In, Sn, Zn

Sensitivity Calibration Curve: In Sn, Zn




CALIBRATIONS



Baseline Stability for Sapphire


Sensitivity Calibration via Cp




EXPERIMENTS


- Isothermal tests
 - 2-3mg HMX samples
 - 280°C,270°C,220°C,210°C,205°C
 - Constant heating rate, 10°C/min
 - Aluminum pans
 - Nitrogen atmosphere

OPEN AND CLOSED ALUMINUM PANS





220°C ISOTHERM







Open and closed pan 2-3mg samples, 220°C isotherm for 24 hours, 10°C/min heating ramp












TG curves of %mass loss over time







Time to ignition and reaction kinetics global chemical HMX model with different observables See: Henson and Smilowitz, IDS, (2010)

CONCLUSION


- Absorbing heat and losing mass for open (210°C-215°C) and closed (206°C-211°C) pans
- Half times 200°C-300°C range for open and closed pans
- Temperature range follows the linear trend for the global chemical kinetic model for HMX

FUTURE WORK

- Repeat experiments to identify impurities
- Preheat the furnace for an "instantaneous heating effect"
- PETN and TATB experiments









REFERENCES

- [1] Henson, B. F., et al. Modeling Thermal Ignition and the Initial Conditions for Internal Burning in PBX 9501.AIP Publishing, American Institute of Physics, 28 Dec. 2009.
- [2] Henson, B.F., et al. Evidence for Thermal Equilibrium in the Detonation of HMX. Shock Wave Science and Technology Reference Library, Vol. 5.
- [3]Burnham, A.K., and R.K. Weese. Thermal Decomposition Kinetics of HMX. Lawrence Livermore National Laboratory, 2005.
- [4]Singh, Arjun. Thermal Decomposition and Kinetics of Plastic Bonded Explosives Based on Mixture of HMX and TATB with Polymer Matrices. Egyptian Journal of Medical Human Genetics, Elsevier, 23 Dec. 2016.
- [5] Peterson, P. and Borovina D. 28 February 2017, *The HE Grand Challenge: LANL Ignition to Detonation Project*, Los Alamos National Laboratory. 19 July 2017.