

LA-UR-21-20554

Approved for public release; distribution is unlimited.

Title: Countroom Gamma Spectroscopy Library Updates


Author(s): James, Michael R.
Chong, Sarah Wing Ting

Intended for: Report

Issued: 2021-01-22

Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.



Countroom Gamma Spectroscopy Library Updates
21-Jan-2021
Michael James and Sarah Chong

The C-NR Countroom analyzes gamma spectroscopy data using three codes, SAP [1](which is the combination of SPECANAL and RAYGUN), GAMANAL [2] and UNISAMPO[3]. The codes use libraries which are chosen by the user of gamma emission data to do the identification and quantification. The Countroom maintains nine unique libraries for use in analysis. Historically, the values contained in these libraries has not been consistently updated and many values are quite old and do not reflect the state-of-the-art.

A comprehensive update of Gamma Spectroscopy Libraries has been undertaken using data from the from IAEA Livechart website [4] which has been confirmed to derive from latest Nuclear Data Sheets. The work to update all nine libraries was carried out in Summer 2020[5]. The SAP code data libraries are maintained in ASCII files which are human-readable and carries a RAY suffix. This file was chosen to the master record for all library updates. From this master file, the RGN program can be used to generated the library binary files which are used as input to SAP processing.

Almost all entries in each library required some change to update the values to modern data. In most cases these changes were minor. Larger changes are called out individually below grouped by library type. A value of 365.2422 was used to convert from years to days from IAEA Livechart website.

Each library was also tested running at least one spectrum with the old and new versions. The counts used are listed by the Countroom serial numbers. The output comparisons are listed in the Appendix.

GAMANAL and UNISAMPO Library Files

For the GAMANAL code, a spreadsheet was used to convert the data from the format in the RAY file to the elibx.in files (where x is the library number). These elibx.in files are then processed by the libgen program to create the GAMANAL libraries.

UNISAMPO uses an ASCII file as its library format with a suffix ILF. In order to update the UNISAMPO libraries, a Python script, sampo_lib.py, was written to create the library in the proper format using the listing which is generated by the RGN program from the original SAP library master file. In this way, the updates were applied universally to all three analysis codes.

Library 1 Changes

K-40 (19040) half-life change from 4.675E+11 to 4.555E+11
Ti-44 (22044) Significant change in half-life from 1.753E+4 to 2.159E+4 days
Cd-109 (48109) Changed 88.03 from 3.9% to 3.644%
Sn-113 (50113) Changed 255.1 line from 2.3% to 2.11%
Sb-124 (51124) 1325.08 line from 1.41% to 1.58%
Sb-126 (51126) Removed 605 – Not in IAEA LiveChart
Sb-126 (51126) Removed 730 - Not in IAEA LiveChart
Te-129m (52429) Half-life change to 33.6 from 4.833e-2 days(reflects metastable half-life)
Nd-147 (60147) Changed 91.1 line from 30.77% to 28.1%
Eu-152 (63152) Changed half-life from 5.128E+3 to 4.937E+3 days
Eu-152 (63152) Removed 443.98 – Not in IAEA LiveChart, merged with adjacent line 0.2 eV away
Eu-154 (63154) Changed half-life from 2.864E+3 to 3.141e+3 days
Eu-154 (63154) Remove 1397.8 line – not in IAEA LiveChart
Eu-155 (63155) Changed half-life from 1.855E+3 to 1.735e+3 days
Eu-156 (63156) Almost every line yield reduced by 10-13%.
Ra-224 (under 90228) 240.99 line increased yield from 2.24% to 4.1%

Baseline count	Test Count
248_509.19	248_510.19
071_486.20	071_487.20
196_487.20	196_488.20

Library 2 Changes

Te-131m (52431) Half-life changed from 1.25d to 1.385d
Te-131m (52431) 452.30 intensity changed from 5.65% to 1.50%
Tc-99m (43399) 140.5 line changed from 80.14% to 89%
I-126 (53126) Removed 656 line – not present on Livechart
I-126 (53126) Removed 664.5 line – not present on Livechart
I-126 (53126) Removed 730 line – not present on Livechart
Xe-133 (54133) 79.6 line changed from 0.22% to 0.44%.
Pu-239 (94239) Removed 42.06 line – not present on Livechart

059_096.20 Baseline count (reprocessed from 059_095.20) using old library 2
059_098.20 Test count with new library 2.

Library 3 Changes

Se-81m (34381) Half changed to isomeric state from 1.278E-02d to 3.978E-02d
102.7 line changed to 103.01 and intensity from 1% to 12.8%
K-42 (19042) 312 line intensity changed from 0.18% to 0.336%
Ca-47 (20047) 1878 line intensity changed from 2.90% to 2.50E-02%
Br-82 (35082) 92 line intensity changed from 0.39% to 0.719%
Ag-110n (47410) changed to 667.9 at 0.75% to 677.62 10.7% (better line to use)

Sb-124 (51124) changes:

400.3 line from 0.530% to 0.139%
444.09 line from 0.350% to 0.189%
525.5 line from 0.310% to 0.138%
1580 line from 0.200% to 0.380%

Eu-154 (63154) changes:

69.67 3 line from 7.80% to 2.42%
97.43 3 line from 56.5% to 29.0%
103.18 3 line from 41.3% to 21.1%

Yb-169 (70169) changes:

515 line from 4.50E-01% to 4.17E-03%
625 line from 4.00E-01% to 4.92E-03%

Yb-175 (70175) changes:

113.81 line from 1.91% to 3.87%
144.86 line from 3.60E-01% to 6.72E-01%
282.52 line from 3.10% to 6.13%
396.33 line from 6.50% to 13.1%

Hf-181 (72181) changes:

136.86 line from 1.8% to 0.86%
475.99 line from 0.43% to 0.703%

Pt-197 (78197) 77.35 line from 80% to 17%

As-76 (33076) Removed 554.45 not present on Livechart
As-76 (33076) Removed 624.00 not present on Livechart
Br-82 (35082) Removed 451.00 not present on Livechart
Br-82 (35082) Removed 1410.00 not present on Livechart
Nb-97 (41097) Removed 480.90 not present on Livechart

Ag-110m (47410) Removed 565.50 not present on Livechart
Ag-110m (47410) Removed 752.70 not present on Livechart
Ag-110m (47410) Removed 867.00 not present on Livechart

010_001.20 Baseline Count old library 3
010_098.20 Test count with new library 3

Library 4 Changes

Set intensity of Rb-82 776 line to 13.4% for Isotope Production program.
Nb-90 (41090) Removed 1522.00 not present on Livechart
Co-55 (27055) 91.9 line from 3.00% to 1.16%
Ni-65 (28065) 1623.42 line from 0.760% to 0.498%
As-76 (35076) 1213.1 line from 0.620% to 1.70%
As-77 (35077) 585.48 line from 1.15% to 1.57%
Kr-79 (36079) 136.09 line from 0.600% to 0.850%
Kr-79 (36079) 217.07 line from 1.80% to 2.37%
Rb-82 (37082) 1395 line from 0.47% to 0.57%
Nb-96 (41096) 350.05 line from 1.11% to 1.54% (added two lines together)
Nb-91m (41391) 1204.67 line from 3.40% to 2.02%

Add Ge-68/Ga-68 parent-daughter with
1077.34 3.22%
1883.16 0.137%

Two tests of Library 4

Baseline Counts	Test Counts
161_081.20	161_087.20 (New LIB4)
210_063.20	210_072.20 (New LIB4)

Library 5 Changes

In previous version of Library 5, numerous half-lives set from long-lived precursors.

Pb-209 (82209) set from Ac-225, 9.92d
At-217 (85217) set from Ac-225, 9.92d
Fr-221 (87221) set from Ac-225, 9.92d
Bi-213 (83213) set from Ac-225, 9.92d

Pb-211 (82211) set from Ra-223, 11.4d
Bi-211 (83211) set from Ra-223, 11.4d
Po-215 (84215) set from Ra-223, 11.4d
Rn-219 (86219) set from Ra-223, 11.4d

Po-215 (84215) 438.9 frequency not listed, kept at 4e-4%

At-217 (85217) 257.9 frequency greatly reduced in new library, 2.87E-4% vs. 5.6e-4%

Ac-225 (89225) 83.23 x-ray line from 1.17% to 0.76%

86.1. x-ray line from 1.94% to 1.26%

These X-ray lines left in the library unchanged:

Th-2289 (90229) 100.13 5.5%

Pa-229 (91229) 105.60 6.91%

104_039.20 Baseline count

104_051.20 Test count (new lib 5)

Library 6 Changes

I-126 (53126) Removed 656 line – not present on Livechart

I-126 (53126) Removed 664.5 line – not present on Livechart

I-126 (53126) Removed 730 line – not present on Livechart

I-135 (53135) 526.56 line has no frequency on Livechart, kept at 13.32%

Replace Sb-127 (51127) 583.4 line at 2.91% with 543.5 line 2.90%

583.4 should have been 0.33%, might have been an error in previous library.

248_040.20 Baseline count

248_045.20 Test count (new lib 6)

Library 7 Changes

Only minor adjustments

196_087.20 Baseline count

196_088.20 Test count (new lib 7)

Library 8 Changes

Only minor adjustments

011_045.20 Baseline count

011_046.20 Test count (new lib 8)

Library 9 Changes

Nb-95m (41395) 204.12 line intensity reduced from 2.34% to 0.960%
Cd-115 (48115) 336.24 line changed from 50.1% to 46.02%
Fix Ru-106/Rh-106 (44106) half-life from 3.212e7 to 371.8 (seconds to days).
Xe-133 (54133) 160.61 line increased from 5.40E-02% to 0.107%
W-187 (74187) Almost all lines frequency increased ~25%
Po-215 (84215) no frequency listed for 438.90, use old value of 4.000E-04
Ac-227 (89227) Yields for all 4 listed lines (70, 99, 147 and 160) reduced by a factor of 4.
Th-232 (90232) 59 line changed to 63.81, yield changed from 0.19% to 0.263%
124 line changed to 140.88, yield changed from 0.043% to 0.021%.
Pa-231 (91231) 303 line increased from 1.71% to 2.30%
Pa-233 (91233) no frequency listed for 111.50 use old value of 6.000E-02
Pa-233 (91233) no frequency listed for 114.90 use old value of 2.000E-02
Pa-234m (91634) 766 line changed from 0.207% to 0.317%
U-235 (92235) 105.36 line not listed in gammas or X-rays but leave in as is (from ThKbeta1) with frequency 1.980E-02
U-237 (92237) no frequency listed for 114.09, use old value of 9.600E-02
U-238 (92238) 110.5 line changed to 113.5. Yield changed from 0.024% to 0.0102%.

056_076.20 Baseline count

056_083.20 Test count (new lib 9)

References

[1] SAP Operations Manual, LANL-INC11-001, April 1993

[2] Gunnink R, Niday JD (1972) UCRL-51061

[3] Aarnio P, Nikkinen M, Routti J (2012) UNISAMPO manual v2.5. Doletum, Helsinki

[4] <https://www-nds.iaea.org/relnsd/vcharthtml/VChartHTML.html>

[5] Chong, Sarah Wing Ting ; James, Michael R., "Gamma Spectroscopy Library Update", Virtual Student Symposium 2020, Los Alamos, New Mexico, United States, LA-UR-20-25368

APPENDIX

Comparisons of SPECANAL outputs using new libraries. The new libraries outputs are listed on the left side and old library outputs are on the right. These comparisons were used to confirm that the library worked and identified the same set of isotopes and computed approximately the same yields.

Library 1 Comparison

248_510_19_specanal.lis										248_509_19_specanal.lis									
Linux-SPECANL Version 1.0 (05-AUG-2016)					DATE:2020/07/15 15:01:02.231					Linux-SPECANL Version 1.0 (05-AUG-2016)					DATE:2020/07/14 23:10:01.895				
*****										*****									
SN: 248_510	SHOT: 9997	JUG: 71	MASS: 2	CTR: 56	ROT: 248.951	ZT: 91.708				SN: 248_509	SHOT: 9997	JUG: 71	MASS: 2	CTR: 56	ROT: 248.951	ZT: 91.708			
NCHAN: 4096 CHAN0: 180.0 CHAN1: 181.0 DELTA TIME: 157.180 GAIN: 0.499950 ZERO: 0.057 CD: 71										NCHAN: 4096 CHAN0: 180.0 CHAN1: 181.0 DELTA TIME: 157.180 GAIN: 0.499950 ZERO: 0.057 CD: 71									
*****										*****									
Nuclide	Energy	DPM at T-zero	Peak wght	Err (%)	Std dev	DPM at count time	Atoms at T-zero	Err (%)		Nuclide	Energy	DPM at T-zero	Peak wght	Err (%)	Std dev	DPM at count time	Atoms at T-zero	Err (%)	
CO 57	-	1.831E+04	-	-	-	1.226E+04	1.033E+10	3.9		CO 57	-	1.836E+04	-	-	-	1.227E+04	1.030E+10	4.0	
	121.95	1.821E+04	1.01	1.0	-0.1						121.95	1.826E+04	1.01	1.0	-0.1				
	136.31	1.965E+04	0.93	3.5	0.4						136.31	1.985E+04	0.93	3.5	0.4				
Y 88	-	7.633E+04	-	-	-	2.747E+04	1.690E+10	1.2		Y 88	-	7.658E+04	-	-	-	2.756E+04	1.696E+10	1.2	
	897.38	7.633E+04	1.00	1.2	0.0						897.38	7.658E+04	1.00	1.2	0.0				
CD109	-	7.819E+05	-	-	-	6.176E+05	7.503E+11	0.8		CD109	-	7.295E+05	-	-	-	5.771E+05	7.044E+11	0.8	
	87.96	7.819E+05	1.00	0.8	0.0						87.96	7.295E+05	1.00	0.8	0.0				
SN113	-	4.555E+04	-	-	-	1.768E+04	1.089E+10	9.6		SN113	-	4.594E+04	-	-	-	1.782E+04	1.098E+10	9.2	
	254.96	5.696E+04	0.80	8.6	0.7						254.96	5.230E+04	0.88	8.6	0.4				
	391.42	4.541E+04	1.00	1.2	-0.0						391.42	4.585E+04	1.00	1.2	-0.0				
TE129M	-	1.565E+02	-	-	-	6.115E+00	1.093E+07	1.7		SB126	-	2.350E+00	-	-	-	4.129E-04	6.152E+04	1.7	
	278.97	1.955E+07	0.00	1.7	0.0						278.97	1.304E+09	0.00	1.7	0.0				
CS137	-	2.235E+04	-	-	-	2.213E+04	5.103E+11	1.0		CS137	-	2.235E+04	-	-	-	2.213E+04	5.089E+11	1.0	
	661.18	2.235E+04	1.00	1.0	0.0						661.18	2.235E+04	1.00	1.0	0.0				
CE139	-	2.626E+04	-	-	-	1.190E+04	7.507E+09	0.9		CE139	-	2.624E+04	-	-	-	1.189E+04	7.501E+09	0.9	
	165.69	2.626E+04	1.00	0.9	0.0						165.69	2.624E+04	1.00	0.9	-0.0				
HG203	-	5.531E+04	-	-	-	5.336E+03	5.354E+09	1.7		EU154	-	3.530E+04	-	-	-	3.398E+04	2.100E+11	12.8	
	278.97	5.531E+04	1.00	1.7	0.0						846.16	3.530E+04	1.00	12.8	0.0				
TH228	-	1.681E+04	-	-	-	1.439E+04	2.439E+10	8.4		HG203	-	5.535E+04	-	-	-	5.340E+03	5.358E+09	1.7	
	510.55	1.681E+04	1.00	8.4	-0.0						278.97	5.535E+04	1.00	1.7	0.0				
U237	-	5.063E-03	-	-	-	4.951E-10	7.100E+01	1.5		TH228	-	1.782E+04	-	-	-	1.525E+04	2.587E+10	8.4	
	59.48	5.727E+11	0.00	1.5	0.0						510.55	1.782E+04	1.00	8.4	0.0				
AM241	-	5.385E+04	-	-	-	5.381E+04	1.768E+13	1.5		U237	-	5.063E-03	-	-	-	4.951E-10	7.100E+01	1.5	
											59.48	5.727E+11	0.00	1.5	0.0				
+-- 2 lines: 59.48 5.385E+04 1.00 1.5 0.0										+-- 2 lines: 59.48 5.385E+04 1.00 1.5 0.0									

071_487_20_specanal.lis										071_486_20_specanal.lis									
Linux-SPECANL Version 1.0 (05-AUG-2016)					DATE:2020/07/15 15:00:02.895					Linux-SPECANL Version 1.0 (05-AUG-2016)					DATE:2020/07/14 23:11:02.147				
*****										*****									
SN: 071_487 SHOT: 9997 JUG: 71 MASS: 205 CTR: 172 ROT: 71.634 ZT: 91.708					NCHAN: 4096 CHAN0: 2485.6 CHAN1: 2514.5 DELTA TIME: 344.052 GAIN: 0.499700 ZERO: 0.436 CD: 41					SN: 071_486 SHOT: 9997 JUG: 71 MASS: 205 CTR: 172 ROT: 71.634 ZT: 91.708					NCHAN: 4096 CHAN0: 2485.6 CHAN1: 2514.5 DELTA TIME: 344.052 GAIN: 0.499700 ZERO: 0.436 CD: 41				
*****										*****									
Nuclide	Energy	DPM at T-zero	Peak wght	Err (%)	Std dev	DPM at count time	Atoms at T-zero	Err (%)		Nuclide	Energy	DPM at T-zero	Peak wght	Err (%)	Std dev	DPM at count time	Atoms at T-zero	Err (%)	
K 40	1460.63	2.184E+03	1.00	9.8	0.0	2.184E+03	2.067E+18	9.8		K 40	1460.63	2.176E+03	1.00	9.8	0.0	2.176E+03	2.113E+18	9.8	
CO 57	122.05	8.232E+04	1.01	0.5	-0.1	3.447E+04	4.681E+10	0.7		CO 57	122.05	8.278E+04	1.01	0.5	-0.1	3.458E+04	4.691E+10	0.7	
	136.45	8.379E+04	0.99	0.6	0.1						136.45	8.489E+04	0.99	0.6	0.2				
CO 60	1173.17	1.671E+05	1.00	0.5	-0.0	1.476E+05	6.681E+11	0.6		CO 60	1173.17	1.666E+05	1.00	0.5	-0.0	1.476E+05	6.681E+11	0.6	
	1332.45	1.675E+05	1.00	0.5	0.0						1332.45	1.675E+05	1.00	0.5	0.0				
Y 88	897.99	3.586E+05	1.00	0.6	-0.0	3.828E+04	7.941E+10	0.7		Y 88	897.99	3.591E+05	1.00	0.6	0.0	3.834E+04	7.952E+10	0.8	
	1836.17	3.591E+05	1.00	0.7	0.0						1836.17	3.585E+05	1.00	0.7	-0.0				
CD109	88.11	4.810E+06	1.00	0.5	0.0	2.870E+06	4.615E+12	0.5		CD109	88.11	4.479E+06	1.00	0.5	0.0	2.682E+06	4.325E+12	0.5	
SN113	255.14	2.448E+05	0.89	2.3	0.7	2.700E+04	5.127E+10	2.8		SN113	255.14	2.222E+05	0.97	2.3	0.2	2.713E+04	5.155E+10	2.5	
	391.69	2.132E+05	1.01	0.6	-0.1						391.69	2.154E+05	1.00	0.6	-0.0				
SB124	1325.08	3.134E+06	1.00	3.3	0.0	5.966E+04	3.920E+11	3.3		SB124	1325.08	3.512E+06	1.00	3.3	-0.0	6.685E+04	4.392E+11	3.3	
CS137	661.62	1.031E+05	1.00	0.5	0.0	1.009E+05	2.355E+12	0.5		CS137	661.62	1.031E+05	1.00	0.5	0.0	1.009E+05	2.348E+12	0.5	
CE139	165.85	1.249E+05	1.00	0.5	0.0	2.208E+04	3.571E+10	0.5		CE139	165.85	1.248E+05	1.00	0.5	0.0	2.206E+04	3.568E+10	0.5	
EU154	691.99	7.034E+03	1.00	24.1	0.0	6.519E+03	4.590E+10	24.1		EU154	691.99	7.631E+03	1.00	24.1	0.0	7.021E+03	4.540E+10	24.1	
HG203	279.21	2.681E+05	1.00	1.3	0.0	1.605E+03	2.595E+10	1.3		HG203	279.21	2.683E+05	1.00	1.3	0.0	1.606E+03	2.597E+10	1.3	
RA226	185.74	2.312E+03	1.00	19.1	0.0	2.311E+03	2.807E+12	19.1		RA226	185.74	2.364E+03	1.00	19.1	0.0	2.363E+03	2.873E+12	19.1	
TH228	238.74	3.042E+02	0.88	15.3	0.3	1.903E+02	3.883E+08	25.7		TH228	238.74	2.960E+02	0.92	15.3	0.2	1.941E+02	3.964E+08	23.9	
	583.06	2.317E+02	1.16	20.0	-0.3						583.06	2.462E+02	1.11	20.0	-0.2				

Left file: C:\Users\121611\Desktop\010_098_specanal.lis
 Right file: C:\Users\121611\Desktop\010_001_specanal.lis

Library 3 Comparison

Show spaces / tabs

Report type: All (Matching and Differences)

Summary: 20 Line(s) match | 14 Block(s) diff | 18 : 23 Line(s) diff

Ignore options: Ignore line terminators

C:\Users\121611\Desktop\010_098_specanal.lis										C:\Users\121611\Desktop\010_001_specanal.lis																							
1	*	FILE:010_098.LIS	SAP Version 4.2 (30-May-2008)				DATE:13-SEP-2020 17:50:07.92				1	*	FILE:010_001.LIS	SAP Version 4.2 (30-May-2008)				DATE:10-JAN-2020 02:01:10.84															
			PAGE: 1											PAGE: 1																			
2	!	>															2																
3			DPM at	Peak	Err	Std	DPM at	Atoms at	Err	3			DPM at	Peak	Err	Std	DPM at	Atoms at	Err														
4		Nuclide	Energy	T-zero	wght	(%)	dev	count	time	T-zero	(%)	4		Nuclide	Energy	T-zero	wght	(%)	dev	count	time	T-zero	(%)										
5		-----											5		-----																		
6												6																					
7	*	MN 54	- - - -	6.987E+00	- - - -	- - - -	- - - -	5.451E+00	4.532E+06	34.8	7																						
8											8	*	MN 54	- - - -	6.985E+00	- - - -	- - - -	- - - -	5.451E+00	4.535E+06	34.8												
9	*		834.84	6.987E+00	1.00	34.8	0.0				9																						
10											10			834.88	6.985E+00	1.00	34.8	0.0															
11	*	CO 58	- - - -	1.519E+05	- - - -	- - - -	- - - -	5.090E+04	2.237E+10	4.0	11																						
12											12	*	CO 58	- - - -	1.519E+05	- - - -	- - - -	- - - -	5.087E+04	2.235E+10	3.9												
13	*		810.72	1.523E+05	1.00	0.5	0.0				13																						
14	*		863.90	1.414E+05	1.07	2.0	-0.5				14			810.76	1.523E+05	1.00	0.5	0.0															
15	*		1674.68	1.810E+05	0.84	3.3	0.9				15			863.94	1.406E+05	1.08	2.0	-0.6															
16											16	*		1674.75	1.766E+05	0.86	3.3	0.8															
17		CO 60	- - - -	3.421E+01	- - - -	- - - -	- - - -	3.286E+01	1.368E+08	12.4	17																						
18											18	*	CO 60	- - - -	3.421E+01	- - - -	- - - -	- - - -	3.286E+01	1.368E+08	12.4												
19	*		1173.10	3.627E+01	0.94	10.9	0.2				19																						
20	*		1332.51	3.288E+01	1.04	9.6	-0.1				20			1173.15	3.627E+01	0.94	10.9	0.2															
21	!	>															21			1332.57	3.288E+01	1.04	9.6	-0.1									
22	!	>															22																
23	!	>															23	*	BA131	- - - -	1.628E+06	- - - -	- - - -	- - - -	2.289E+03	3.991E+10	9.5						
24	!	>															24																
25	!	>															25			92.57	1.628E+06	1.00	9.5	0.0									
26											26																						
27	*	EU152	- - - -	1.811E+02	- - - -	- - - -	- - - -	1.783E+02	1.858E+09	3.6	27	*	EU152	- - - -	1.822E+02	- - - -	- - - -	- - - -	1.793E+02	1.839E+09	3.6												
28											28																						
29	*		122.01	1.811E+02	1.00	3.6	0.0				29			122.02	1.822E+02	1.00	3.6	0.0															
30	*	TM170	- - - -	4.199E+02	- - - -	- - - -	- - - -	2.299E+02	1.122E+08	27.6	30																						
31											31	*	TM170	- - - -	3.230E+02	- - - -	- - - -	- - - -	1.749E+02	8.474E+07	27.6												
32	*		84.15	4.199E+02	1.00	27.6	0.0				32																						
33											33			84.16	3.230E+02	1.00	27.6	0.0															
34	*	YB169	- - - -	4.217E+02	- - - -	- - - -	- - - -	3.750E+01	2.805E+07	38.8	34																						
35											35	*	YB169	- - - -	4.522E+02	- - - -	- - - -	- - - -	3.623E+01	2.884E+07	38.9												
36	*		63.31	3.684E+02	1.14	14.8	-0.4				36																						
37	*		197.73	3.400E+03	0.12	12.0	2.5				37			63.31	3.962E+02	1.14	14.8	-0.4															
38											38			197.74	3.766E+03	0.12	12.0	2.5															
39											39																						
40	*	HF181	- - - -	1.127E+03	- - - -	- - - -	- - - -	1.811E+02	9.926E+07	19.5	40	*	HF181	- - - -	1.080E+03	- - - -	- - - -	- - - -	1.737E+02	9.517E+07	19.5												
41											41																						
42	*		136.57	1.127E+03	1.00	19.5	0.0				42			136.58	1.080E+03	1.00	19.5	0.0															
43											43																						

Left file: C:\Users\121611\Desktop\161_087_specanal.lis
 Right file: C:\Users\121611\Desktop\161_081_specanal.lis

Show spaces / tabs

Library 4 Comparison

Report type: All (Matching and Differences)

Summary: 14 Line(s) match | 5 Block(s) diff | 7 : 6 Line(s) diff

Ignore options: Ignore line terminators

C:\Users\121611\Desktop\161_087_specanal.lis										C:\Users\121611\Desktop\161_081_specanal.lis									
1	*	FILE:161_087.LIS	SAP Version 4.2 (30-May-2008)	DATE:13-SEP-2020	17:24:07.76					1	FILE:161_081.LIS	SAP Version 4.2 (30-May-2008)	DATE: 9-JUN-2020	20:16:09.50					
			PAGE: 1										PAGE: 1						
2	<!									2									
3										3									
4		ANALYTICAL RESULTS USING LIBRARY #4								4		ANALYTICAL RESULTS USING LIBRARY #4							
5	0	NUCLIDE	DPM AT	PEAK	PEAK	DPM AT	ATOMS AT			5	0	NUCLIDE	DPM AT	PEAK	PEAK	DPM AT	ATOMS AT		
6		ENERGY	ZERO TIME	FRAC	ERR	COUNT TIME	ZERO TIME	ERROR		6		ENERGY	ZERO TIME	FRAC	ERR	COUNT TIME	ZERO TIME	ERROR	
7										7									
8	*	1 AS 71	- - - -	4.361E+08	- - - -	7.945E+03	2.465E+12	21.2		8	*	1 AS 71	- - - -	4.294E+08	- - - -	7.186E+03	2.408E+12	21.2	
9										9									
10	*		175.	4.361E+08	1.00	21.2				10	*		175.	4.294E+08	1.00	21.2			
11										11									
12	*	2 SR 82	- - - -	5.136E+06	- - - -	1.592E+06	2.704E+11	2.5		12	*	2 SR 82	- - - -	5.163E+06	- - - -	1.599E+06	2.717E+11	2.5	
13										13									
14	*		777.	5.156E+06	1.00	2.4				14	*		777.	5.158E+06	1.00	2.4			
15	*		1396.	4.436E+06	1.00	17.6				15	*		1396.	5.383E+06	1.00	17.6			
16										16									
17		3 SR 85	- - - -	1.489E+06	- - - -	9.422E+05	2.006E+11	1.7		17		3 SR 85	- - - -	1.489E+06	- - - -	9.422E+05	2.006E+11	1.7	
18										18									
19			514.	1.489E+06	1.00	1.5				19			514.	1.489E+06	1.00	1.5			
20										20									
21										21									

Left file: C:\Users\121611\Desktop\210_072_specanal.lis

Right file: C:\Users\121611\Desktop\210_063_specanal.lis

Report type: All (Matching and Differences)

Summary: 9 Line(s) match | 3 Block(s) diff | 12 : 3 Line(s) diff

Ignore options: Ignore line terminators

C:\Users\121611\Desktop\210_072_specanal.lis									C:\Users\121611\Desktop\210_063_specanal.lis														
1	*	FILE:210_072.LIS	SAP Version 4.2 (30-May-2008)			DATE:13-SEP-2020 17:18:08.34			1	FILE:210_063.LIS	SAP Version 4.2 (30-May-2008)			DATE:28-JUL-2020 18:41:07.94									
			PAGE: 1								PAGE: 1												
2	<!								2														
3									3														
4									4														
5			DPM at	Peak	Err	Std	DPM at	Atoms at	Err	5			DPM at	Peak	Err	Std	DPM at	Atoms at	Err				
6		Nuclide	Energy	T-zero	wght	(%)	dev	count	time	T-zero	(%)	6		Nuclide	Energy	T-zero	wght	(%)	dev	count	time	T-zero	(%)
7		-----										7		-----									
8												8											
9	<!	GE 68	- - - -	5.300E+06	- - - -	- - - -	- - - -	4.751E+06	2.983E+12	2.5	9												
10	<!										10												
11	<!	1077.75	7.273E+06	0.73	2.5	0.0					11												
12	<!										12												
13	<!	RB 86	- - - -	3.234E+06	- - - -	- - - -	- - - -	6.591E+05	1.252E+11	2.5	13												
14	<!										14												
15	<!	1077.75	1.192E+07	0.27	2.5	0.0					15												
16	<!										16												
17	*	TC 96	- - - -	4.079E+08	- - - -	- - - -	- - - -	4.000E+05	3.626E+12	29.7	17	*	TC 96	- - - -	4.062E+08	- - - -	- - - -	- - - -	3.984E+05	3.612E+12	29.7		
18											18												
19	*	314.29	4.079E+08	1.00	29.7	0.0					19	*	314.29	4.062E+08	1.00	29.7	0.0						
20											20												
21											21												

Left file: C:\Users\121611\Desktop\196_088_specanal.LIS
 Right file: C:\Users\121611\Desktop\196_087_specanal.LIS

Show spaces / tabs

Library 7 Comparison

Report type: All (Matching and Differences)

Summary: 10 Line(s) match | 4 Block(s) diff | 16 : 14 Line(s) diff

Ignore options: Ignore line terminators

C:\Users\121611\Desktop\196_088_specanal.LIS										C:\Users\121611\Desktop\196_087_specanal.LIS																
1	*	FILE:196_088.LIS	SAP Version 4.2 (30-May-2008)			DATE:22-SEP-2020 20:08:08.32				1	*	FILE:196_087.LIS	SAP Version 4.2 (30-May-2008)			DATE:21-SEP-2020 18:44:08.27										
2	<!									2																
3										3																
4										4																
5			DPM at	Peak	Err	Std	DPM at	Atoms at	Err	5			DPM at	Peak	Err	Std	DPM at	Atoms at	Err							
6		Nuclide	Energy	T-zero	wght	(%)	dev	count time	T-zero	(%)	6		Nuclide	Energy	T-zero	wght	(%)	dev	count time	T-zero	(%)					
7											7															
8											8															
9	*	EU152	-	-	-	1.172E+06	-	-	-	-	4.305E+05	1.202E+13	24.9	9	*	EU152	-	-	-	1.169E+06	-	-	-	4.300E+05	1.201E+13	24.8
10											10															
11	*		121.85	1.176E+06	1.00	0.7	0.0				11	*		121.85	1.172E+06	1.00	0.7	0.0								
12	*		244.75	1.236E+06	0.95	1.7	0.4				12	*		244.75	1.229E+06	0.95	1.7	0.4								
13	*		344.30	1.137E+06	1.03	1.3	-0.3				13	*		344.30	1.139E+06	1.03	1.3	-0.2								
14			411.15	1.155E+06	1.01	5.7	-0.1				14			411.15	1.155E+06	1.01	5.7	-0.1								
15	*		443.98	1.117E+06	1.05	6.2	-0.2				15	*		443.98	1.107E+06	1.06	6.2	-0.2								
16	*		778.85	1.132E+06	1.04	4.9	-0.2				16	*		778.85	1.130E+06	1.04	4.9	-0.2								
17	*		867.31	1.034E+06	1.13	10.2	-0.4				17	*		867.31	1.029E+06	1.14	10.2	-0.4								
18	*		964.03	1.184E+06	0.99	5.1	0.0				18	*		964.03	1.174E+06	1.00	5.1	0.0								
19	*		1085.86	1.203E+06	0.97	7.0	0.1				19	*		1085.86	1.189E+06	0.98	7.0	0.1								
20	*		1112.03	1.174E+06	1.00	6.0	0.0				20	*		1112.03	1.175E+06	1.00	6.0	0.0								
21	*		1213.08	1.172E+06	1.00	24.2	0.0				21	*		1213.08	1.164E+06	1.00	24.2	0.0								
22	*		1299.18	1.262E+06	0.93	19.2	0.2				22	*		1299.18	1.268E+06	0.92	19.2	0.2								
23	*		1407.96	1.134E+06	1.03	6.2	-0.1				23	*		1407.96	1.125E+06	1.04	6.2	-0.2								
24	<!									24	<!															
25										25																
26										26																

Left file: C:\Users\121611\Desktop\011_046_specanal.LIS

Right file: C:\Users\121611\Desktop\011_045_specanal.LIS

Report type: All (Matching and Differences)

Summary: 14 Line(s) match | 10 Block(s) diff | 22 : 22 Line(s) diff

Ignore options: Ignore line terminators

Library 8 Comparison

C:\Users\121611\Desktop\011_046_specanal.LIS										C:\Users\121611\Desktop\011_045_specanal.LIS											
1	*	Analysis Code: specan1		Peaks: 36	Overall QA: 0.00	Library: 8	CD: 71	Timestamp: 2020-09-22 20:13:02.0		1	*	Analysis Code: specan1		Peaks: 36	Overall QA: 0.00	Library: 8	CD: 71	Timestamp: 2020-03-10 3:18:07.0	1		
2		Analysis ZT: 008.0000		AnalysisZY: 2020	Original: Yes	Comment:				2		Analysis ZT: 008.0000		AnalysisZY: 2020	Original: Yes	Comment:					
3										3											
4		Line	Nuclide	Isotope	Meta	DPM@ZT	DPM@CT	Atoms@ZT	Error	Accep	4		Line	Nuclide	Isotope	Meta	DPM@ZT	DPM@CT	Atoms@ZT	Error	Accep
5		Energy				Fraction	Error	StdDev			5		Energy				Fraction	Error	StdDev		
6											6										
7	*	1	CS	134		1.7070E+02	1.7010E+02	2.6740E+08	13.9		7	*	1	CS	134		1.7060E+02	1.7010E+02	2.6740E+08	13.9	
8											8										
9	*	802.51				1.00	13.90	0.00			9	*	802.54				1.00	13.90	0.00		
10	*	2	EU	152		2.1960E+02	2.1950E+02	2.2520E+09	27.6		10	*	2	EU	152		2.1960E+02	2.1950E+02	2.2550E+09	27.6	
11											11										
12	*	121.71				0.98	2.70	0.10			12	*	121.72				0.99	2.70	0.10		
13	*	244.68				0.56	9.00	1.50			13	*	244.69				0.56	9.00	1.50		
14	*	344.01				1.06	3.60	-0.30			14	*	344.03				1.05	3.60	-0.30		
15	*	443.58				1.06	18.10	-0.10			15	*	443.60				1.06	18.10	-0.10		
16	*	778.26				1.03	8.80	-0.10			16	*	778.29				1.03	8.80	-0.10		
17	*	866.84				0.97	23.80	0.10			17	*	866.87				0.97	23.80	0.10		
18	*	3	EU	154		4.9190E+02	4.9150E+02	3.2100E+09	33.5		18	*	3	EU	154		4.8510E+02	4.8470E+02	3.1630E+09	33.5	
19											19										
20	*	692.88				1.00	33.50	0.00			20	*	692.91				1.00	33.50	0.00		
21	*	4	PB	212		6.2920E+03	2.8230E+01	5.7950E+06	9.5		21	*	4	PB	212		6.3360E+03	2.8430E+01	5.8350E+06	9.5	
22											22										
23	*	238.35				1.00	9.50	0.00			23	*	238.36				1.00	9.50	0.00		
24	*	5	AC	228		4.7280E+05	4.1060E+01	2.5170E+08	18.2		24	*	5	AC	228		4.7280E+05	4.1060E+01	2.5170E+08	18.2	
25											25										
26	*	910.51				1.00	18.20	0.00			26	*	910.55				1.00	18.20	0.00		
27	*	6	TH	232		1.7860E+03	1.7860E+03	1.8970E+19	14.5		27	*	6	TH	232		2.0380E+03	2.0380E+03	2.1720E+19	14.5	
28											28										
29	*	63.37				0.10	14.50	0.00			29	*	63.37				0.12	14.50	0.00		
30	*	7	TH	234		1.1970E+03	1.0840E+03	5.9940E+07	14.5		30	*	7	TH	234		9.0740E+02	8.2150E+02	4.5430E+07	14.5	
31											31										
32	*	63.37				0.90	14.50	0.00			32	*	63.37				0.88	14.50	0.00		
33	*	92.64				1.00	7.90	0.00			33	*	92.64				1.00	7.90	0.00		
34	*	8	U	235		1.7670E+01	1.7670E+01	9.4370E+15	10.70	u	34	*	8	U	235		1.7610E+01	1.7610E+01	9.4040E+15	10.70	u
35											35										
36	*	185.50				1.00	10.70	0.00			36	*	185.51				1.00	10.70	0.00		

Left file: C:\Users\121611\Desktop\056_083_specanal.LIS
Right file: C:\Users\121611\Desktop\056_076_specanal.LIS

Show spaces / tabs

Library 9 Comparison

Report type: All (Matching and Differences)

Summary: 18 Line(s) match | 10 Block(s) diff | 145 : 142 Line(s) diff

Ignore options: Ignore line terminators

C:\Users\121611\Desktop\056_083_specanal.LIS										C:\Users\121611\Desktop\056_076_specanal.LIS											
1										1											
2	* Analysis Code: specanl Peaks: 102 Overall QA: 0.00 Library: 9 CD: 71 Timestamp:2020-11-18 21:05:02.0									2	Analysis Code: specanl Peaks: 102 Overall QA: 0.00 Library: 9 CD: 71 Timestamp:2020-10-26 03:23:31.0										
3	Analysis ZT: 052.7597		AnalysisZY: 2020		Original: Yes		Comment:			3	Analysis ZT: 052.7597		AnalysisZY: 2020		Original: Yes		Comment:				
4										4											
5	Line	Nuclide	Isotope	Meta	DPM@ZT	DPM@CT	Atoms@ZT	Error	Accept	Comment	5	Line	Nuclide	Isotope	Meta	DPM@ZT	DPM@CT	Atoms@ZT	Error	Accept	Comment
6	Energy	DPM@ZT			Fraction	Error	StdDev				6	Energy	DPM@ZT			Fraction	Error	StdDev			
7										7											
8	* 1	CO	60		8.2060E+06	8.1940E+06	3.2820E+13	11.70	u		8	1	CO	60		8.2020E+06	8.1900E+06	3.2800E+13	11.70	u	
9										9											
10	* 1174.01		8.2060E+06		1.00	11.70	0.00				10	1174.03		8.2020E+06		1.00	11.70	0.00			
11	* 2	ZR	95		1.3380E+08	1.2800E+08	1.7790E+13	2.60	u		11	2	ZR	95		1.3350E+08	1.2770E+08	1.7750E+13	2.60	u	
12										12											
13	* 723.62		2.1090E+08		0.63	2.60	0.00				13	723.64		2.1170E+08		0.63	2.60	0.00			
14	* 757.12		1.3380E+08		1.00	2.00	0.00				14	757.13		1.3350E+08		1.00	2.00	0.00			
15	* 3	ZR	97		1.0350E+10	1.7810E+08	1.5010E+13	22.30	u		15	3	ZR	97		1.0010E+10	1.7860E+08	1.4640E+13	22.70	u	
16										16											
17	* 355.78		9.8680E+09		1.05	18.60	-0.10				17	355.79		8.7620E+09		1.14	18.60	-0.30			
18	* 602.72		1.3750E+10		0.75	22.10	0.00				18	602.73		1.3160E+10		0.76	22.10	0.00			
19	* 743.77		1.0350E+10		1.00	1.30	0.00				19	743.78		1.0010E+10		1.00	1.30	0.00			
20	* 4	NB	95		8.3330E+06	7.6850E+06	6.0580E+11	8.80	u		20	4	NB	95		8.3360E+06	7.6860E+06	6.0560E+11	8.80	u	
21										21											
22	* 766.12		8.3330E+06		1.00	8.80	0.00				22	766.14		8.3360E+06		1.00	8.80	0.00			
23	* 5	MO	99		1.4960E+09	5.3310E+08	8.5390E+12	6.50	u		23	5	MO	99		1.4830E+09	5.2850E+08	8.4660E+12	6.50	u	
24										24											
25	* 181.11		1.3400E+09		1.12	2.20	-0.80				25	181.11		1.3350E+09		1.11	2.20	-0.70			
26	* 739.86		2.7890E+09		0.54	1.50	-0.10				26	739.87		2.8020E+09		0.53	1.50	-0.20			
27	* 778.29		3.0010E+09		0.50	2.90	2.90				27	778.30		2.9720E+09		0.50	2.90	2.90			
28	* 6	TC	99	M	7.2480E+10	8.7350E+05	3.7690E+13	0.60	u		28	6	TC	99	M	2.8440E+13	3.4270E+08	1.4790E+16	0.60	u	
29										29											
30	* 140.49		2.7860E+13		0.00	0.60	0.00				30	140.50		2.8440E+13		1.00	0.60	0.00			
31	* 7	RU	103		9.1310E+07	8.4940E+07	7.4450E+12	19.10	u		31	7	RU	103		9.3670E+07	8.7150E+07	7.6380E+12	19.00	u	
32										32											
33	* 497.28		9.1210E+07		1.00	1.60	0.00				33	497.29		9.3580E+07		1.00	1.60	0.00			
34	* 610.70		1.0830E+08		0.84	18.30	0.40				34	610.72		1.1070E+08		0.85	18.30	0.40			
35	* 8	RU	106		1.4350E+08	1.4240E+08	1.1090E+14	6.10	u		35	8	RU	106		1.4250E+08	1.4250E+08	9.5060E+18	6.10	u	
36										36											
37	* 621.31		1.4350E+08		1.00	6.10	0.00				37	621.33		1.4250E+08		1.00	6.10	0.00			
38	* 9	RH	105		7.7210E+08	1.1270E+08	2.3630E+12	3.20	u		38	9	RH	105		7.8260E+08	1.1420E+08	2.3950E+12	3.20	u	
39										39											
40	* 319.36		8.5110E+08		0.91	3.20	0.00				40	319.37		8.5560E+08		0.91	3.20	0.00			
41	* 10	SB	127		3.6190E+07	1.7330E+07	2.8950E+11	30.10	u		41	10	SB	127		3.3940E+07	1.6250E+07	2.7150E+11	32.30	u	
42										42											
43	* 473.49		5.0570E+07		0.72	21.90	0.60				43	473.50		5.2820E+07		0.64	21.90	0.80			
44	* 685.38		3.3220E+07		1.09	15.20	-0.20				44	685.39		3.4640E+07		0.98	15.20	-0.30			
45	* 11	TE	129		1.4890E+08	1.3690E+08	1.0400E+13	1.00	u		45	11	TE	129		2.8040E+08	2.5770E+08	1.9570E+13	1.00	u	
46										46											
47	* 487.24		1.4800E+10		0.01	1.00	0.00				47	487.24		1.4800E+10		0.02	1.00	0.00			
48	* 12	TE	132		1.1400E+09	4.7080E+08	7.5910E+12	0.60	u		48	12	TE	132		1.1380E+09	4.6970E+08	7.5740E+12	0.60	u	
49										49											
50	* 228.37		1.1400E+09		1.00	0.60	0.00				50	228.38		1.1380E+09		1.00	0.60	0.00			
51	* 13	TE	131	M	3.5530E+08	4.5880E+07	1.0220E+12	29.60	u		51	13	TE	131	M	4.4330E+08	4.5890E+07	1.1510E+12	31.00	u	
52										52											
53	* 149.50		5.9120E+08		0.60	17.70	0.90				53	149.50		7.0850E+08		0.63	17.70	0.90			
54	* 241.35		5.3660E+08		0.66	18.20	0.00				54	241.36		6.4640E+08		0.69	18.20	0.70			
55	* 334.48		3.3720E+08		1.05	23.50	-0.10				55	334.49		4.0610E+08		1.09	23.50	-0.20			
56	* 794.23		3.4460E+08		1.03	12.80	-0.10				56	794.25		4.1740E+08		1.06	12.80	-0.20			
57	* 852.68		3.3100E+08		1.07	10.50	-0.20				57	852.70		4.0630E+08		1.09	10.50	-0.30			
58	* 1126.10		3.7400E+08		0.95	14.10	0.10				58	1126.12		4.4970E+08		0.99	14.10	0.00			
59	* 14	I	124		2.5280E+06	1.2820E+06	2.1930E+10	22.10	u		59	14	I	124		2.5260E+06	1.2820E+06	2.1940E+10	22.10	u	
60										60											

61	*	602.72	1.0230E+07	0.25	22.10	0.00				61	602.73	1.0550E+07	0.24	22.10	0.00			
62	*	15	I	130	1.7380E+10	7.0700E+07	1.8600E+13	1.60	u	62	15	I	130	1.7800E+10	7.2390E+07	1.9040E+13	15.50	u
63	*									63								
64	*	668.06	2.0590E+11	0.08	0.70	0.20				64	668.08	2.0570E+11	0.09	0.70	0.30			
65	*	739.86	3.6360E+10	0.48	1.50	-0.10				65	685.39	1.3470E+11	0.13	15.20	-0.30			
66	*	16	I	131	4.2010E+08	2.9510E+08	7.0040E+12	5.80	u	66	739.87	3.6230E+10	0.49	1.50	-0.20			
67	*									67	16	I	131	4.2120E+08	2.9600E+08	7.0350E+12	5.80	u
68	*	284.54	5.4770E+08	0.77	3.90	1.20				68								
69	*	364.60	4.1700E+08	1.01	0.80	-0.10				69	284.55	5.5270E+08	0.76	3.90	1.20			
70	*	637.30	4.1110E+08	1.02	4.60	-0.10				70	364.61	4.1830E+08	1.01	0.80	-0.10			
71	*	17	I	132	1.7730E+09	7.3210E+08	1.1800E+13	4.20	u	71	637.31	4.0460E+08	1.04	4.60	-0.20			
72	*									72	17	I	132	1.7550E+09	7.2430E+08	1.1680E+13	4.20	u
73	*	506.12	1.6080E+09	1.10	3.50	-0.60				73								
74	*	522.91	1.5870E+09	1.12	1.30	-1.00				74	506.13	1.5820E+09	1.11	3.50	-0.60			
75	*	630.52	1.7890E+09	0.99	1.40	0.10				75	522.92	1.5770E+09	1.11	1.30	-1.00			
76	*	668.06	1.9740E+09	0.90	0.70	0.20				76	630.54	1.7250E+09	1.02	1.40	-0.10			
77	*	773.10	1.8170E+09	0.98	0.70	0.30				77	668.08	1.9740E+09	0.89	0.70	0.30			
78	*	812.61	2.0780E+09	0.85	2.70	0.90				78	773.11	1.8030E+09	0.97	0.70	0.30			
79	*	955.18	1.7400E+09	1.02	1.50	-0.20				79	812.63	2.0400E+09	0.86	2.70	0.90			
80	*	18	I	133	9.3810E+09	3.5780E+08	1.6910E+13	20.70	u	80	955.20	1.6910E+09	1.04	1.50	-0.30			
81	*									81	18	I	133	9.5330E+09	3.6190E+08	1.7160E+13	20.70	u
82	*	530.18	9.3590E+09	1.00	0.90	0.00				82								
83	*	707.06	1.0650E+10	0.88	16.70	0.30				83	530.19	9.5110E+09	1.00	0.90	0.00			
84	*	856.71	8.1150E+09	1.16	18.80	-0.40				84	707.08	1.0850E+10	0.88	16.70	0.30			
85	*	1298.88	1.4610E+10	0.64	8.40	1.20				85	856.73	8.2190E+09	1.16	18.80	-0.40			
86	*	19	XE	133	5.7570E+08	3.3540E+08	6.2760E+12	0.90	u	86	1298.91	1.4800E+10	0.64	8.40	1.20			
87	*									87	19	XE	133	5.5890E+08	3.2570E+08	6.0970E+12	0.90	u
88	*	80.92	5.7770E+08	1.00	0.90	0.00				88								
89	*	20	XE	135	4.1600E+10	2.4310E+07	3.2910E+13	2.80	u	89	80.92	5.6080E+08	1.00	0.90	0.00			
90	*									90	20	XE	135	4.2840E+10	2.4310E+07	3.3760E+13	2.80	u
91	*	249.79	4.1600E+10	1.00	2.80	0.00				91								
92	*	21	CS	136	1.2450E+07	1.0010E+07	3.3650E+11	12.20	u	92	249.80	4.2840E+10	1.00	2.80	0.00			
93	*									93	21	CS	136	1.1960E+07	9.6440E+06	3.2700E+11	12.20	u
94	*	340.24	1.2490E+07	1.00	12.20	0.00				94								
95	*	22	BA	140	5.2020E+08	4.1650E+08	1.3780E+13	16.10	u	95	340.25	1.1990E+07	1.00	12.20	0.00			
96	*									96	22	BA	140	5.2040E+08	4.1670E+08	1.3790E+13	16.10	u
97	*	162.70	2.8640E+08	1.82	4.10	-4.00				97								
98	*	305.11	6.5120E+08	0.80	4.80	0.90				98	162.71	2.8690E+08	1.81	4.10	-4.00			
99	*	423.85	6.3800E+08	0.82	5.50	0.80				99	305.12	6.4970E+08	0.80	4.80	0.90			
100	*	437.67	5.8470E+08	0.89	9.80	0.40				100	423.86	6.3400E+08	0.82	5.50	0.80			
101	*	537.54	5.9530E+08	0.87	1.20	1.10				101	437.67	5.8440E+08	0.89	9.80	0.30			
102	*	23	LA	140	2.2740E+09	4.2030E+08	7.9330E+12	8.90	u	102	537.55	5.9530E+08	0.87	1.20	1.10			
103	*									103	23	LA	140	2.2370E+09	4.1290E+08	7.7970E+12	8.90	u
104	*	328.88	2.0590E+09	1.10	1.30	-0.90				104								
105	*	432.62	3.3340E+09	0.68	5.60	1.30				105	328.88	2.0220E+09	1.11	1.30	-0.90			
106	*	487.24	2.2970E+09	0.99	1.00	0.00				106	432.63	3.2370E+09	0.69	5.60	1.30			
107	*	751.95	2.4800E+09	0.92	5.30	0.40				107	487.24	2.2800E+09	0.98	1.00	0.00			
108	*	816.23	2.7330E+09	0.83	1.70	1.30				108	751.96	2.5000E+09	0.89	5.30	0.50			
109	*	868.30	2.5610E+09	0.89	4.20	0.50				109	816.24	2.6990E+09	0.83	1.70	1.30			
110	*	920.14	2.5070E+09	0.91	6.90	0.40				110	868.31	2.5230E+09	0.89	4.20	0.60			
111	*	925.76	2.6790E+09	0.85	3.80	0.80				111	920.16	2.4910E+09	0.90	6.90	0.40			
112	*	24	CE	141	7.8260E+07	7.1730E+07	5.2860E+12	1.70	u	112	925.78	2.6330E+09	0.85	3.80	0.80			
113	*									113	24	CE	141	7.8260E+07	7.1730E+07	5.2840E+12	1.70	u
114	*	145.44	7.8260E+07	1.00	1.70	0.00				114								
115	*	25	CE	143	4.0280E+09	5.1400E+08	1.1520E+13	6.90	u	115	145.45	7.8260E+07	1.00	1.70	0.00			
116	*									116	25	CE	143	4.1100E+09	5.2290E+08	1.1740E+13	7.10	u
117	*	293.35	4.0090E+09	1.00	0.70	-0.10				117								
118	*	350.77	4.2260E+09	0.95	4.10	0.20				118	293.36	4.0980E+09	1.00	0.70	0.00			
119	*	490.46	5.8110E+09	0.69	5.60	1.30				119	350.77	4.0270E+09	1.02	4.10	-0.10			
120	*	26	ND	147	1.3390E+08	1.0340E+08	3.0540E+12	3.20	u	120	490.47	6.3910E+09	0.64	5.60	1.50			
121	*									121	26	ND	147	1.3440E+08	1.0380E+08	3.0650E+12	3.20	u
122	*	91.09	1.3390E+08	1.00	2.90	0.00				122								
123	*	319.36	1.4420E+09	0.09	3.20	0.00				123	91.09	1.3440E+08	1.00	2.90	0.00			
124	*	27	PM	151	6.7750E+05	6.1680E+04	1.6650E+09	12.20	u	124	319.37	1.5750E+09	0.09	3.20	0.00			

125 *										125	27	PM	151	6.1150E+05	5.5670E+04	1.5030E+09	12.20	u	
126 *	340.24	2.2940E+08	0.00	12.20	0.00					126									
127 *	28	SM	153	6.8900E+07	1.5850E+07	2.7610E+11	15.20	u		127	340.25		2.3460E+08	0.00	12.20	0.00			
128 *										128	28	SM	153	7.0300E+07	1.6380E+07	2.8420E+11	15.20	u	
129 *	102.82	6.8900E+07	1.00	15.20	0.00					129									
130 *	29	EU	156	7.2870E+08	6.0470E+08	2.3000E+13	2.60	u		130	102.82		7.0300E+07	1.00	15.20	0.00			
131 *										131	29	EU	156	6.7280E+08	5.5830E+08	2.1230E+13	2.60	u	
132 *	723.62	1.9930E+09	0.37	2.60	0.00					132									
133 *	30	TB	161	4.1240E+08	2.7330E+08	5.9030E+12	24.30	u		133	723.64		1.8210E+09	0.37	2.60	0.00			
134 *										134	30	TB	161	4.2830E+08	2.8420E+08	6.1490E+12	24.70	u	
135 *	57.37	1.1770E+09	0.35	15.10	1.70					135									
136 *	74.92	4.0170E+08	1.03	5.20	-0.10					136	57.37		1.3150E+09	0.33	15.10	1.70			
137 *	31	AC	228	6.6920E+10	1.0290E+06	3.5560E+13	9.40	u		137	74.92		4.1760E+08	1.03	5.20	-0.10			
138 *										138	31	AC	228	5.8030E+02	8.7670E-03	3.0790E+05	9.40	u	
139 <!	910.74	2.7850E+12	0.02	9.40	0.00					139									
140 <!	32	TH	228	4.1960E+07	4.1790E+07	6.0860E+13	18.20	u		140	910.76		2.5210E+12	0.00	9.40	0.00			
141 <!										141	32	TH	228	3.8250E+07	3.8090E+07	5.5520E+13	9.40	u	
142 <!	241.35	1.2420E+08	0.34	18.20	0.00					142									
143 <!	910.74	4.2990E+07	0.98	9.40	0.00					143	910.76		3.8250E+07	1.00	9.40	0.00			
144 <!	33	TH	232	1.4190E+12	1.4190E+12	1.5080E+28	0.60	u		144	33	TH	234	1.6670E+04	1.4820E+04	8.3470E+08	8.80	u	
145										145									
146 *	140.49	1.4230E+12	1.00	0.60	0.00					146	766.14		4.1680E+09	0.00	8.80	0.00			
147 *	34	TH	234	6.4390E+03	5.7250E+03	3.2240E+08	8.80	u		147	34	PA	231	2.2710E+07	2.1140E+07	5.6480E+17	0.90	u	
148										148									
	!>									149	80.92		1.4830E+10	0.00	0.90	0.00			
	!>									150	35	U	234	1.8750E+04	1.8750E+04	3.4920E+15	8.80	u	
	!>									151									
149 *	766.12	1.0780E+10	0.00	8.80	0.00					151									
150 *	35	PA	231	2.4100E+07	2.2430E+07	5.9930E+17	0.90	u		152	766.14		3.7050E+09	0.00	8.80	0.00			
151										153	36	U	235	2.2710E+07	2.2710E+07	1.2130E+22	0.90	u	
152	80.92	1.4830E+10	0.00	0.90	0.00					154									
153 *	36	U	234	7.2430E+03	7.2430E+03	1.3490E+15	8.80	u		155	80.92		1.3800E+10	0.00	0.90	0.00			
154										156	37	U	235	2.2710E+07	2.2710E+07	1.2130E+22	0.90	u	
155 *	766.12	9.5880E+09	0.00	8.80	0.00					157									
156 *	37	U	235	2.4100E+07	2.4100E+07	1.2870E+22	0.90	u		158	80.92		1.3800E+10	0.00	0.90	0.00			
157										159	38	U	238	2.0810E+03	2.0810E+03	7.0550E+18	8.80	u	
158	80.92	1.3800E+10	0.00	0.90	0.00					160									
159 *	38	U	238	8.0380E+02	8.0380E+02	2.7250E+18	8.80	u		161	766.12		9.5880E+09	0.00	8.80	0.00			
160										162									
161 *	766.12	9.5880E+09	0.00	8.80	0.00					163									
162																			
163																			