

## LA-UR-21-21117

Approved for public release; distribution is unlimited.

Title: Medium Resolution Gamma Spectrometers for Isotopic Analysis of Uranium and Plutonium at Los Alamos National Laboratory

Author(s): Vo, Duc Ta

Intended for: International Workshop on Isotopic Analysis of Uranium and Plutonium by Nondestructive Assay Techniques for Nuclear Safeguards

Issued: 2021-02-08

---

**Disclaimer:**


Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

# Medium Resolution Gamma Spectrometers for Isotopic Analysis of Uranium and Plutonium at Los Alamos National Laboratory

Duc Vo

# FRAM v.6.1

- FRAM v.6.1 (released March 2020) can analyze the MRGS spectra.
- FRAM uses a nonlinear least squares fit technique, combining the Powell's minimization method with the linear least squares fit to fit the peaks of the LaBr3 and CZT spectra.



# Uranium ANOVA analysis


Parameter	Global biases	<sup>234</sup> U	<sup>235</sup> U	<sup>238</sup> U
GePInr_ULEU_060-250 GePInr_UHEU_060-250	Estimate of $\delta_R$	9.0	2.1	8.4
	Estimate of $\delta_S$	1.4	0.9	5.8
	STDEV of Relative Bias	3.8	0.9	6.7
GeCoax_ULEU_120-1010 GeCoax_UHEU_120-1010	Estimate of $\delta_R$	16.0	2.8	2.6
	Estimate of $\delta_S$	1.2	0.7	0.9
	STDEV of Relative Bias	4.3	1.1	1.2
CZT500_U_120-1010	Estimate of $\delta_R$	33.3	32.4	13.6
	Estimate of $\delta_S$	11.3	6.7	23.6
	STDEV of Relative Bias	13.4	10.1	23.6
LaBr_U_120-1010	Estimate of $\delta_R$	9.1	8.8	10.1
	Estimate of $\delta_S$	18.8	9.1	50.9
	STDEV of Relative Bias	18.4	9.3	49.1

- The random uncertainties  $\delta_R$ , item-specific systematic uncertainties  $\delta_S$ , and the STDEV of the relative biases of the MRGS data are about a factor of ten worse than those of HPGe



# Plutonium ANOVA analysis

Parameter	Global biases	<sup>238</sup> Pu	<sup>239</sup> Pu	<sup>240</sup> Pu	<sup>241</sup> Pu	<sup>241</sup> Am
GePInr_Pu_060-230	Estimate of $\delta_R$	6.0	0.2	0.8	0.5	0.7
	Estimate of $\delta_S$	2.0	0.1	0.7	0.4	1.2
	STDEV of Relative Bias	2.4	0.1	0.8	0.4	1.2
GePInr_Pu_120-420	Estimate of $\delta_R$	6.0	0.4	2.1	0.8	1.4
	Estimate of $\delta_S$	1.9	0.2	0.8	0.4	1.6
	STDEV of Relative Bias	2.4	0.1	0.8	0.4	1.2
GeCoax_Pu_120-420	Estimate of $\delta_R$	6.1	0.3	2.1	0.4	0.8
	Estimate of $\delta_S$	0.8*	0.1	0.7	0.4	1.0
	STDEV of Relative Bias	1.8	0.1	0.9	0.4	0.9
GeCoax_Pu_180-1010	Estimate of $\delta_R$	10.9	0.6	3.3	0.7	1.0
	Estimate of $\delta_S$	7.7	0.4	1.7	0.8	0.9
	STDEV of Relative Bias	7.9	0.4	1.8	0.8	0.9
CZT500_Pu_120-500	Estimate of $\delta_R$	59.8	2.0	26.2	20.8	48.2
	Estimate of $\delta_S$	35.1	0.9	10.2	13.1	17.2
	STDEV of Relative Bias	37.9	0.9	11.0	13.5	25.8
LaBr_Pu_200-750	Estimate of $\delta_R$	83.7	3.7	34.0	11.4	5.3
	Estimate of $\delta_S$	95.2	1.5	12.7	38.1	3.4
	STDEV of Relative Bias	106.0	1.9	16.2	37.9	3.9

# MARIA (MARIA Analyses Ratios of Isotope Activities)



Regions for: LABR\_Pu\_A\_16.CHN  
Reduced Chi-square: 90651.63



- MARIA analyzes LaBr3 spectra of plutonium and MOX.
- It works better with MOX than plutonium.


# Differences between FRAM and MARIA

---

- FRAM steps:
  - Fit peaks in different regions one region at a time
  - Determine the relative efficiency curve
  - Calculate the relative activities of the isotopes
  - Adjust the background of the analysis regions
  - Repeat the above steps several times
  
- MARIA fits all the peaks in all the regions, determines the efficiency curve, and calculates the relative activities of the isotopes in one single step.
  
- MARIA has the potential to be an excellent tool for isotopic analysis of plutonium and MOX spectra acquired by LaBr3 detectors.



# H3D CZT peak fitting (high energy tail)



- Left: 186-keV peak; right: 1001-keV peak.
- Top: FRAM v.6.1; bottom: updated FRAM