

LA-UR-21-21931

Approved for public release; distribution is unlimited.

Title: Exploring Neutron Imaging for Non-DE Analysis of corrosion in ICCWR of 3013 Containers

Author(s): Kruk, Izabela Iwona

Intended for: Corrosion Working Group Meeting, 2021-02-23/2021-02-24 (Los Alamos, New Mexico, United States)

Issued: 2021-02-25

Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

Exploring Neutron Imaging for Non-DE Analysis of corrosion in ICCWR of 3013 Containers

Acknowledgments:

MIS team (AMPP4): Juan G. Duque, Daniel Rios, Joshua E. Narlesky, Ed Romero,
Mary Ann Stroud, Douglas Veirs, Elizabeth Kelly, Kim Kaufeld, Laura Worl, John Berg (C-PCS)
Erik B. Watkins (MPA-11), Markus Strobl and Jacopo Valescchi (PSI), Switzerland, El'ad Caspi (NRCN), Israel

Packaging style

DOE stores excess weapons plutonium in nested, welded stainless steel containers (3013 Containers)

3013 containers: 3 layers

1. Convenience container (different features depending on origin, some have filters for outgassing, ie: Hanford-Livermore)-CC. In contact with material.
2. Inner container (welded)-IC. In contact with gases generated by the material.
3. Outer container (welded)-OC.

Expected lifetime of packaged 3013 containers is 50 years.



Corrosion Observations at DE

DE of some 3013 container packaged with chloride-bearing PuO_2 and moisture have shown corrosion on the convenience and inner containers.

Small scale corrosion experiments suggest chlorine-containing gases (HCl and Cl_2) may have caused the observed headspace corrosion in 3013 containers.


CaCl_2 supports release of Cl_2 and HCl

MgCl_2 supports release of HCl only

Convenience Can



Inner Can Lid



Most Common Corrosion Analysis Techniques used in 3013 containers

Laser confocal microscope (LCM)


Identification of pits, cracks and surface corrosion.

X-ray tomography


Identification of pits and cracks.

Metallographic cross-sections

Identification of pit and crack depths



Cross-section



All great techniques for corrosion analysis but they require small sections of the 3013 container which means the container has to be DE.

Can we perform this level of corrosion analysis without sectioning the 3013 container?

Why is it important to maintain the 3013 container intact?

Being able to do corrosion analysis in the ICCWR of intact 3013 containers allow us to study the progression of corrosion in these container over time.

Benefits:

1. Monitoring corrosion over extended period of time without interrupting the experiment
2. Observing changes in corrosion rate for better understanding the reaction kinetics
3. Measuring crack growth rate over time
4. Inspection of pits and cracks at the moment of formation

Why neutrons?



1. Higher resolution
2. Deeper penetration of material
3. Better contrast among elements with similar Z
4. Neutrons activate materials, so built-in infrastructure to mitigate radioactive hazards
5. Neutrons are isotope selective

Attenuation coefficients for thermal neutrons [cm⁻¹]



1a	2a	3b	4b	5b	6b	7b	8	1b	2b	3a	4a	5a	6a	7a	0
H															He
3.44															0.02
Li	Be									B	C	N	O	F	Ne
3.30	0.79									101.60	0.56	0.43	0.17	0.20	0.10
Na	Mg									Al	Si	P	S	Cl	Ar
0.09	0.15									0.10	0.11	0.12	0.06	1.33	0.03
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se
0.06	0.08	2.00	0.60	0.72	0.54	1.21	1.19	3.92	2.05	1.07	0.35	0.49	0.47	0.67	0.73
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te
0.08	0.14	0.27	0.29	0.40	0.52	1.76	0.58	10.88	0.78	4.04	115.11	7.58	0.21	0.30	0.25
Cs	Ba	La	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po
0.29	0.07	0.52	4.99	1.49	1.47	6.85	2.24	30.46	1.46	6.23	16.21	0.47	0.38	0.27	
Fr	Ra	Ac	Rf	Ha											
0.34															
	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu	
	0.14	0.41	1.87	5.72	171.47	94.58	1479.04	0.93	32.42	2.25	5.48	3.53	1.40	2.75	
	Th	Pa	U	Np	Pu	Am	Cm	Bk	Cf	Es	Fm	Md	No	Lr	
	0.59	8.46	0.82	9.80	50.20	2.86									

Attenuation coefficients for X-ray [cm⁻¹] (150kV)


1a	2a	3b	4b	5b	6b	7b	8	1b	2b	3a	4a	5a	6a	7a	0
H															He
0.02															0.02
Li	Be									B	C	N	O	F	Ne
0.06	0.22									0.28	0.27	0.11	0.16	0.14	0.17
Na	Mg									Al	Si	P	S	Cl	Ar
0.13	0.24									0.38	0.33	0.25	0.30	0.23	0.20
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se
0.14	0.26	0.48	0.73	1.04	1.29	1.32	1.57	1.78	1.96	1.97	1.64	1.42	1.33	1.50	1.23
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te
0.47	0.86	1.61	2.47	3.43	4.29	5.06	5.71	6.08	6.13	5.67	4.84	4.31	3.98	4.28	4.06
Cs	Ba	La	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po
1.42	2.73	5.04	19.70	25.47	30.49	34.47	37.92	39.01	38.61	35.94	25.88	23.23	22.81	20.28	20.22
Fr	Ra	Ac	Rf	Ha											
	11.80	24.47													
	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu	
	5.79	6.23	6.46	7.33	7.68	5.66	8.69	9.46	10.17	10.91	11.70	12.49	9.32	14.07	
	Th	Pa	U	Np	Pu	Am	Cm	Bk	Vf	Es	Fm	Md	No	Lr	
	28.95	39.65	49.08												



Cold neutron phase contrast imaging (PCI) at LANSCE: in-line propagation method (resolution $>50\text{ }\mu\text{m}$)



Needs high spatial coherence, results in low flux and long count times ($\sim 24\text{ hrs}$)




PCI of Pb fishing weights:


Due to the low absorption contrast, a radiographs contain little information (center).

Phase contrast imaging ($\sim 50\mu\text{m}$ resolution) was used to enhance the visibility of interfaces and interior defects in the Pb (right).


More capabilities of neutron PCI



Photograph of
cracked 50x50 mm Al plate
Crack width: 150 μm







Phase contrast
image using the
propagation method,
10 minutes exposure




photos of the cracked plate between
20 cm of Al shielding


Article

Neutron Imaging at LANSCE—From Cold to Ultrafast

Ronald O. Nelson ^{1,*}, Sven C. Vogel ² , James E. Hunter ³, Erik B. Watkins ⁴, Adrian S. Losko ², Anton S. Tremsin ⁵ , Nicholas P. Borges ^{2,6}, Theresa E. Cutler ⁷, Lee T. Dickman ⁸ , Michelle A. Espy ³, Donald Cort Gautier ³, Amanda C. Madden ⁹, Jaroslaw Majewski ^{10,11}, Michael W. Malone ¹², Douglas R. Mayo ¹³, Kenneth J. McClellan ², David S. Montgomery ¹⁴, Shea M. Mosby ¹, Andrew T. Nelson ¹⁵ , Kyle J. Ramos ¹⁶, Richard C. Schirato ⁹, Katlin Schroeder ¹⁷, Sanna A. Sevanto ⁸, Alicia L. Swift ¹⁸, Long K. Vo ^{2,19}, Thomas E. Williamson ²⁰ and Nicola M. Winch ⁷




PCI of the plate
inside 8 cm of Al
shielding, 20 minutes
exposure




PCI of the plate inside
20 cm of Al shielding,
25 minutes exposure

More capabilities of neutron imaging: Isotope separation




0.2% ^{235}U U_3Si_5



8.8% ^{235}U U_3Si_5

Examples of transmission data analysis for two fuel samples with different enrichment

Sample data collected with a liquid scintillator detector.






UN/ U_3Si_5 and U_3Si_5 fuel pellets with different enrichment levels

U-235 0.34 g/cc
in U-238 12.43 g/cc

U-235 0.61g/cc
in U-238 6.91 g/cc



U-235 0.02g/cc
in U-238 7.56 g/cc



tomography




Neutron computed tomography

Can we perform this level of corrosion analysis without sectioning the 3013 container?



- Exposure to corrosive conditions for 109 days

- 4.5 mL of 5.4 M HCl
- 2.39 g of TCCA




In-line propagation neutron imaging of 3013 container fragment


Transmission based neutron radiography does not detect any inhomogeneity in the fragment. Using in-line propagation imaging, ~100-500 micron features appear at the top of the sample. Increasing propagation distance increases contrast but decreases resolution.




250mm propagation



400mm propagation



650mm propagation



The H6284 piece of IC was imaged at following conditions:

Neutron wavelengths:
4-13 Angstroms

Source: 2mm diameter
pinhole

Source-to-detector
distance: 5.75 m

Camera: ATIK 490 EX CCD


Scintillator: 10x10cm, 35
micron resolution,

ZnS:Cu/6LiF (RCTritec)

Proof of Concept experiments

Obtain neutron PC Images of ICCWR samples and Compare with LCM results:

- 1/8 ICCWR exposed to HCl gas and 1/8 ICCWR exposed to Cl_2 gas
- ICCWR of pristine full scale container to establish baseline and neutrons transmission



Examination of root growth of lupin over 4 weeks as an example of an in situ, continuous experiments


Achievements: whole piece of the sample imaged; corrosion features identified also sub-surface

Experiments *in situ*

Based on proof of concept results, neutron images of ICCWR could be collected on a single 3013 inner container being exposed to Cl_2 gas at 50% RH:

- Baseline measurement on a pristine 3013 inner container
- ICCWR area of an intact IC exposed to HCl gas
- ICCWR area of an intact IC exposed to Cl_2 gas

Further plans: Checking whether the corrosion effects continue developing under conditions as observed in DE 3013 (lower RH)



Concluding remarks

1. Goal for Non-Destructive Analysis: Ability to study crack growth over time and evaluate impact of various parameters on crack growth
2. Potential of Neutron Imaging for corrosion analysis of 3013 containers
3. Neutron source available
4. No reports of in situ corrosion monitoring inside objects
5. DOE and NNSA mission relevance: defining safety envelope of storage of the plutonium oxide inventory that contains chloride salts
6. Challenges: new technique under development, relying on the accelerator cycle (6 months not operating)