

## LA-UR-21-22328

Approved for public release; distribution is unlimited.

Title: Construction, characterization and crystal structure of a fluorescent single-chain Fv chimera

Author(s): Velappan, Nileena

Intended for: Webinar organized by Journal and Publisher for an Editor's choice article

Issued: 2021-03-09

---

**Disclaimer:**

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.


# Construction, characterization and crystal structure of a fluorescent single-chain Fv chimera

Nileena Velappan, Biosciences Division

March 24, 2021

# Fluorescent Antibodies: Workhorse of Biomedical Research



- Antibodies (IgG) are large Y shaped proteins of our immune system that identify and neutralize pathogenic bacteria and viruses
- Monoclonal antibodies are effective therapeutic agents as well as a highly valuable and ubiquitous biomedical research reagent
- Fluorescently labeled antibodies are extensively used in diagnostic and research methods such as immunofluorescent microscopy and flow cytometry.



R. Rouet et al. / FEBS Letters 588 (2014) 269–277



N. VELAPPAN ET AL MABS 2020 e1843754-2.



# Generating Fluorescent Antibodies: Current techniques and Concerns

- Fluorescent IgG are typically produced by chemical reaction, e.g. a succinimidyl ester functional group attached to a fluorophore core reacts with primary amines to label the antibody.



## Potential Problems

- Labeling of primary amines in the binding region may block antigen recognition
- Different antibodies react with different fluorophores at different rates, creating variability in experimental settings
- Labeling reaction is stochastic, some antibodies will have large numbers of fluorophores and some very few
- Over conjugation can result in quenching
- Batch-to-batch variations may result in variations in antigen recognition signal




# *In vitro* Generation of Highly Specific Antibodies

- *In vitro* antibody selection utilizing large combinatorial libraries displayed on phage and/or yeast can generate highly specific antibodies with high affinity
- These display technologies often use scFv -V<sub>H</sub> and V<sub>L</sub> regions with an unstructured linker peptide
- However, scFvs tend to
  - aggregate
  - low expression levels
  - unstable during long-term storage
  - Binding signal detection requires a labeled secondary antibody that recognizes an expression tag





# Strategies to improve scFv functionality

- scFv-Fc improved stability and allow recognition by secondary antibodies
- scFv-Enzyme (e.g. alkaline phosphatase)
- Single chain-fluorescent protein (scFP) chimera
  - FP expression on C-terminus
  - FP protein as linker between VH and VL regions of the scFv



# Thermal Green Protein (TGP) based scFPs


- Extremely stable
- Highly soluble
- Non-aggregating green fluorescent protein
- TGP is monomeric
- TGP was engineered by disrupting crystal lattice contacts and introducing high-entropy glutamate residues to improve crystallization and prevent oligomerization.
- Potentially suitable for insertion into scFv genes, substituting the linker originally used to tether the VL and VH.





# Construction of scTGPs


- Inverse PCR amplified scFv gene and vector sequence
- TGP was amplified with overlapping primers
- Circular polymerase extension cloning (CPEC) assembly was used to create scTGP constructs
  - Yeast display vectors
  - Protein expression vectors



# Evaluating scTGP functionality using yeast display

scFvs p1-1 and p13B1 recognizes phospho-tyrosines on the Fc $\epsilon$ R1 receptor

Antigen recognition



Display level (PE or Green fluorescence)




- Yeast can display folded scTGP
- Display level and antigen recognition signal for scTGPs is lower than scFv
- ScTGPs recognize specific antigen with similar affinity (~3x range)

	scFv p1-1	scTGP p1-1	scFv p13B1	scTGP p13B1
<b>Display level (MFI) (SV5-PE)</b>	9111	1280	10,215	5168
<b>Recognition signal for specific target (MFI) (biotinylated target-streptavidin A633)</b>	4696	283	4882	893
<b>Recognition signal for non-specific target M2 (MFI) (biotinylated target-streptavidin A633)</b>	160	127	160	118
<b>K<sub>D</sub> values for specific antigen (nM)</b>	2.43±0.22	8.55±2.33	1.40±0.65	0.51±0.56




# Determining protein expression system for scTGPs

- scFvs require periplasmic expression
- Fluorescent proteins express in cytoplasm, which usually produces higher quantities of protein
- Evaluated scTGP expression in both periplasm and cytoplasm and compared functionality using FLISA
- Cytoplasmic expression (pET CK3) produced statistically higher amounts of scTGPs and they also recognized the specific antigen in FLISA



# scFv- scTGP Protein Expression Comparison

- Three scFv and scTGP pairs were compared for protein yield
- scFvs with good expression levels (p1-1 and p13B1), had comparable expression and yield as scTGPs
- For  $\alpha$ TyrS antibody a problem scFv, conversion to scTGP increased purified protein yield by 10x
- Monitoring protein expression and purification is considerably easier for scFPs





Antibody name	scFv (mg)	scTGP (mg)
p1-1	4.3	6.9
p13B1	5.5	6.9
$\alpha$ TyrS	0.7	6.0




# scFv – scTGP Protein Functionality Comparison

- Equal concentration of scFv and scTGP used to access antigen recognition of three antibodies
- ELISA data showed both antibody formats have equivalent functionality for antigen recognition
- FLISA data shows that antigen recognition by scTGPs can be measured using fluorescence in a single step assay format
- Additional experiments also showed that at equimolar concentrations both antibody format give equivalent binding signal



# Crystal Structure of scTGP p1-1

- Bipyramid shaped crystals of scTGP p1-1 formed with 1.6 M citric acid pH 6.2 as crystallization buffer after 2 months.
- Diffracted to better than 2.5 Å resolution
- The final  $R_{\text{work}}$  and  $R_{\text{free}}$  values were 16.3% and 21.1%
- The refined structure contains 1 scTGP p1-1, 2 glycerols, and 136 water molecules.




PDB ID:6WZN



# Epitope binning of $\alpha$ *Yersinia pestis* antibodies

- The effect of unlabeled IgGs on scTGP's binding was used to determine whether the two antibody formats bind to same or different epitopes, in a one-step assay (A)
- An scTGP-F1V binding curve (B) was used to estimate the half-saturating concentration of scTGP, optimal for detection of scTGP displacement
- Binding of scTGP was detected by fluorescence in the presence of non-fluorescent IgGs
- Panel C shows that IgG YP3 and YP8 does not displace scTGP 2 (no drop in signal) and IgG YP2 displaces scTGP2 as expected.
- Panel D shows antibodies YP3 and YP8 compete for the same epitope



# Conclusions

- The scTGPs are recombinant intrinsically fluorescent antibodies, ready for use in traditional fluorescent antibody based assays
- The scFPs can be constructed as vL-FP-vH format, TGP an extremely stable monomeric fluorescent protein is suitable FP
- The scFPs can be effectively displayed on yeast, where they function similarly to their corresponding scFv counterparts, recognize specific antigen similar affinity
- Functional expression levels of scFPs in the bacterial cytoplasm are relatively high and can rescue some scFvs otherwise expressed at low levels in the periplasm
- Straightforward assessment of expression levels and monitoring of protein purification steps. Antigen recognition in one step binding assays, novel assays such as epitope binning
- TGP is a monomeric fluorescent protein amenable to crystallization. These features provide unique and valuable characteristics to scTGP molecules for their use in protein chemistry
- We present unique antibody format with wide variety of use in research and clinical setting





# Acknowledgments

- **Andrew R.M. Bradbury**
- **Devin Close**
- **Li-Wei Hung**
- **Leslie Naranjo**
- **Colin Hemez**
- **Natasha DeVore**
- **Donna K. McCullough**
- **Antonietta M. Lillo**
- **Geoffrey S. Waldo**
- National Institutes of Health
- Foundation for the National Institutes of Health
- Los Alamos National Laboratory's Laboratory Directed Research & Development
- DOE Office of Science through the National Virtual Biotechnology Laboratory, a consortium of DOE national laboratories focused on response to COVID-19

