

LA-UR-21-23589

Approved for public release; distribution is unlimited.

Title: Agnostic Immunity

Author(s): Mukundan, Harshini

Intended for: CNES CAPABILITY REVIEW

Issued: 2021-04-14

Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.



Agnostic Immunity


Harshini Mukundan

May 2021


Replace and add LA-UR number




Managed by Triad National Security, LLC., for the U.S. Department of Energy's NNSA.




There is a need for universal, agnostic, layered strategies to combat unknown unknowns – especially at the point of need





Innate Immunity




Elegant pattern
recognition network



Early (pre-symptomatic)
recognition system

Universal

Rapid




Bacterial Pathogen Associated Molecular Patterns (PAMPs)




Association of Bacterial Amphiphiles with Host Carriers


Amphiphilic PAMPS universally bind to host lipoprotein carriers (HDL, LDL), and occur associated with them in blood



Copyright © 2005 by Elsevier Inc.







Universal *Bacterial* Biosensor






Clinical Validation

Lipopolysaccharides (Gram-negative) and Lipoteichoic acids (Gram-positive) bacterial pathogens






Clinical Validation

Lipoarabinomannans
(Gram-indeterminate)



Engineering Lipoprotein Nanodiscs

Reproducible, Controlled Lipoprotein Nanodiscs can be potential ligands for lipids, concentration matrices or early therapeutics




MSP	Lipid	Size (nm)	Zeta (mV)*
1D1	100% DMPC	9.9 ± 2.5	-3.9
	20% DMTAP	10.9 ± 2.8	4.2
	30% DMTAP	8.6 ± 1.5 Multiple peaks	7.9
	40% DMTAP	Multiple large (>100) peaks	N/A
	50% DMTAP	1.591 ± 0.33	N/A

Mycolactone, Small Molecule Toxin associated with Buruli Ulcer




MD simulations and all-atomistic models to generate ND


Lipopolysaccharides from *F. tularensis*




Lipopolysaccharides



Host
Lipoproteins



LPS binds to host lipoproteins





Presentation of LPS impacts assay performance and antibody selection




Select agent PAMPs associate with host lipoprotein, and this presentation impacts assay performance

Presentation Matters



Condition 1: LPS in buffer
Condition 2: LPS in serum

Lipopolysaccharides




Lipoarabinomannans

Association with carriers completely changes cytokine profiles – specifically impacting pro-inflammatory cytokines and chemokines


Intelligent Immunity

Laboratory Data




20,000 data points/biomarker combination every 2 weeks

Model Development



Model Development


Clinical Validation





3 time pts/patient, 2200 patients, 18 years,
1600 parameters/collection
= 18,625,200 data points

Separation of Proteins and lipids from Blood – The Traditional Way

Step 1. Separation of serum from blood





Step 2: Separation of amphiphilic biomarkers



Catch-All, Integrated Sample Processing at the point of need



- Two-tiered novel cross-flow filtration system for separation of aqueous and organic components
- Membrane interfaces compatible with lipid extraction – use of chloroform and methanol in a microfluidics cassette



First microfluidic system for separation of lipids and proteins




Performance of the Chip is Comparable to the Laboratory Method




Comparison of Cell Counts (Cells/mL)

Comparison of Lipoarabinomannan Extraction (nM)

Field Ready Waveguide-based Optical Biosensor Platform





Bench-top Instrument



Field-Ready Instrument




App-based read-out on phone



Comparison of Performance of the Old Vs. New Sensor

- Enhanced dynamic range
- Better sensitivity
- Adapted for amplification-free nucleic acid detection



Universal Bacterial Biosensor

Harshini Munkundan and team

- **Comprehensive:** Identifies any bacteria from a small sample of water or blood
- **Robust:** Cannot be tricked by bacterial evolution or engineering
- **Fast:** Runs in 30 minutes from sample collection to result
- **Easy:** Can be performed by a non-expert at the point of need with no laboratory equipment
- **Flexible:** Interfaces with other detection systems, can be adapted to detect other oily biomarkers for diseases like cancer



2018
**R&D
100**
WINNER

2021 R&D 100 ENTRY

PEGASUS

Portable EnGineered Analytic Sensor with **Automated Sampling**
Novel analytical tool brings the lab to your fingertips

Innovative: Supplies point-of-need biochemical analysis

Rapid: Delivers accurate results in less than 30 minutes

Versatile: Detects a variety of biomarkers and biochemicals

Unique: Processes samples with a centrifugal microfluidic device

Robust: Provides a rugged, portable unit to resource-poor settings



 **Los Alamos**
NATIONAL LABORATORY



BIOLOGY

Jessica Kubicek-Sutherland
Loreen Stromberg
Zachary Stromberg
Alexis Bitzer
Shailja Jakhar



CHEMISTRY

Laura Lilley
Aaron Anderson
Joanna Casson
Stosh Kozimor



ENGINEERING

Kiersten Lenz
Aaron Anderson
Jennifer Harris
David Mascarenas
John Morales



THEORY


Sara DelValle
Nick Hengartner
Ben McMahon
Carrie Manore



COLLABORATORS

University of New Mexico
National Institutes of Allergy and Infectious Disease
Johns Hopkins University
Univ Nebraska Medical Center
Foundation for Innovative New Diagnostics
Medical University of South Carolina





Nature's response to infectious disease is agnostic

Why isn't ours?