

LA-UR-21-25684

Approved for public release; distribution is unlimited.

Title: Nondestructive Assay for Nuclear Safeguards

Author(s): Trahan, Alexis Chanel

Intended for: MTV Summer School webinar

Issued: 2021-06-16

Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

Nondestructive Assay for Nuclear Safeguards

LA-UR-21-

Alexis Trahan


June 17, 2021




Safeguards and the International Atomic Energy Agency

International Nuclear Safeguards

- Set of technical measures applied by the International Atomic Energy Agency (IAEA) to independently verify that nuclear materials are not being diverted to illicit purposes
- Safeguards play a central role in international nonproliferation efforts, i.e., in preventing the spread of nuclear weapons
- Tools and methods for implementing safeguards at nuclear facilities include:
 - Nuclear Material Accountancy
 - Nondestructive Assay (NDA)
 - Destructive Analysis (DA)
 - Containment & Surveillance
 - Environmental Sampling
 - Unattended and Remote Monitoring



The International Atomic Energy Agency



The International Atomic Energy Agency


- Currently, the IAEA is working to advance the following initiatives (among many others...)
 - Universal acceptance of the Additional Protocol
 - Safeguards-by-design
 - Integrated within a facility's design, covering safeguards and security
 - Unattended monitoring & data integration
 - Robust data management systems to reduce on-site inspector presence
 - State-level Concept
 - Assessing each State as a whole
 - Developing unified and consistent State-Level Approaches
 - Establishing safeguards measures based on path attractiveness rather than simply material attractiveness



Nondestructive Assay

Nondestructive Assay

- NDA is the most commonly employed technique for material accountancy
- A series of gamma or neutron detectors are typically used to measure radiation emitted from the sample of interest
- Energy, timing, and intensity of radiation may be correlated to isotope type and quantity in the sample



- Passive interrogation requires good signal intrinsic to sample (^{240}Pu , ^{252}Cf)
- Active interrogation requires fissile material or material prime for gamma interactions (^{235}U , ^{239}Pu)

Neutrons and Photons as NDA Signatures



Neutrons



Photons

Origins

- Spontaneous and induced fission
- (α ,n)
- Cosmic rays
- (p,n)
- (n,2n)
- (γ ,n)

Less common

- Nucleus (gamma-ray)
- Nuclear collision (gamma-ray)
- Electron cloud (x-ray)

Time and correlations

Energy

Low Z material

High Z material

^3He , Scintillators, fission chambers

HPGe, Scintillators, NaI, CZT, LaBr

Signal

Shielding


Detectors

Neutrons

History of Neutron Counting for NDA


- **TOTAL NEUTRON**
 - Record the total number of neutrons detected in a certain amount of time
 - Accurate assays can be obtained only for very few types of SNM
- **COINCIDENCE COUNTING**
 - Record the number of times two neutrons arrive within a set time window (gate)
 - Wide application for international safeguards
 - focused on verifying declared materials
- **NEUTRON MULTIPLICITY COUNTING**
 - Extension of neutron coincidence counting
 - Record the number of times we detect 2, 3, 4, etc. neutrons within a gate
 - It improves neutron assay accuracy dramatically by adding more measured information

Neutron Coincidence Counter




- ^3He neutron detectors
- Fission source (Pu) surrounded by neutron detectors
- Emission of multiple **prompt** neutrons from fission detected as coincident neutron events
- Multiplicity information is used to calculate the mass of fissile isotopes

Neutron Coincidence Counting




Neutron Coincidence Counting

Spontaneous
Fission




Detected
Neutrons


Rossi-Alpha Distribution




Epithermal Neutron Multiplicity Counter (ENMC)



- $\varepsilon = 65.0\%$
- $\tau = 22.0 \mu\text{sec}$
- 121 tubes
- 27 preamplifier channels




Advanced Experimental Fuel Counter



Advanced Experimental Fuel Counter

Active Doubles – Passive Doubles – Cf Doubles = Net Active Doubles



Gamma Rays

Photons

Generic Assay Equation

$$M_{SNM} = \frac{R_{Rad} \times CF}{Cal}$$

M_{SNM} = Mass of special nuclear material

R_{Rad} = Measured radiation rate (counts per unit time) from SNM item

CF = Correction for losses due to:

- item self absorption
- container absorption
- measurement system electronics

Cal = Calibration constant


FRAM

- FRAM is an isotopic analysis code nominally designed for plutonium and uranium.
- Fixed-energy Response-function Analysis with Multiple efficiencies.
- Self-calibration using several gamma-ray peaks.
- User-editable analysis parameters.
- Analyze gamma ray data from 30keV to >1MeV of HPGe, CdTe, CZT, and LaBr3 detector.


Peak Fitting

- FRAM uses linear least squares to fit the peaks of the HPGe spectra.
- FRAM uses a nonlinear least squares fit technique, combining the Powell's minimization method with the linear least squares fit to fit the peaks of the LaBr3 and CZT spectra.

Pu, 285-481 keV, LaBr3




U, 152-226 keV, CZT




Passive Gamma Emission Tomography (PGET)

- Three simultaneous measurements: gross neutron, gamma spectroscopy, and 2D emission tomography
- Create an axial image of emission locations to detect pin-level diversions
- Measurements take 3-5 minutes



Mayorov et al., IEEE, 2017




Miller et al., PNNL, 2017

- Neutron data are used for BU, spectroscopy data for CT or to verify non-fuel items
- Has been tested for burnups from 5.7-58 GWd/tU and cooling times from 1.9-27 years

Heat


Calorimetry

- Well-established, precise method of NDA
- Uses thermal power generated by radioactive decay in the sample to determine the mass of special nuclear material
- Heat flow calorimetry is most commonly used for materials control and accounting
- 60 Wheatstone bridge calorimeters currently being used for Pu and tritium measurements at LANL
- Bulk measurements can be taken without issues from absorption or self-shielding
- Takes much longer than other NDA techniques




Microcalorimetry


- Ultra-high energy resolution microcalorimeter technology offers a path to overcome NDA performance limits
 - *10-50x better energy resolution than semiconductor detectors*
- Improve economics and performance of safeguards and material accounting approaches




Decay energy spectroscopy



X-ray spectroscopy





Microcalorimeters measure the heat energy of individual photons or nuclear decays




$$E = C \times \Delta T$$

$$\Delta E \approx \sqrt{4kT^2C}$$



Potential to improve nondestructive measurements of ^{235}U enrichment



Safeguards Research at LANL


LANL Support

Over 50 years of support for the IAEA
through...

Technology development




Training



Expertise

Updated PANDA Manual

- An updated version of the PANDA manual (published in 1990) will be released soon
- New addenda
- New technologies
- New characterization methods
- New electronics
- ...much more!



Thank you!

Questions?