

LA-UR-21-26759

Approved for public release; distribution is unlimited.

Thermal Analysis and Flowability of Biomass Feedstock Title:

Semelsberger, Troy Allen Powers, Fiona Mae Author(s):

Intended for: Student Symposium (LANL)

2021-07-30 (rev.1) Issued:





Rheological Characterization of Biomass Feedstock for Alternative Energy Applications

8/4/2021

Fiona Powers

Mentor: Troy A. Semelsberger

Department: MPA-11

LAUR 21-26759

Background: Introduction to Bulk Solids


What are bulk solids?
 Loose solid material, e.g. sugar, salt, coal powder, ground corn stover.

Biomass Reusable Material





- What is biomass feedstock?
 Organic wastes, agriculture residues that can be sustainable for reusable material.
- What is a biorefinery?
 Converts biomass feedstock to useful fuels and chemicals



Corn Stover

Objectives

Short-term:

 Characterize the flow behavior of corn stover bulk solids with different moisture contents

Long-term:

 Use the measured flow behavior to re-design storage/handling equipment to mitigate the impact of moisture



Challenges

Flowability:

- Dependent on many environmental factors
- Particle size/geometry
- Variability in biomass feedstock (moisture content, chemical composition, storage conditions, etc.)

Problem with Moisture:


- Clogging of the processing/handling equipment
- Long down-time, high operational cost

How do we maximize flowability to transport bulk solids?


Need to characterize the flow (rheological) behavior as a function of moisture content



Flowability vs Moisture Content



What percentage of moisture yields the greatest flowability in bulk solids?





Sample Information:

- Low ash corn stover, ground to 2 mm in diameter
- Moisture: manually add DI water to corn stover
- Moisture content=mass of water/total mass

Instrument: Freeman Technology Powder Rheometer



Powder rheometer









Compression piston



I. Stability and Variable Flow Rate Test

Procedure

- 11 alternating conditioning/testing cycles with variable blade speed
- Tests 1~8: 100 mm/s
- Test 9, 10, 11: 70, 40 and 10 mm/s

Measurements

- Flow energy (mJ): measurement of amount of work required to move blade though powder
- Basic flowability energy=flow energy of Test 7

Purpose

- Access the stability in flow energy: tests
 1~8
- Access the dependence of flow energy on blade speed: tests 8~11







II. Compressibility Test

Procedure

- 1 conditioning cycle
- Increasing normal stresses are applied

Measurement

% Compressibility =
$$\frac{V_0 - V}{V} \times 100\%$$

Purpose

Obtain % compressibility vs. applied normal stress



III. Wall Friction Test

- Procedure
 - Conditioning Cycle
 - Sample Compression
 - Normal stress is applied with rotation (shearing) disks
 - Materials tested: Stainless steel, PTFE
- Measurements
 - Wall friction angle: amount of friction/resistance against material.
- Purpose
 - What material causes the least amount of friction between corn stover and transport material?









Images: Courtesy of Freeman Technology Ltd.


Stability & Variable Flow Rate Test





- 0%~15% moisture: increasing amount of energy needed to turn the blade
- 15% moisture: maximum BFE
- 15% ~75% moisture: decreasing amount of energy needed to turn the blade


Compressibility Test



- As normal stress increases, compressibility increases, and at a slower rate.
- Greater moisture content → greater compressibility
- Larger clusters form with increasing moisture → looser packing




Wall Friction Test



- Stainless Steel:
- (1.2 um surface roughness) Wall friction first increases and then reaches a plateau when exposed to 25% or greater moisture.
- PTFE: (Polytetrafluoroethylene): wall friction first decreases, plateaus after 15% moisture, and increases at 75% moisture exposure.



Summary and Conclusions



3D plot for the stability and variable flow rate test

- Based on the stability and variable flow rate test, the best and worst moisture content is 75% and 15%, respectively.
- The compressibility continuously increases with moisture content
- The relationship between wall friction angle and moisture are different on hydrophilic (stainless steel) and hydrophobic (PTFE) surfaces
- The data collected in this study can help with designing handling/ processing equipment for continuous flow and improved process economics



Future Work

 Compare chemical compositions and rheological properties of corn stover based on anatomical fraction.

How do the leaf, stalk and cob parts of the corn stover differ in terms of rheological properties?

Compare the rheological properties of high ash vs. low ash corn stover

How does high ash corn stover (biomass with high inorganic salt content) differ in rheological characterization and handling needs?

Use existing data for hopper design

Can we use existing handling equipment to economically transport and handle biomass feedstock?

Rheological data collected will be used in an upcoming publication



Acknowledgments

- Mentor: Dr. Troy A. Semelsberger
- Colleagues: Dr. Juan Hilario Leal, Dr. Lily Cheng, David Gao, Jarrod Ronquillo
- This research was supported in part by an appointment with the Energy Efficiency & Renewable Energy (EERE) Energy Storage Internship Program sponsored by the U.S. Department of Energy (DOE), EERE Advanced Manufacturing Office (AMO). This program is administered by the Oak Ridge Institute for Science and Education (ORISE) for the DOE. ORISE is managed by ORAU. All opinions expressed in this paper are the author's and do not necessarily reflect the policies and views of DOE, ORAU, or ORISE.

