

LA-UR-21-27425

Approved for public release; distribution is unlimited.

Title: Improved Temperature Control for Measuring the Humidity Dependence of Aerosol Optical Properties

Author(s): Benedict, Katherine Beem
Moosmuller, Michael Urs
Gorkowski, Kyle Joseph
Dubey, Manvendra Krishna

Intended for: Summer Internship Presentation

Issued: 2021-07-28

Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

Improved Temperature Control for Measuring the Humidity Dependence of Aerosol Optical Properties


EES14 - ACES
Center for Aerosol Forensic Experiments
(CAFE)

By: Urs Moosmuller
Mentors: Kyle Gorkowski,
Katie Benedict,
Manvendra Dubey




Importance of Measuring Optical Properties of Humidified Aerosols

- **Climate Modeling** – quantitative methods to simulate interactions of the important drivers of climate, including atmosphere, oceans, land surface, and ice.
- **Climate Forcing** – the physical process of affecting the climate on Earth through different forcing factors. (variations in solar radiation levels, volcanic eruptions, changing albedo, changing greenhouse gases, etc.).
- **Cloud Formation** – water vapor in the air condenses into visible water droplets or ice crystals through condensation of saturated parcels of air.




Optical Properties and Humidified Aerosols

Optical Properties




Humidified Aerosols



Relative Humidity and Temperature


- As the temperature of water increases, the kinetic energy of its molecules increases and the number of molecules transitioning into water vapor also increases, thereby increasing the vapor pressure.

$$\text{Relative Humidity (RH)} = \frac{\text{partial pressure of water vapor}}{\text{equilibrium vapor pressure of water}}$$



CAPS PM_{ssa} Monitor


- Provides simultaneous measurement of both particle-based total optical extinction and scattering.
- Produces a direct measurement of the particle single scattering albedo (SSA), the ratio of scattering to total extinction.
- Laser wavelength blue light 450nm
- Measure ambient aerosol directly in the CAPS monitor but also simulate different levels of humidity to characterize the roll of water uptake on optical properties.



DUAL CAPS DESIGN




1. Temperature Controller and water cooler/heater added to water system.
2. Temperature Controller and air cooler/heater added to the aerosol system.
3. Second CAPS instrument.



Temperature Control to Humidify Aerosol

Thermocouple Temperature Sensor



Temperature Controller




Water Reservoir



Water



Aerosol Humidifier




Water Heater/Cooler

Temperature Control for Aerosols


Thermocouple Temperature Sensor



Air Cooler




Thermocouple Temperature Sensor




CAPS Comparison

Fullerene Soot (Black Carbon Surrogate)
Time Period: 7/22/2021 8:00 – 10:30am

Direct Comparison



Distribution of the Difference (CAPS2-CAPS1)



Percentage Error:

Extinction: 4.9%

Scattering: 6.2%

Slope:

1.054

1.072



CONCLUSION

- There aren't significant biases between the CAPS systems so we can directly compare the dry and humidified measurements.
 - Temperature control on the dry system creates a more robust and stable measurement.
 - On the wet system we now can more precisely control the different levels of humidification.
 - Now we can start looking at ambient aerosols to look at their behavior and changes to the optical properties with water uptake.
- 