

LA-UR-21-27755

Approved for public release; distribution is unlimited.

Title: A new approach to measure the light absorption of aerosols

Author(s): Gorkowski, Kyle Joseph
Jordan, Spencer Halloran
Benedict, Katherine Beem

Intended for: SULI Student Presentation

Issued: 2021-08-04

Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.




A new approach to measure the light absorption of aerosols

Spencer Jordan

Mentored by: Kyle Gorkowski and Katie Benedict


Solar Absorption by Atmospheric Aerosols



- Influence how much solar radiation gets absorbed/reflected
 - Controls the warming and cooling of the planet
- Understanding the wavelengths at which aerosols absorb light is important to climate models
 - Better parameterization
- In this study – used organic dyes to validate the methods
 - Sunset yellow, nigrosin, fluorescein-sodium-salt
 - Molecular similarities to common Organic Aerosols






Current Methods of Aerosol Absorption Measurement

- Aerosol instruments typically measure absorption and scattering of light
 - Inversion algorithms are used to retrieve the refractive index
 - This process has a high uncertainty
- Photoacoustic methods
 - Measure the absorption of single particles with laser beam
 - Single wavelength instruments
- Filter based methods
 - Collects samples on filter then measures attenuation
 - Biases exist due to light scattering off the filter paper
 - Can get overloaded with aerosol




Bringing Thin-film Techniques to Aerosol Measurement

- We used a spin coater to create thin, uniform layers of solutions onto glass slides
- Slides could then be scanned with a spectrometer to measure reflection and transmission of light
- Spectrum is used to calculate both the thickness and refractive index




Spectrometer Analysis




Direct Chemical Doping Method


Injection of solution onto
slide - while spinning at
3000 rpm



Spun for an additional
55 seconds at 3000 rpm



Resultant
thin-film





Poly-acrylic Acid Doped with Dye

Fused Silica Glass Slide


- Poly-acrylic acid (PAA) and water solutions doped with several concentrations of sample dyes
- Very low non-uniformity (< 10 nm), consistent thicknesses

Deposition of Aerosolized Dyes




- Aerosols were collected in small areas on the surface of the slides
 - 'Layers' were less uniform than the doping method
- Future work will look at improved impaction methods

Sunset Yellow: All Collection Methods Compare Well




Challenges of Fluorescein Sodium Salt



- Not fully miscible in water
 - Acetone was used to bring it into solution
- New layering process was used
 - PAA layer immediately followed by dye solution


Aerosolized Nigrosin Samples Agree well with Literature



Doped and Aerosolized Slides Generally Agreed

Poly-acrylic Acid Doped with Dye

Fused Silica Glass Slide



- Future Works
 - Develop aerosol collection process
 - 3D printed slide holders
 - Aerosol Impaction Simulations
 - Deploy in the field
 - Urban pollution
 - Wildfire observations