

#### LA-UR-21-29673

Approved for public release; distribution is unlimited.

Title: Exploring Interactions among Nanocrystals with Multidimensional

Nanoscopy


Author(s): Ryan, Duncan Patrick

Intended for: Idaho State University physics colloquium. 10/4/21. Virtual.

Issued: 2021-09-29



# Exploring Interactions among Nanocrystals with Multidimensional Nanoscopy




Duncan P Ryan, Los Alamos National Laboratory









# **Quantum Dot Nanocrystals**




#### **Fluorescent Particles**

2 – 10 nm diameter
Tunable emission in optical region
High quantum yield (bright)



Halber, Spectrum (2011)

#### **Structural components:**




Core/shell are semiconductor material (ex. CdSe/ZnS), define basic emission properties through intrinsic bulk material bandgap.

Ligands protect and stabilize NCs from environment.

# **Quantum Dot Nanocrystals**




The structure of NCs energy levels are the result of quantum confinement.



# **Quantum Dot Spectral Properties**





An **exciton** is the combined pair of an electron and hole.

Shells confine excitons to core of the nanoparticle.

# **Quantum Dot Nanocrystals**




#### **Displays**



Campbell, Apple Insider (2014)


#### Research & Imaging



**Life Technologies** 

# 45 mm






#### **Properties of NPs:**

Consist of a core/shell structure, often made up of two different semiconductor materials.

Strongly absorbing in the UV and high quantum yields in the visible (for a CdSe system).

The size of the NP determines the energy levels: smaller cores have larger bandgaps.


#### **Photovoltaics & Solar**



Los Alamos Science & Tech Magazine


# **Networks of Nanoparticles**









Talapin, Nano Letters (2007)



Small clusters of NPs are a model system to study interactions with single-molecule techniques.






Xu, JACS (2011)


# **Fluorescence from Interacting NP Clusters**






Clusters are formed by slowly stripping ligands in solution to form surface-to-surface attachment.

CdSe/CdS or CdSe/ZnS semiconductor NPs




A distribution of cluster sizes and geometries are produced

#### **Emission fluctuations: blinking**



Individual NPs exhibit binary or telegraph blinking




Clusters of NPs will interact, resulting in **enhanced** blinking with rapid intensity modulations that typically lack discrete intensity levels.

# **Optical Super-resolution Imaging of NP Clusters**



#### Where in a NP cluster is the emission coming from?

To study collective fluorescence behaviors in NP clusters, turn to super-resolution imaging.



NP clusters are smaller than the diffraction-limit of their fluorescence: ~8nm. Resolving individual NPs in a cluster requires nanometer resolution.

# **Super-resolution Microscopy**



Particles can be resolved below the diffraction-limit of visible light. From 250  $\mu m$  to 5 nm.


Super-resolution microscopy is a collection of techniques:

- Fluorescent probe physics/dynamics
- Instrumentation
- Analysis algorithms


Many super-resolution methods are **single-particle** approaches: particles are examined individually

Super-resolution microscopy earned the 2014 Nobel Prize in Chemistry for Eric Betzig, Stefan Hell, and William Moerner.

Wide-field image





Super-resolution image (dSTORM)




# **Super-resolution Imaging of NP Clusters**










Correlated fluorescence super-resolution and SEM imaging compared implied cluster geometry with the true geometry.

"Hot-spots" emerge as dominant emission regions within NP clusters

### **Correlation with Cluster Size**






Total cluster emission intensity does not scale simply with cluster size. Additional NPs do not increase emission as significantly as independent NPs.

Exciton loss mechanisms suppress emission intensity.

# **Energy Transfer**





#### **Model for energy transfer:**

An excited donor can (a) emit a photon or (b) transfer the exciton to a nearby acceptor.

Because fluorescence and energy transfer are competing processes, the direct donor emission will appear quenched with a shorter lifetime and at a lower intensity.

Emission from an acceptor will be enhanced due to the additional excitons funneled through energy transfer and will have the intrinsic lifetime of an individual NP.


#### **Energy transfer has signatures in:**

- Intensity (enhanced blinking)
- Spatial (emission center moves)
- Lifetime
- Spectral
- Orientation

Combining localization information with lifetime information can reveal dynamics and relationships among energy transferring NPs.

# **Orientation Components of Energy Transfer**






Orientation factor  $\kappa^2$  strongly dictates coupling efficiency in FRET.

Other coupling strength parameters: spectral overlap and separation distance

#### NPs emit from a 2D dipole




Bene, L.; Bagdány, M.; Damjanovich, L. "Checkpoint for Helicity Conservation in Fluorescence at the Nanoscale: Energy and Helicity Transfer (HFRET) from a Rotating Donor Dipole." Biophysical Chemistry 2018, 239, 38–53.

As 2D emitters, NPs require additional considerations to determine  $\kappa^2$ .

In fixed nanostructures, using the typical averaged  $\kappa^2$  value is not appropriate.

# **Orientation Mapping of NP Clusters**





#### Different orientations from diffraction patterns.

Irregular shapes and shell thicknesses from contrast imaging. Not ideal NPs for energy transfer studies.


Never use commercial NPs: reviewers will hate it.

# **Polarization Imaging for Orientation Measurement**



#### **Conventional confocal setup (continuous monitoring)**


- Established method for non-imaging based polarization and orientation measurements.
- Low throughput
- Does not super-resolve emitter locations



Lethiec, C.; et al. "Measurement of Three-Dimensional Dipole Orientation of a Single Fluorescent Nanoemitter by Emission Polarization Analysis." Phys. Rev. X 2014, 4 (2), 021037.

#### Wide-field imaging setup (fixed polarization angles)

- High throughput & intrinsic drift correction
- Super-resolves emitter locations
- Orientation information from intensity distribution into 4 polarization channels




# **Polarization Imaging for Orientation Measurement**



#### **Super-resolved orientation**

4 images


Calibration corrects for optical path efficiency differences




#### **Super-resolved orientation + spectral**

8 images

Calibration also considerations of polarization responses of dichroic mirrors








Ryan, D. P. *et al.* Quantum dot clusters visualized by polarizationand wavelength-resolved nanoscopy. *In prep*.

# **Orientation Imaging**










Anisotropy Orientation angles

Intensity anisotropy between orthogonal polarization pairs measures the shape of the projected ellipse.


The orientation of a NP also requires an appropriate model of imaging system to determine the polar angle  $\theta$ .


The azimuthal angle  $\phi$  is two-fold degenerate.




# **Orientation of a Single NP**











Isolated NPs localize to a single position.

A single NP has a well-defined wavelength at a given intensity (bright state).

#### **Orientation of a Small NP Cluster**









Because of PSF asymmetry, localization information from orthogonal pairs of polarization channels are required to generate a single super-resolution image.

Super-resolved images from the final two channel pairs produce similar shapes.







Spectral-intensity distribution follows energy transfer model: highest emission is from longer wavelength NPs and lower intensity donor emission is shorter wavelength.


Orientation signature indicates a single orientation.


# **Orientation of a Large NP Cluster**











This larger NP cluster shows structural features, but a common alignment among the NPs, indicating oriented attachment.


The spectral signature of larger NP clusters tends to not exhibit the energy transfer relationship between wavelength and intensity.

The mild orientation variation is only present in the polar angle  $\theta$ .

Color coding (blue/red) the orientation relative to the mean emission wavelength shows shorter wavelength emission originates from NPs with small tilt angles.

#### **Conclusions**







- Demonstrated orientation imaging using a wide-field polarization microscope.
- NP clusters have a narrow distribution of orientations, indicating particle alignment within cluster.

#### **Acknowledgements:**



Somak Majumder Jennifer A Hollingsworth Chris J Sheehan James H Werner Peter M Goodwin





Megan Dunlap Martin Gelfand Alan Van Orden



Research presented in this presentation was supported by the Laboratory Directed Research and Development program of Los Alamos National Laboratory under project number 20180189ER.

LA-UR #21-xxxxx