

LA-UR-21-30551

Approved for public release; distribution is unlimited.

Title: NHMFL Pulsed Field Facility Capabilities Overview

Author(s): Wampler, James Paris
Zapf, Vivien

Intended for: For use in conversations with other scientists

Issued: 2021-10-22

Disclaimer:

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.



NHMFL Pulsed Field Facility Capabilites Overview

James Wampler, Vivien Zapf

Los Alamos National Laboratory (LANL)



Outline:

1. Magnets
2. Probes
3. Process

3N070A

1. Four (4) 65 T Magnets
2. 75 T Magnet
3. Mid Pulse magnet
4. 100 T Magnet in... soon?

@Los Alamos National Lab

A User Facility:

Magnetization

Resistivity

Electric polarization/dielectric

Magnetostriction

Magnetocaloric

dHv α , sdH: fermi surface


Tunnel-diode oscillators: rf

penetration depth

Optical absorption

Credit : Vivien Zapf

- Include 73 T and mid pulse
- include pic of generator in the warehouse





NATIONAL HIGH

MAGNETIC FIELD LABORATORY

Four (4) 65 T Magnets

@Los Alamos National Lab



Capacitor bank ~ 4 MJ


A User Facility:

- Magnetization**
- Resistivity
- Electric polarization/dielectric
- Magnetostriction**
- Magnetocaloric**
- dHv α , sdH
- Tunnel-diode oscillators
- Optical absorption

Credit : Vivien Zapf

73 T Magnet

@Los Alamos National Lab



A User Facility:

Magnetization

Resistivity

Electric polarization/dielectric

Magnetostriction

Magnetocaloric

dHv α , sdH


Tunnel-diode oscillators

Optical absorption


Credit : Vivien Zapf

@Los Alamos National Lab

101 Tesla multishot (*world record*)



60 Tesla long pulse (*world record*)



A User Facility:

Magnetization

Resistivity

Electric polarization/dielectric

Magnetostriction

Magnetocaloric



dHv α , sdH

Tunnel-diode oscillators

Optical absorption

Credit : Vivien Zapf

High Field Magnetization



- Measurement : voltage across the pickup coil $\propto \frac{dM}{dH}$.
- Integrate $\frac{dM}{dH}$ with field to obtain the magnetization.


Black: **Sample Cylinder – 1 mm inner diameter, 5 mm length.**

- Sample must be held still through N Grease or similar


Red: Pickup coil – Cu - 1000 turns (CCW) + 500 (CW) + one turn of compensation coil

Minimum temperature: ~0.6 K

High Field Magnetization




Raw data



Integrated magnetization data

High field Polarization




Samples must be insulating

We measure I_{me} and integrate with time to obtain the change of polarization


$$Q = \int I_{me} dt$$

$$\frac{dQ}{dt} = \frac{dQ}{dB} \frac{dB}{dt} = I_{me}$$

Current is proportional to dB/dt -> faster field sweep rate increases the signal size!




High field Polarization



Raw data

Polarization



Integrated delta P data

Polarization Sample Prep

- Sample must be insulating ($> 10^9 \Omega\text{cm}$)
- ~ 1 pA minimum signal
 - Current depends on size of surface, so big pads is good.
- Flat bulk or thin film – two well defined surfaces
 - Multilayer for bottom lead or some other way of accomplishing that
 - Electrodes on both sides
- Au or Pt wires
- Can measure 2 samples simultaneously; 1 in plane, 1 out of plane
- 5 mm^2 max
- $< 2 \text{ mm}$


Magnetocaloric effect

- Measuring temperature change of sample as a function of field in the adiabatic condition.



AuGe: Semiconductor – Thermometer.

Au: Electrical Contact




Measure resistance of AuGe - > Convert it to temperature.

(Can also measure resistance and not convert to temperature!)

- Can sputter on the contacts at LANL or at UCSD
- Same technique used to measure AC resistivity
 - AuGe not necessary
- 1 M Ω practical limitation

Picture of α -RuCl₃




Magnetocaloric effect – example data



PHYSICAL REVIEW B 102, 214432 (2020)



High Field Magnetostriction

Limitation: Mostly only useful for bulk samples!



Credit : fbgs.com

High Field Magnetostriction – example data



Applying for magnet time

- The Pulsed Field Facility (part of the NHMFL) is a user lab
- 3 page long proposal covering the proposed science and impact
- 1 page long experimental proposal covering the experiment and prior work
- Deadlines are rolling – announced as we approach
 - Submit any time
 - Next deadline Nov 12, 2021
- I am happy to support your magnet times

<https://nationalmaglab.org/user-resources/request-magnet-time>

Maglab in the time of Covid

- Only vaccinated individuals may visit the lab
 - Non-US residents require additional testing upon arrival
- NM very safe. Los Alamos is particularly safe.
 - LANL has a vaccine mandate effective Oct 15
- Flying still risky. 14 hour drive still lengthy.
- Since you are applying for time well in advance, there's always unknowns.