

LA-UR-21-31467

Approved for public release; distribution is unlimited.

Title: Development of a pixelated BaF2 test bed for timing applications

Author(s): Jordan, Tyler Alexander

Intended for: Slides for PhD qualifying exam

Issued: 2021-11-18







Development of a pixelated BaF2 test bed for timing applications

Tyler Jordan

November 19, 2021

Berkeley UNIVERSITY OF CALIFORNIA

LA-UR-xx-xxxxx

Outline

- Motivation: What can BaF₂ do?
- Background: How does BaF₂ do it?
- Focus: Time-of-flight positron emission tomography (TOF-PET)
- Prototype development: Design considerations
- Prototype development: Components and assembly
- Next Steps: Characterizing detector performance




Introduction



Motivation

- BaF₂ one of the fastest known scintillators, emitting at 195 & 220 nm with subnanosecond decay constant
- Fast emission only comprises ~13% of total light
 - Most light is emitted at 310 nm, 600 ns decay constant
 - Count rate permitting → fast component for timing, slow component for energy resolution



Motivation

Medical imaging

- Time-of-flight positron emission tomography (TOF-PET)
- Hadron therapy range monitoring

Low energy physics

- Neutron capture cross section measurements
- Nuclear lifetime measurements

Timing applications

- High time resolution
- High count rate

Materials science

 Positron annihilation spectroscopy

High energy physics

Calorimetry



Motivation

- Ongoing BaF₂ research largely focuses on material and/ or photodetector characterization → small, single crystals
- Different applications require different detection system geometries
- Some applications require slow component rejection

Goal: Develop a pixelated BaF₂ test bed to characterize the performance of various configurations for timing applications.

- Focus: TOF-PET → design inspiration
- Deployable form factor
- Variable system geometry



Background




BaF₂ scintillation mechanisms

Slow component

- Self-trapped exciton (STE) luminescence
- Self-trapped hole: two adjacent F⁻ ions share a hole, effectively forming F₂⁻
- A self-trapped hole captures an electron, forming a self-trapped exciton which subsequently decays


Fast component

- Core-valence (CV) luminescence
- An electron from the 5pBa²⁺ band is excited to the conduction band
- The hole left behind is quickly annihilated by an electron from the 2pF⁻ band






 Heating BaF₂ activates thermal dissociation of STEs






- Doping BaF₂ with rare earths preferentially quenches STE luminescence
 - Inhibits STE formation
 - Provides STE dissociation channels





 Filtering BaF₂ emission allows selective collection of the fast component







Example: HEP colliders – High bunch crossing rate and intense product showers make slow component suppression ideal.




An overview of TOF-PET



Conventional PET – Principles

 Antiparticle pair: electrons and positrons annihilate one another, converting the rest mass of both particles into electromagnetic energy




- PET leverages these kinematics to localize a positron emitter
 - Near-coincident detection of back-to-back gammas defines a line of response (LOR)
 - Multiple LORs indicate source position via their intersection
- Conventional PET uses TOF only for event identification



TOF-PET – Principles


- TOF can localize source position along each line of response
 - The difference in time of arrival of annihilation quanta yields the distance from field of view (FOV) center to the source
- Ideally, enables direct source localization without image reconstruction
 - Obviates the need for measurements at multiple detector angles
 - Requires 1.5 mm resolution along the line of response → 10 ps FWHM coincident time resolution (CTR)





TOF-PET – Principles

TOF Improves image signal-to-noise ratio (SNR)



Gaussian probability kernel FWHM determined by coincidence time resolution (CTR)





TOF-PET – Materials


- Historically, BGO is the material of choice for conventional PET
 - TOF-PET not possible
- LSO is considered the best scintillator overall for PET
 - 300 ps CTR for two crystals → TOF-PET feasible

Property	BaF2	LSO:Ce	BGO
Density ρ (g cm ⁻³)	4.88	7.4	7.13
Effective atomic number, $Z_{\rm eff}$	53	66	75.2
Photon absorption α @511 KeV (cm ⁻¹)	0.085	0.28	0.336
Radiation length X_0 (cm)	2	1.1	1.12
Intrinsic light yield, LY_{intr} (ph/MeV)	$1400^{a\ b\ h} - 7000^{b\ i}$	$40\ 000^{\mathrm{b}}$	10 000 ^b
Decay time τ (ns)	$0.6-0.8^{\rm h}/620^{\rm i}$	$22/44^{\rm b}$	$46/365^{b}$
Photon fraction @ 0.5 MeV	0.19°)	0.34°	0.43 ^c
Emission peak(s) λ_{max} (nm)	$195^{\rm i}220^{\rm i}310^{\rm h}$	420^{b}	480^{b}
Refractive index (RI) @ λ_{max}	$1.56^{\mathrm{d}}1.55^{\mathrm{d}}1.50^{\mathrm{d}}$	1.82 ^b	2.1 ^b
Melting point (°C)	1280 ^e	2150 ^f	1050^{f}
$Cost (\$ cm^{-3})$	15 ^g	60 ^g	35 ^g



TOF-PET – Detector structure

- Block detector
 - Tightly packed array of crystals, separated by reflective material
 - Monolithic light guide
 - Array of PMTs
 - Interaction position reconstructed from light intensity ratios across PMT array





TOF-PET – Problems with BaF₂

- Earliest TOF-PET systems utilized BaF₂, but material property difficulties and the advent of BGO all but halted development
- Low density, low effective Z
 - Poor efficiency
- Low light output
 - Poor energy resolution (photopeak identification)
 - Poor position reconstruction assuming standard block detector
- UV emission
 - Light collection challenges



Prototype development

Design considerations

Goal: Develop a pixelated BaF₂ test bed to characterize the performance of various configurations for timing applications.

- Focus: TOF-PET → design inspiration
- Deployable form factor
- Variable system geometry

- Low density, low effective Z
 - Larger pixels
- Low light output
 - Read out each pixel individually
- UV emission
 - Total internal reflection
 - UV sensitive photodetector


- High time resolution
 - Total internal reflection
 - Low transit time spread photodetector
- Modular
 - Easily configurable experimental geometries



Photonis Planacon


- 59 mm x 59 mm x 21 mm
- Sapphire window → Good UV transmission
- MCP gain stage
 - 35 ps transit time spread
- 8x8 grid of anodes
 - 5.9 mm square anodes, 6.5 mm pitch
- Two Samtec QRM8 connectors for signal transmission
- Potted HV divider (SHV to multi-pin LEMO)








BaF₂ array – Pure BaF₂ pixels, 4.6 mm x 4.6 mm x 30 mm




- Stack plates to form a 3D pixel frame
- Rounded corners on pixel slots enforce an air gap around each pixel
- Align pixels over anodes
- Laser-cut black acetal





Slow component suppression

- Slow component suppression does not necessarily improve TOF-PET performance
- Instead, facilitates system characterization for high count rate applications
 - Slow component suppression via heating isn't practical for instrumentation
 - Doping can also degrade fast component intensity
 - Optical filtering is relatively easy to implement
- Interference filters can be manufactured to match pixel size, preserving optical segmentation
- Complication: filter transmission is angle-ofincidence (AOI) dependent







Filter optimization

- Choice of long wavelength cutoff must account for AOI distribution of collected light
- Simplified 2D light transport, AOI-dependent filter transmission, and QE calculations employed to optimize long wavelength cutoff for fast photoelectron production

$$\lambda_{\theta} = \lambda \sqrt{1 - \left(\frac{n_0}{n_{eff}} \sin \theta\right)^2}$$








Filter optimization

- Increasing the filter cutoff from 227 nm to 242 nm preserves more off-angle fast light, at the cost of slightly more slow component transmission
- Slow component suppression via doping would allow for an even higher cutoff





Long wavelength cutoff	Number of fast photoelectrons (pe)	[Fast pe] / [Total pe]
227nm	1.74 / MeVee	0.999
242nm	8.17 / MeVee	0.985



DAQ Electronics

- Timing: CAEN VX1290N-2eSST TDC
 - 16 channels
 - 25 ps resolution
 - Timing signal can be picked off from the HV divider
 - Modify HV divider provided by Photonis (if possible)
 - Design custom HV divider
 - A 10 ps resolution model coming soon...
- Spectroscopy: CAEN VX2740 digitizer
 - 64 channels
 - 125 MS/s \rightarrow 8 ns samples
 - 2 V full scale range, 16-bit
 - Spectroscopy helps with event identification in TOF-PET








Detector assembly




BaF₂ array construction


~60 plate stack forms pixel frame



BaF₂ array construction

- Air gap around each pixel enables total internal reflection
 - The top of each pixel shows the color of the surface the array sits on








Detector assembly

Stack up (single pixel)




 Filtered BaF₂ array sits at the end of an aluminum housing





Detector assembly


 Planacon is air-coupled to the array





BaF₂ array construction

- Four corner screws for mounting
- 4 detectors assembled!







Next steps



Simulations

- Toy Monte Carlo models are useful for system performance estimation and building intuition, but they are not sufficient
- GEANT4: G4Scintillation and G4OpticalPhoton
 - One simulation spanning from radiation incident on a BaF₂ crystal to digitized waveforms
- Validation via comparison to experimental data enables design optimization in simulation space
 - Optimization parameters
 - Crystal dimensions
 - Crystal emission spectrum (i.e., doping)
 - Filter wavelength cutoff
 - Optical coupling




Data acquisition software

- VX2740 digitizers represent a significant departure from previous CAEN product design
 - Hardware: direct communication via USB/ ethernet rather than VME backplane
 - Software: string-based register interface rather than hexadecimal
- This is an opportunity to build new DAQ control software from scratch, rather than making discontinuous modifications to existing code
 - Must incorporate CAEN TDC, which isn't supported by the current software
- Set up time signal pickoff from HV divider
 - Requires either existing circuit modification or custom circuit design



Experiments


- Filter characterization
 - 1" x 1" cylindrical pure BaF₂ crystal
 - Simultaneous collection of filtered and unfiltered light leads to filtered photopeak identification
 - Feasibility study for filtered spectroscopy






Experiments

- PET-like system performance
 - Mounting hardware to secure detectors to ring stands
 - 4 detectors have been constructed
 - ²²Na source localization







Final thoughts



Acknowledgements and Disclaimer

I offer my sincere thanks to Madison Andrews, Edward McKigney, Michael McCumber, Sy Stange, and Cameron Bates for their invaluable guidance and continual assistance throughout this work.

Los Alamos National Laboratory has provided all of the resources that made this work possible, for which I am exceedingly grateful.

This material is based on work supported by the Department of Energy National Nuclear Security Administration through the Nuclear Science and Security Consortium under Award Number DE-NA0003180.

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor any of their employees, makes any warranty, express or limited, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.



References

- [1] M. Laval et al., Nucl. Instrum. Methods Phys. Res. 206, 169 (1983)
- [2] C.W.E. van Eijk, *Journal of Luminescence* **60**, 936 (1994)
- [3] P. Dorenbos et al., *IEEE Trans. On Nucl. Sci.* **40**, 424 (1993)
- [4] P. Schotanus et al., Nucl. Instrum. Methods Phys. Res. A 238, 564 (1985)
- [5] P. Schotanus et al., Nucl. Instrum. Methods Phys. Res. A 281, 162 (1989)
- [6] S. Diehl et al., J. Phys.: Conf. Ser. 587 012044 (2015)
- [7] M. Conti, Physica Medica 25, 1 (2009)
- [8] S. Gundacker et al., Phys. Med. Biol., 66 114002 (2021)
- [9] M. Andrews et al., LA-UR-20-23526 (2020)
- [10] RoentDek, Manual for MCP ToF detector DET25/40/75/100/120(P)
- [11] T.A. Chubb, J. Opt. Soc. Am. 46, 362 (1956)
- [12] M.R. Farukhi et al., IEEE Trans. On Nucl. Sci. NS-18, 200 (1971)



Questions?



Toy Monte Carlo models



Light transport

- N events with E energy deposited
- For each event
 - Position sampled uniformly throughout crystal
 - Number of photons sampled from Poisson distribution with mean E*TLY
 - For each photon
 - Emission component sampled according to FLY/TLY
 - Emission time sampled from exponential distribution
 - Wavelength sampled from digitized spectrum
 - Direction sampled isotropically
 - Photons moving the wrong way are thrown out
 - Time to hit crystal face in each direction computed from position and direction
 - Indices of refraction and critical angles computed from wavelength



Light transport

- Based on times and TIR criteria, the following are calculated
 - Collected (bool)
 - TIR at collection face (bool)
 - AOI on collection face
 - Refracted AOI after escaping through collection face
 - Transit time in crystal
 - Collection time (emission time + transit time)
 - Number of reflections



Filter transmission

- For each photon
 - Uncollected photons are thrown out
 - Effective photon wavelength calculated from wavelength and RCAOI, according to transmission AOI dependence
 - Photons with effective wavelength outside of digitized transmission range are thrown out
 - Transmitted (bool) sampled according to probability at effective wavelength




Quantum efficiency

- Runs off of either light transport or filter transmission output
- For each photon
 - Uncollected/ untransmitted photons are thrown out
 - Photons with wavelength outside of digitized QE range are thrown out
 - Converted (bool) sampled according to probability at wavelength




Schematics

Signal cabling – QRM8 to JTAG






Signal cabling – QRM8 to MCX





Corner radius

Filleted corners maintain an airgap around each pixel



$$r = x\sqrt{2} \Longrightarrow x = \frac{r}{\sqrt{2}}$$


$$a + x = r \Rightarrow a + \frac{r}{\sqrt{2}} = r$$

 $\Rightarrow a = r\left(1 - \frac{1}{\sqrt{2}}\right)$

- A slightly more complicated equation accounts for machining tolerance
- Air gap: 0.05mm
- Corner radius: 0.11mm
- Pixel size: 4.6mm
- Hole width: 4.7mm

BaF₂ array structural components – pixel plates


- Stack up to form a 3D pixel frame
- Rounded corners on pixel slots enforce an air gap around each pixel
- Align pixels over anodes
- Laser-cut black acetal
 - Nonconductive material avoids potential interference with electron optics





BaF₂ array structural components – filter plates


- Align filters under pixels
- Laser-cut black acetal
 - Nonconductive material avoids potential interference with electron optics





BaF₂ array structural components – spacer plates


- Diaphragm blocking unfiltered light to account for pixel-filter size mismatch
- Laser-cut black acetal
 - Nonconductive material avoids potential interference with electron optics





Detector housing – cap


- Contains the pixel-filter-Planacon stack-up
- Wall thickness set such that two adjacent detectors will have 3 columns of "dead pixels" between them
- Aluminum
 - Low Z material for minimal gamma attenuation





Detector housing – back plate

- Ports for signal and HV cables
- 4 mounting holes
 - #10-32 screws
- Aluminum
 - Low Z material for minimal gamma attenuation





Principles of radiation detection

What is a radiation detector?

- A transducer that converts radiation. into electric current
 - The output current encodes information about radiation interactions within the detector
- An analogy: cryptography
 - Goal: to translate information's syntax while preserving its substance
 - Physics handles the encryption
 - We are responsible for decryption, so we must understand the physics




By Hubert Berberich (HubiB) - Own work, Public Domain, https://commons.wikimedia.org/w/index.php?curid=25964875



Two stages of general detector response


- 1. Excitation: a free, energetic charged particle interacts Coulombically with the electrons in the detector active volume
 - Direct, e.g., a decay alpha in a gas detector
 - Indirect, e.g., a photoelectron in a semiconductor detector
- 2. Collection: the intervening process(es) between excitation and the final output current
 - Intrinsic, e.g., scintillation
 - Extrinsic, e.g., pulse shaping





Two stages of general detector response – timing

- Extrinsic collection processes can be optimized for timing
 - Ex: high time resolution photodetectors in a scintillation detector
- Excitation time and intrinsic collection time are determined by detector material properties
 - Ex: stopping power, electron/ hole mobilities (semiconductors), fluorescence decay time (scintillators)
 - Define the fundamental limits of detector performance
- New material properties → new detectors






Two stages of general detector response – timing

- Extrinsic collection processes can be optimized for timing
 - Ex: high time resolution photodetectors in a scintillation detector
- Excitation time and intrinsic collection time are determined by detector material properties
 - Ex: stopping time (scintill
 Define the fill

 The discovery of the BaF₂ fast scintillation component prompted a wave of BaF₂ R&D
- New material properties → new detectors





BaF₂ Timeline


A brief history of BaF₂ R&D


 BaF₂ was originally investigated for vacuum ultraviolet optical elements

Non-hygroscopic

Very slightly water soluble









- 1971: Farukhi et al. observe scintillation under ionizing irradiation
 - Emission peaked at 325 nm
 - 630 ns decay time constant
 - Scintillation mechanism not well understood, attributed to modification of hole centers



A brief history of BaF₂ R&D

- 1971-1982: very little progress, still mostly focused on UV transmission
- 1982: Ershov et al. and Laval et al. independently observe the fast scintillation component
 - Laval et al. propose TOF-PET as a potential application








A brief history of BaF₂ R&D

- 1971-1982: very little progress, still mostly focused on UV transmission
- 1982: Ershov et al. and Laval et al. independently observe the fast scintillation component





BaF₂ in the 21st century

Photodetectors

- Solar-blind PMTs
- **VUV SiPMs**

Scintillation characterization

Discovery of "very fast" emission below 240 nm with 100 ps decay time

- **Medical imaging**
- **High energy physics**
- Low energy physics
- **Materials science**

Other applications

- Oil well logging
- Thermal neutron detection
- Double beta decay detection
- Dark matter detection

Doping

- La, Y, Lu, Sc
- Tm, Ce, Nd

BaF₂-based materials

- Nanoparticles
- Ceramics
- Composites

