

LA-UR-21-31602

Approved for public release; distribution is unlimited.

Title: Radiochemical Separation and Spectrometric Analyses of Americium and Plutonium, Thesis Prospectus

Author(s): Mueller, Rebecca Jane

Intended for: Report

Issued: 2022-02-04 (rev.1)



Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

RADIOCHEMICAL SEPARATION AND SPECTROMETRIC ANALYSES OF AMERICIUM AND PLUTONIUM, Thesis Prospectus

Rebecca Mueller

February 11, 2022

LA-UR

Welcome




- Committee Members
 - Dr. Ralf Sudowe (chair)
 - Dr. John Auxier, II
 - Dr. Thomas Johnson
 - Dr. Sherrod Maxwell
 - Dr. Alexander Brandl
 - Dr. Thomas Borch
- Special thanks
 - Other Sudowite graduate students
 - Kim Bonilla
 - Dr. Stephanie Velardo

Outline

- Why now? (Current Views)
- History
- Experimental & Results
- Issues/Future Work

Why now? An Intro to Am

- Discovered in 1944 by Seaborg, et al. at the University of Chicago (schema to right)
- Importance
 - Am-241 is a by-product of Pu-239 production and decay, and is removed during Pu purification
 - Very useful due to its low cost, convenience, spectral purity and half-life
 - Largest use as a low-energy gamma-ray source



Why now? An Intro to Am (cont.)



Image: EPA, "Radionuclide Basics: Americium-241"

- Uses
 - Fire-detection applications
 - Important in energy security as it can be paired with Be to make neutron sources for well-logging
 - Gauging thickness, density, and radiographic measurements
 - Production of Cm-242, for the start-up of nuclear reactors

Why now? Why Am and Pu?

- C-AAC performs analyses for programmatic support at TA-55
 - Certain analysis are performed to measure concentrations of Pu and Am-241 in process solutions
 - Solutions contain high levels of Pu (e.g., g/L) and Am-241, with significant amounts of stable isotopes
- Pu-238, 239, 240, 242 contribute alpha activity.
 - Am-241, alpha emitter, overlaps with Pu-238 emission
 - Pu-239, a major component of weapons grade Pu, has a half-life of 24,000 yrs, in contrast, Am-241 has a shorter half-life of 400 yrs. This creates the situation where a small mass fraction of Am-241 can contribute a majority of the overall alpha activity.
 - For example 5.0% Am-241 by mass contributes 74.4% of the overall alpha activity.
- Currently two techniques are employed, depending on the level of Am-241 in the sample.

Why now? Separation Schemes

- Less Am-241 (Am-241 α activity $< 90\%$ overall α activity):
 - Dilute sample and divide into two portions
 - One portion, place on glass slip, to be counted in gas-proportional counter to determine total alpha
 - Second portion counted on an NaI(Tl) counter to determine Am-241 (59 keV gamma)
- More Am-241 (Am-241 α activity $> 90\%$ overall α activity):
 - Liquid-liquid separation to separate Pu and Am-241
 - Counts are done of the two resulting fractions
 - Avoids the situation where the difference between large alpha counts and NaI(Tl) count is dominated by relative uncertainties of both counts
- Liquid-Liquid separation generates a waste stream with limited disposition options.
 - Flammability issues with the organic phase (TTA/xylene)


Outline

- Why now?
- History (A Review of the Separation)
- Experimental & Results
- Issues/Future Work

History

- Thenoyltrifluoroacetone (TTA) has been used for the separation of rare earth metals since the late 1940s (1)
- TTA use for extracting actinides via liquid-liquid extraction has been used since the 1940s/1950s (2)

Group	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
Period	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
1	1 H																	2 He
2	3 Li	4 Be											5 B	6 C	7 N	8 O	9 F	10 Ne
3	11 Na	12 Mg											13 Al	14 Si	15 P	16 S	17 Cl	18 Ar
4	19 K	20 Ca	21 Sc	22 Ti	23 V	24 Cr	25 Mn	26 Fe	27 Co	28 Ni	29 Cu	30 Zn	31 Ga	32 Ge	33 As	34 Se	35 Br	36 Kr
5	37 Rb	38 Sr	39 Y	40 Zr	41 Nb	42 Mo	43 Tc	44 Ru	45 Rh	46 Pd	47 Ag	48 Cd	49 In	50 Sn	51 Sb	52 Te	53 I	54 Xe
6	55 Cs	56 Ba	* 71 Lu	72 Hf	73 Ta	74 W	75 Re	76 Os	77 Ir	78 Pt	79 Au	80 Hg	81 Tl	82 Pb	83 Bi	84 Po	85 At	86 Rn
7	87 Fr	88 Ra	* 103 Lr	104 Rf	105 Db	106 Sg	107 Bh	108 Hs	109 Mt	110 Ds	111 Rg	112 Cn	113 Nh	114 Fl	115 Mc	116 Lv	117 Ts	118 Og
			* 57 La	58 Ce	59 Pr	60 Nd	61 Pm	62 Sm	63 Eu	64 Gd	65 Tb	66 Dy	67 Ho	68 Er	69 Tm	70 Yb		
			* 89 Ac	90 Th	91 Pa	92 U	93 Np	94 Pu	95 Am	96 Cm	97 Bk	98 Cf	99 Es	100 Fm	101 Md	102 No		




(1) *Anal. Chem.* Vol. 29, No. 12, pp. 1767-1770 (1957)
 (2) *JACS.* Vol. 76, pp. 1982-1984 (1954)

History

- Solvent extraction
 - Distribution of a chemical species between two immiscible phases—typically aqueous and organic
 - Wide range of available reagent concentrations
 - Rapid and specialized separations
 - Distribution ratio given by $D^{SE} = \frac{c_{org}}{c_{aq}}$
 - Separation factor given by $SF_{AB} = \frac{D_A}{D_B}$


Modern Applications, edited by Frank Rösch, De Gruyter, Inc., 2016.
ProQuest Ebook Central,
<https://ebookcentral.proquest.com/lib/alamos/detail.action?docID=4644589>.
Friedlander, Gerhart, J.W. Kennedy, E.S. Macias, J.M. Miller. Nuclear and
Radiochemistry, 3rd Ed. 1981.

History




- TTA is selective for +4 oxidation state Pu
- Was originally done in benzene, known carcinogen

$$D_w = \frac{\text{Activity}_{\text{organic phase}}}{\text{Activity}_{\text{aqueous phase}}}$$



History

- Extraction Chromatography
 - Extractant material mounted on inert substrate
 - Aqueous solution passed through column material, typically inorganic acids or mineral acids
 - Effective for the separation of metal ions
 - Advantages of good selectiveness of extractants and multistages of chromatography




Conceptual section through an extraction chromatography (EXC) porous bead impregnated with the stationary, organic extracting phase, and the surrounding aqueous mobile phase.
E.P. Horwitz,
<http://www.eichrom.com/products/extraction.cfm>.

Outline

- Why now?
- History
- Experimental (Procedure) & Results
- Issues/Future Work

Experimental

1. Adjust pH of solution
2. Add $\text{NH}_2\text{OH HCl}$ for Pu(III)
3. 10 minutes in hot water bath, if needed
 1. Allow to cool (if hot water bath was used)
1. Add 2 M NaNO_2 for Pu(IV)
2. Add 0.5 mL TTA in xylene
3. Shake
4. Transfer aliquot to cover slip
5. Count when dry



FLOW DIAGRAM FOR DILUTION AND EXTRACTION PROCEDURES




Figure 1


Results

TABLE IV
EFFECT OF TEMPERATURE ON RATE CONSTANTS

Temperature, °C	k_1 , mol l. ⁻¹ sec ⁻¹	k_2 , sec ⁻¹	k_3 (K ₁)
25	7.04 ± 1.14	48.2 ± 9.1	2.14 ± 0.49
35	26.6 ± 3.5	135. ± 19.	2.10 ± 0.28
40	38.1 ± 2.9	347. ± 32.	2.49 ± 0.11
45	68.6 ± 7.0	720. ± 67.	2.82 ± 0.16

Barney, G S. KINETICS AND MECHANISM OF PLUTONIUM(IV)
REDUCTION BY HYDROXYLAMINE.. United States: N. p.,
1971. Web. doi:10.2172/4591822.

10 min in hot water bath v. 20 min
sitting



Results

- Rate limiting step is transfer of chelated Pu across aqueous-organic boundary layer
- Increase in stirring speed does lead to faster extraction
- Amount of stirring when hand shaking does agree with this (vortex mixing preferred)

Shake Time v. Recovery of Pu in TTA Xylene






FIG. 4: THE EXTRACTION OF Pu(IV) FROM THE AQUEOUS TO THE ORGANIC PHASE AS A FUNCTION OF TIME AT TWO SHAKING SPEEDS

Experimental


1. Control redox state of Pu-239 with hydroxylamine hydrochloride and sodium nitrite
2. Add 50 Bq Pu-239 to each sample
3. Contact 2 mL aqueous (of various acid concentrations) with 2 mL 0.5 M TTA in xylene
4. Allow to settle and remove 1 mL from each phase
5. Add scintillation cocktail and count via LSC

Results


Pu-239 avg Dw v. Nitric Acid Concentration




EFFECT OF HNO_3 CONCENTRATION ON THE EXTRACTION COEFFICIENT OF PU(IV) INTO 0.2 M TTA IN CCL_4 AT ROOM TEMPERATURE



Results




Results



Separation Factors

- Optimal separation is at 1M nitric acid
- Separation factors are higher at 0.5 and 3.5 M
- Likely due to the very very small Dw of Am at these concentrations



Nitric Acid Concentration / M	Dw: Am(standard deviation)	Dw: Pu (standard deviation)	SF
0.1	0.22(0.21)	2(0.74)	9.2
0.35	0.20(0.11)	2.88(0.63)	14.51
0.5	0.14(0.06)	7.98(1.64)	56.2
1	0.24(0.29)	9.68(0.47)	40.97
3.5	0.07(0.01)	4.97(0.71)	66.49
5	0.08(0.01)	0.94(0.23)	11.6


Conclusion

1. 1 M Nitric Acid worked best
2. TTA-xylene separation has separation factor >9 for Pu relative to Am at all concentrations tested
3. Effective separation

HOWEVER!....

Issues and Further Work

1. TTA has 3 fluorides
2. Low atomic numbers have higher cross-sections for (α, n) reaction
3. Fluorine undergoes (α, n) reaction leading to large neutron dose of the waste drum
4. Potential neutron dose to workers should be avoided when ever possible





Issues and Further Work

1. Xylene is a problematic waste constituent
2. Other solvents with TTA
 1. AmeriClear
 2. HistoClear
 3. Dodecane
 4. 1,3-diethylbenzene

Column Separation


- When Am is major constituent, use separation medium that selects for tetravalent ions, like TEVA (trialkyl, methylammonium nitrate or chloride)
- When Pu is major constituent, use separation medium that selects for trivalent ions, like DGA



Column Separation

- TEVA extractant structure shown (below)
- Experimental flow diagram shown (right)
- Control redox state of Pu with ascorbic acid and sodium nitrite

Trialkyl, methylammonium
nitrate (or chloride)



$\text{R} = \text{C}_8\text{H}_{17} \text{ and } \text{C}_{10}\text{H}_{21}$






Figure 1

TEVA Column Extraction



5. Add 3mL 0.02M HNO_3 -0.005M HF
4. Rinse column with 3M Nitric Acid x2
3. Rinse sample tube with 3M Nitric Acid
2. Load Sample (Pu & Am) onto Column
1. Nitric Acid precondition




1. Collection vessel for disposal
2. Collect (Am) in solution, not retained
3. Collection vessel for disposal
4. Collection vessel for disposal
5. Collect (Pu) in new vial, plate fraction for α counting

Results: Initial TTA/TEVA Comparison

- Same sample, two analyses
- Am determined concentration was equal in both
- Pu determined concentration was higher with TEVA

Comparison of Recovered Activity for TTA and TEVA Extractions



Outline





- History
- Why now?
- Experimental
- Results
- Issues/Future Work (Moving Forward)

Issues and Further Work

1. Xylene is a problematic waste constituent
 - Other solvents???
2. Fluorine undergoes (α , n) reaction leading to large neutron dose of the waste drum
3. PuTTA complex on UV-Vis (next slide)
4. Other extractants
 - Dihexyl Hexanamide (DHHA)
 - Dihexyl Octanamide (DHOA)
 - Tetra- or tri-alkylcarbamides

Issues and Further Work: Other Extractants

- Dihexyl Hexanamide (DHHA) (left)
- Dihexyl Octanamide (DHOA) (middle)
- Tetra- or tri-alkylcarbamides (basic carbamide structure shown right)



Other Extractants, cont.

Figure 1. Distribution ratios as functions of acidity for hexavalent actinides Pu(VI) (), Np(VI) (), and U(VI) () using 1 M monoamide/Exxsol D60. Lines are plotted as a visual guide




Figure 4. Distribution ratios as functions of acidity for tetravalent actinides Pu(IV) (), Np(IV) (), and Th(IV) () using 1 M monoamide/Exxsol D60. Lines are plotted as a visual guide




Figure 3. Distribution ratios as functions of acidity for Pa(V) 1 M DEHBA/dodecane () vs. Np(V) DEHBA (), DEHBA (), DHOA (), and DHHA () using 1 M monoamide/Exxsol D60. Lines are plotted as a visual guide





Figure 5. Distribution ratios as functions of acidity for extracting Am(III) (left) and Pu(III) with 0.43 M hydroxylamine (right). DEHBA (), DEHBA (), DHOA (), and DHHA () using 1 M monoamide/Exxsol D60. Lines are plotted as a visual guide



Gogolski, Jarrod M, and Mark P Jensen. "Using N,N-dialkylamides for Neptunium Purification from Other Actinides for Space Applications." *Separation Science and Technology* 56.16 (2021): 2775-788. Web.



Americium Isotopics

- TIMS: Thermal Ionization Mass Spectroscopy
 - Very Expensive but provides
 - Pu Assay values
 - Am isotopic ratio values (^{241}Am and ^{243}Am)
 - Requires chemical separation similar to previous work
 - Also requires a more lengthy turn-a-round time.




Americium Isotopics

Nuclide	Decay Method	Half-life	Daughter	Emax, β (MeV)	E, α (MeV)	E, γ (MeV)
Sr-89	β, γ	50.5 d	Y-89	1.5	N/A	N/A
Sr-90	β	28.8 y	Y-90	0.546	N/A	N/A
Y-89	STABLE	N/A	N/A	N/A	N/A	N/A
Y-90	β	64.0 h	Zr-90	2.28	N/A	N/A
Am-241	α, γ	432.6 y	Np-237	N/A	5.486	0.0595
Am-243	α, γ	7364 y	Np-239	N/A	5.275	0.0747
Np-237	α	2.144×10^6 y	Pa-233	N/A	4.788	0.0294
Np-239	β, γ	2.356 d	Pu-239	0.438	N/A	0.1061



Proposed Detection Methods

- ORTEC Octete Spectrometer



Silicon Detectors

450 mm²: 17-21 keV, FWHM

50 mm²: 12-14 keV, FWHM resolution



Photos courtesy of ORTEC, Inc.

References

1. Modern Applications, edited by Frank Rösch, De Gruyter, Inc., 2016. ProQuest Ebook Central, <https://ebookcentral.proquest.com/lib/alamos/detail.action?docID=4644589>.
2. GREGORY R. CHOPPIN, JAN-OLOV LILJENZIN, JAN RYDBERG, Editor(s): GREGORY R. CHOPPIN, JAN-OLOV LILJENZIN, JAN RYDBERG, Radiochemistry and Nuclear Chemistry (Second Edition), Butterworth-Heinemann, 1995, ISBN 9780750623001, <https://doi.org/10.1016/B978-0-7506-2300-1.50006-3>. (<https://www.sciencedirect.com/science/article/pii/B9780750623001500063>)
3. Ambrosi, R.M., Williams, H., Watkinson, E.J. et al. European Radioisotope Thermoelectric Generators (RTGs) and Radioisotope Heater Units (RHUs) for Space Science and Exploration. *Space Sci Rev* 215, 55 (2019). <https://doi.org/10.1007/s11214-019-0623-9>
4. Annika Tovedal, U. Nygren, H. Ramebäck. Determination of ^{90}Sr in preparedness: Optimization of total analysis time for multiple samples. *Journal of Radioanalytical and Nuclear Chemistry* (2008) 276:357-362.
5. Annika Tovedal, U. Nygren, H. Ramebäck. Methodology for determination of ^{89}Sr and ^{90}Sr in radiological emergency: I. Scenario dependent evaluation of potentially interfering radionuclides. *Journal of Radioanalytical and Nuclear Chemistry* (2009) 282:455-459.
6. Annika Tovedal, U. Nygren, P. Lagerkvist, A. Vesterlind, H. Ramebäck. Methodology for determination of ^{89}Sr and ^{90}Sr in radiological emergency: II. Method development and evaluation. *Journal of Radioanalytical and Nuclear Chemistry* (2009) 282:461-466.
7. U.S. Department of Energy Isotope Program <https://isotopes.gov>
8. "Considerations for Use of Am-241 for Heat Source Material for Radioisotope Power Systems" Stephen G. Johnson. January 2017. INL/EXT-16-40336-Revision-1
9. Plutonium Handbook, 2nd Edition. Editors: Davic L. Clark, David A. Geeson, Robert J. Hanrahan, Jr. American Nuclear Society, 2019.
10. "Exploration of Pu and Am extraction studies with TTA/xylene" RJ Mueller, et al. Fall American Chemical Society Meeting, Conference Proceedings. August 25, 2021. LA-UR-21-28404.
11. Irwin L. Hinds. "A Comparison of Three Xylene Substitutes." *Laboratory Medicine*: Vol. 17 (12). December 1986.
12. "Hazardous Waste Management" from the United States Environmental Protection Agency. www.epa.gov/hw/defining-hazardous-waste-listed-characteristic-and-mixed-radiological-wastes#PandU
13. "DGA Resins" from Eichrom Technologies, LLC. www.eichrom.com/eichrom/products/dga-resins/
14. E. P. Horwitz, D. R. McAlister, A. H. Bond & R. E. Barrans Jr (2005) Novel Extraction of Chromatographic Resins Based on Tetraalkyldiglycolamides: Characterization and Potential Applications, *Solvent Extraction and Ion Exchange*, 23:3, 319-344, DOI: 10.1081/SEI-200049898
15. "TEVA Resin" from Eichrom Technologies, LLC. www.eichrom.com/eichrom/products/teva-resin/

References, cont.

16. Horwitz, E.; Dietz, M.; Chiarizia, R.; Diamond, H.; Maxwell, S.; and Nelson, M.; Separation and Preconcentration of Actinides by Extraction Chromatography Using a Supported Liquid Anion-Exchanger—Application to the Characterization of High-Level Nuclear Waste Solutions, *Analytica Chimica Acta*, 1995, 310 (1), 63-78.
17. Maxwell, S., Rapid Actinide Separation Methods, *Radioactivity and Radiochemistry*, 1998, 8, 36-44.
18. Samantha K. Cary, Maksim Livshits, Justin N. Cross, Maryline G. Ferrier, Veronika Mocko, Benjamin W. Stein, Stosh A. Kozimor, Brian L. Scott, and Jeffrey J. Rack. Advancing Understanding of the +4 Metal Extractant Thenoyltrifluoroacetate (TTA-); Synthesis and Structure of MIVTTA4 (MIV = Zr, Hf, Ce, Th, U, Np, Pu) and MIII (TTA)4- (MIII = Ce, Nd, Sm, Yb). *Inorg. Chem.* 2018, 57, 7, 3782–3797.
19. Kurosaki, Hiromu, Rebecca J Mueller, Susan B Lambert, and Govind R Rao. "Alternate Method of Source Preparation for Alpha Spectrometry: No Electrodeposition, No Hydrofluoric Acid." *Journal of Radioanalytical and Nuclear Chemistry* 311.1 (2017): 323-29. Web.
20. Live Chart of Nuclides, IAEA. Web.
21. Krmpotić, Matea, Martina Rožmarić, and Ljudmila Benedik. "Investigation of Key Factors in Preparation of Alpha Sources by Electrodeposition." *Applied Radiation and Isotopes* 136 (2018): 37-44. Web.
22. Gupta, K. K, V. K Manchanda, M. S Subramanian, and R. K Singh. "N,N-Dihexyl Hexanamide: A Promising Extractant for Nuclear Fuel Reprocessing." *Separation Science and Technology* 35.10 (2000): 1603-617. Web.
23. McCann, Kevin, Sergey I Sinkov, Gregg J Lumetta, and Jenifer C Shafer. "Inner: Versus Outer Sphere Metal-monoamide Complexation: Ramifications for Tetravalent & Hexavalent Actinide Selectivity." *New Journal of Chemistry* 42.7 (2018): 5415-424. Web.
24. K.K. Gupta , V.K. Manchanda , S. Sriram , George Thomas , P.G. Kulkarni & R.K. Singh (2000) THIRD PHASE FORMATION IN THE EXTRACTION OF URANYL NITRATE BY N,N-DIALKYL ALIPHATIC AMIDES, SOLVENT EXTRACTION AND ION EXCHANGE, 18:3,421-439, DOI: 10.1080/07366290008934690
25. Gupta, K. K, V. K Manchanda, M. S Subramanian, and R. K Singh. "Thermodynamics of Extraction of Uranium(VI) and Plutonium(IV) with Some Long-chain Aliphatic Amides." *Radiochimica Acta* 85.3 (1999): 103-06. Web.
26. Gogolski, Jarrod M, and Mark P Jensen. "Using N,N-dialkylamides for Neptunium Purification from Other Actinides for Space Applications." *Separation Science and Technology* 56.16 (2021): 2775-788. Web.
27. Friedlander, Gerhart, J.W. Kennedy, E.S. Macias, J.M. Miller. *Nuclear and Radiochemistry*, 3rd Ed. 1981.
28. Johnson, Thomas. *Introduction to Health Physics*, 5th Ed. 2017.
29. Nelson S. Stalnaker & TJ Beugelsdijk
30. Nelson S. Stalnaker & TJ Beugelsdijk. Robotic Sample Preparation for Radiochemical Plutonium and Americium Analyses. LA-UR—85-3444 DE86 000801.
31. "Applications of Liquid Scintillation Counting." Ralf Sudowe, Webinar by National Analytical Management Program (NAMP), Series 2: Environmental Radiochemistry/Bioassay Series, Number 11. Recorded May 22, 2014.


Acknowledgements

- Colorado State University
- Los Alamos National Laboratory
- Dr. John Auxier II & Dr. Ralf Sudowe
- Dr. Stephanie Velardo, Kimberly Bonilla, Nell Carver, Robert Jump
- Funding:
 - National Isotope Development Center for funding on the TIMS analysis
 - MR&R (Material Recycle and Recovery), NASA, and Pu Sustainment for funding on radiochemical methods.
 - Pu Sustainment (NA-10)
 - Seaborg Fellowship

Questions?


Why now?

- Purity
 - For the majority of these applications, Am purity needs to have > 99% purity (e.g. < 1% Pu and other elements)
 - The feedstock for Am typically doesn't contain other elements so Pu is the major impurity
 - Isotopically needs to be greater than 99% Am-241.
 - Needs to be greater than 95% AmO₂ by weight.
- Radiochemical separation?
 - Current is TTA/xylene



UV-Vis Results

• UV-Vis Spectrometer for Am/Pu Assay





Am Transitions

- (${}^7F_0 \rightarrow {}^5L_6$) @ 503 nm
- (${}^7F_0 \rightarrow {}^5L_6$) @ 806 nm

Pu⁴⁺ Transitions

- 470 nm
- 653 nm
- Thus no interference with Am³⁺



Single-crystal UV-vis-NIR absorption spectra from UIV(TTA)4 (top, red), NpIV(TTA)4 (middle, green), and PuIV(TTA)4 (bottom, orange) crystals.