

LA-UR-22-20888

Approved for public release; distribution is unlimited.

Title: P. putida Learn Experiments – TrRNA library – LANL & ANL

Author(s): Davison, Steffi Angela
Kern, Theresa Lynn
Dale, Taraka T.

Intended for: Meeting with collaborators

Issued: 2022-02-01



Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by Triad National Security, LLC for the National Nuclear Security Administration of U.S. Department of Energy under contract 89233218CNA000001. By approving this article, the publisher recognizes that the U.S. Government retains nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this article as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

P. putida Learn Experiments – TrRNA library – LANL & ANL

Steffi Davison, Theresa Kern, Naresh Panday, Claire Sanders, Bert Huttanus, Gyorgy Babnigg, Phil Laible, Rose Wilton, Debbie Hanson, Ramesh Jha, Scott Hennelly, Taraka Dale

31 January, 2022





Background and context

ABF libraries in past

- *Cis*-repressor library (LANL, Scott Hennelly and Naresh Pandey)
- Biosensor libraries (LANL Ramesh)
- RB-Tn-Seq (ORNL, Adam Guss)
 - Adam's TnSeq library was sorted a couple of times
 - Conceptually this should work, but there were technical issues that we didn't completely work through
 - low [DNA] led to insufficient sequence depth

In progress ABF Libraries + point of contact


- IMRE library (LANL, Theresa Kern)
- RB-Tn-Seq library (NREL, Caroline Amendola)
- CRISPRi library (Oakridge, Carrie Eckert)



How does IMRE library differ from other *putida* work?

- How is this different from the ongoing work?
 - No ALE or random mutagenesis
 - So far, sensor/sorting work has focused only on ‘top performers’
 - Here, we will also collect ‘poor and intermediate’ performance’ data: important for LEARN to bound what doesn’t work
 - Ideally this work will also give strains that increase TRY, but the main goal is to give **LEARN a richer dataset to improve Learn outcomes for increasing TRY**
- How is this different from using a transposon or other KO library?
 - Conceptually the same! But practically it’s different
 - Here we will use a **plasmid-based approach**, which is easier for DNA isolation & sequencing, AND we’ve already done it before
 - Also interested in seeing how **knockdowns perform versus knockouts**


IMRE Riboregulator RNA structures (Naresh Pandey and Scott Hennelly)



- TrRNA, a Twist library of all *P. putida* genes assembled
- Knockdown all ~5500 genes in *P. putida*

Riboregulators controls RBS accessibility, in turn, allows tunable protein expression control of a target mRNA '*in trans*'

Screening with biosensor work (Ramesh Jha and Taraka Dale)



Use fluorescent biosensors & flow cytometry to rapidly screen and isolate populations of strains with varied phenotypes




Hypothesis: Adding “negative”/”poor performance” data to the Learn models will better inform what makes a good performer, leading to improved predictions for increased TRY.


The Learn team has been saying this for a while, but experimentally the *Putida* team has thus far only focused on generating and characterizing and better strains, we have not placed resources on looking at poor performers

Goal: Provide the Learn team with a richer dataset, resulting from isolating and characterizing *P. putida* strains with a range of phenotypes for muconate production, including poor performers.


General experimental workflow




General experimental workflow plan



Lag in GFP signal
vs
Drop out rate



Define populations
to be collected by
“gating”
= 2 mil. cells per
population




GFP = phenotype

General experimental workflow plan

NGS sample preparation



Deep sequences of sorted populations




Enriched sequence analysis




- Plasmid sequencing of each sorted population to look for enriched sequences that can be correlated to each phenotype
- Identify genes of interest from enriched pools
- Finally, validate the genes in 'clean' background

Outgrowth and PCR Bias Investigation

1 step barcode/adaptor pcr





Library Deep Sequencing using Illumina MiSeq Technologies



Challenges with barcoding, library pool prep

- Plasmid 1 Step PCR: Meant to simplify and expedite process


- Long primers : Fwd_illumina_ppIMRE_A501 AATGATACGGCGACCACCGAGATCTACAC_{TGAACCTT}
ACACTCTTTCCCTACACGACGCTCTTCCGATCT_{CGC GGC ACA CTC TAG AGT TT}
- **Presence of primer dimers/nonspecific bands and low target amplicon yields:** time consuming and tedious library prep.
 - Multiple PCRs with high cycle number
 - Gel band excision cleanup



Troubleshooting/protocol discussion with LBL

- Preferential formation of dimers/concatamers
- Low complexity of amplicons
 - Difficult for MiSeq analysis

Index (Variable region)




Next steps for IMRE library


- Repeat Library growth/sorting/ pcr bias with “more sensitive” muconate sensor. -This should reduce GFP signal lag time
- Improved Miseq library pool prep with guidance from LBL:
 - Several options being explored:
 - 1 step pcr with altered long primers: removal of extra “G: residue in P7 adapter
- i7 index - Read 2 primer region. (Can inhibit template binding)

5' CAAGCAGAAGACGGCATAACGAGAT- index - **G** TGAAGTGG
AGTTCAGACG TGTGCTCTTCCGATCT---AAC AGA TAA AAC
GAA AGG CCC 3'
 - 2 step PCR with short primers/limited cycle numbers
 - 10 cycle of PCR + Illumina adapter
 - Or instead of 2nd pcr, adapter ligation starting with 100 ng from 1st pcr
 - Improved Thermocycling protocols
 - Magnetic bead cleanup


General experimental workflow plan



Lag in GFP signal
vs
Drop out rate



Define populations
to be collected by
“gating”
= 2 mil. cells per
population



GFP = phenotype

RB_Tn-Seq library (NREL)


Muconate biosensors:

- **CA089 and CA090** - CJ781 (muconate from aromatics)
 - M9 minimal media with 20 mM p-coumarate and 10 mM glucose
- **CA091 and CA092** – KH083 (muconate from sugars)
 - M9 minimal media with 30 mM glucose and 15 mM xylose

BKA biosensors:


- **CA093 and CA094** – AW311 (BKA from aromatics)
 - M9 minimal media with 20 mM p-coumarate and 10 mM glucose
- **CA095 and CA096** – GR038 (BKA from sugars)
 - M9 minimal media with 30 mM glucose and 15 mM xylose

Biosensors: more sensitive vs less sensitive




Monoclonal strain growth curves

10mM glucose/20mM p-coumarate




30mM glucose/15mM xylose



KT2440 Library collapsing experiment


KT2440B library collapse




Purpose: to track library changes throughout growth.

- M9, 50mM glucose
- Performed in triplicate
- Samples taken for each time point shown.
- *MiSeq analysis pending protocol revisions.*

General experimental workflow plan




Cell sorting plan, added detail



- Sort all 3 triplicates 4 ways in parallel = 15 sorts & 60 samples - 2 full days of sorts
- ~1 x 10⁶ cells collected per sort


transrepressor length = 100 (200+ with terminator)	transrepressors -- 1	Actuals
population desired	6000designed per gene	5500
sequencing coverage	300x	182
need to collect	1.80E+06cells	1.00E+06
FACS optimal event rate	8000cells/sec	8000
# of parallel sorts	4	4
Fraction of population gated	50%population gated	64%
cells per sort (cells/m)	1000tube	1280
	1800sec	781
	30min of sorting	13

FACS of trans repressor library SD3




16h

- Aria data with four gates
- Actual 1M cells collected per gate
- How collected (tubes)
- 55% of population collected




40h

FACS of trans repressor library LC236




16h

- Aria data with four gates
- Actual 1M cells collected per gate
- How collected (tubes)
- 57% of population collected



40h

Example: analytical flow cytometry of sorted populations




- Overnight re-growth and analysis shows sorted populations remain (mostly) distinctive
- Plasmid miniprep
- PCR barcode+ illumina adapters
- Sequencing

Example: sequence selection in sorted populations

- Sequences selected “by eye/hand”

Deep-sequence sorted pools



Comments from meeting

- LC224-B FACs – re-draw the gates!
- Sort NP162-B when library is completed (with and without gates)
- Amplicon sequencing vs plasmid sequencing
 - Go ahead with amplicon sequencing for the samples
 - Send pooled samples (SD3-B, LC224-B, LC236-B and LC238-B) and some separate samples (NP162 samples)
 - Send PCR cycles (10-40 cycles) for sequencing (determine PCR bias)
 - LBL can do the bulk 100+ sequencing samples, Argonne can do a few sequencing samples
 - Work with Josh McCaully, not Garima (left LBL)
 - Phil will check with Gyorgy about cycle number (usually less than 15 cycles, no PCR reach exponential phase yet)
 - Nathan thinks they will sequence whole 450bp
- Productivity vs yield muconate (perform a time-course with FI, muconate yield and OD of the selected strains)
- Share on a Google Drive (Jeremy will send a place to port if over) - old sequencing data to either Jeremy Zucker and Hector
- Jeremy will set up a meeting with Craig Bakker early Oct for stats analysis on old data)
- Send samples mid-October to LBL
- If we need a re-sort, then we can quickly get through it as we have a setup pipeline (milestone)
- Sort with gaps between the gates (maybe do NP162-B)
- Hector needs to get together Ramesh to get the background for the experiment
- Next meeting will revolved around Hector's understandings and aims of the project (maybe combination of genes in the future, maybe knockout strains and in combination)